

Introduction

There are many challenges associated with metagenome assembly:

- the presence of multiple species
- uneven and unknown species abundances
- conserved genomic regions shared across species
- strain-level variation within species

PacBio HiFi sequencing produces highly accurate long reads (>Q20, >99% accuracy) which provide major advantages for metagenome assembly. New metagenome assembly algorithms have been developed specifically for HiFi reads, including **hifiasm-meta**¹ and **metaMDBG**². These methods can reconstruct full metagenome-assembled genomes (**MAGs**) for many higher abundance species.

Metagenome assembly of soil microbiomes has been historically difficult using short reads. The combination of high species diversity and ultra-low relative abundances poses a challenge, and requires a higher sequencing depth to achieve success with long reads. Here, we demonstrate that the amount of HiFi data from the high-throughput PacBio Revo system is sufficient to assemble high-quality MAGs in complex microbiomes, including soil.

Methods

PacBio HiFi sequencing

A soil core was obtained from a northern California wetland and sampled along six depths. The six samples were prepped and sequenced on a PacBio Revo system using three 25M SMRT Cells. Sequencing resulted in a total of 20.6 million HiFi reads, 196 Gb total data, and a median read QV of 41 (representing >99.99% accuracy).

Sequencing	HiFi reads (million)	Total data (gigabases)	Average read length (kb)	Median QV
SMRT Cell 1	8.21 M	79.83 Gb	9.7 kb	Q40
SMRT Cell 2	6.65 M	57.15 Gb	8.6 kb	Q43
SMRT Cell 3	5.76 M	58.90 Gb	10.2 kb	Q41

Metagenome assembly and postprocessing

The combined sequencing dataset of 20.6M HiFi reads was assembled using **hifiasm-meta** and **metaMDBG**, and each contig set was processed using the PacBio **HiFi-MAG-Pipeline**. The complete analysis workflow is shown visually in Figure 1.




Figure 1. Visual overview of methods used for assembly and post-processing.

MAGs were categorized based on quality scores from CheckM2:

- Medium-quality (MQ): >50% completeness and <10% contamination
- High-quality (HQ): >90% completeness and <5% contamination
- Single-contig high-quality (SC-HQ): same criteria as high-quality, but the MAG also consists of one contig

Results

HiFi metagenomics routinely produces fully resolved MAGs

Both **hifiasm-meta** and **metaMDBG** produced over 350 single-contig, HQ-MAGs directly from the assembly step, which often appeared as circular genomes in the assembly graphs (Fig. 2).




Figure 2. A partial view of the hifiasm-meta assembly graph for the combined soil dataset. The graph reveals many large circular contigs (0.5–13 Mb) produced directly from assembly.

Over 500 high-quality MAGs assembled from a soil sample

- **hifiasm-meta** produced over 1,200 total MAGs, with 550 HQ-MAGs (391 are single-contig; Fig. 3)
- **metaMDBG** resulted in 1,000 total MAGs, with over 500 HQ-MAGs (355 are single-contig; Fig. 3)




Figure 3. Counts of MAGs assigned to different quality categories.


Assembling genomes for hundreds of uncultured species

None of the MAGs could be assigned to the species level by GTDB-Tk. Recovered **1,097 bacterial** MAGs from 51 phyla (513 HQ-MAGs; Fig. 4). Recovered **135 archaeal** MAGs from 8 phyla (50 HQ-MAGs; Fig. 5).

Figure 4. Phylogeny of 1,097 bacterial MAGs that were assembled.


The most abundant phyla include Acidobacteriota, Proteobacteria, Patescibacteria, and Chloroflexota.

Outside edges of the phylogeny show whether a MAG is classified as HQ (navy), SC-HQ (red), and the relative depth of coverage (light blue).



Archaea

- Nanoarchaeota
- Aenigmataarchaeota
- Iainarchaeota
- Micrarchaeota
- Hadarchaeota
- Thermoplasmata
- Halobacteriota
- Thermoproteota



The 50 HQ-MAGs here are in many cases the first high-quality genomes available for a given taxonomic group.

Outside edges of the phylogeny show whether a MAG is classified as HQ (navy), SC-HQ (red), and the relative depth of coverage (light blue).

Long-read sequencing recovers genomes with diverse properties

- Size range across HQ-MAGs was 0.5–12Mb (median 4Mb)
- HQ-MAGs display a range of 35–73% GC content (Fig. 6)
- HQ-MAGs displayed 4X–400X depth of coverage (median 17X)




Figure 6. Characteristics of the 1,232 MAGs that were recovered by hifiasm-meta, including the genome size (x-axis) and average depth of coverage per genome (y-axis). Each point represents an individual MAG, and they are color-coded based on their estimated percent GC content. Three high coverage genomes were excluded from the plot (400-600X).

Conclusions

- PacBio HiFi sequencing offers major advantages for metagenome assembly, particularly for difficult environmental samples.
- Single-contig HQ-MAGs are routinely assembled with HiFi reads.
- HiFi sequencing is effective for obtaining large numbers of high-quality MAGs from uncultured species in complex microbiomes.

All PacBio metagenomics workflows are open-source and publicly available on **GitHub**:



References

1. Feng et al. 2022. Metagenome assembly of high-fidelity long reads with hifiasm-meta. *Nature Methods*, 19: 671–674.
2. Benoit et al. 2024. High-quality metagenome assembly from long accurate reads with metaMDBG. *Nature Biotechnology*, <https://doi.org/10.1038/s41587-023-01983-6>
3. Chklovskii et al. 2023. CheckM2: a rapid, scalable and accurate tool for assessing microbial genome quality using machine learning. *bioRxiv*, <https://doi.org/10.1101/2022.07.11.499243>
4. Kang et al. 2019. MetaBAT 2: an adaptive binning algorithm for robust and efficient genome reconstruction from metagenome assemblies. *PeerJ*, 7: e7359.
5. Pan et al. 2023. SemiBin2: self-supervised contrastive learning leads to better MAGs for short- and long-read sequencing. *Bioinformatics*, 39: 121–129.
6. Sieber et al. 2018. Recovery of genomes from metagenomes via a dereplication, aggregation and scoring strategy. *Nature Microbiology*, 3: 836–843.
7. Chaumeil et al. 2019. GTDB-Tk: a toolkit to classify genomes with the Genome Taxonomy Database. *Bioinformatics*, 35: 1925–1927.