

Targeted long-read sequencing of native DNA for genetic disease diagnostic and screening research

Sarah B. Kingan, Guilherme De Sena Brandine, Jeff Zhou, Tom Mokveld, Jocelyne Bruand, Egor Dolzhenko, Michael A. Eberle, Duncan Kilburn

PacBio, 1305 O'Brien Drive, Menlo Park, CA 94025

Background

Despite advances in DNA sequencing, the causal mutations of many human diseases remain challenging to characterize with standard methods. The adoption of long-read sequencing can improve and simplify sequencing and analysis of these complex regions. The ability to detect methylation enables simultaneous generation of sequencelevel and epigenetic data. We present new panel content for the PureTarget assay, which uses the CRISPR/Cas9 system to generate targeted sequencing libraries of native DNA, sequenced with long and accurate HiFi sequencing.

Scalable amplification-free workflow



Figure 1. The PureTarget technology is a robust amplification-free approach to generate long-read HiFi sequencing libraries for panel targets. Starting with high molecular weight DNA, the workflow employs Cas3 and pairs of guide RNAs to target loci ranging in size from 3 – 14 kb.

Automated and manual workflows

	Automated	Manual
Input DNA	1 – 1.5 µg	1 – 4 µg
Sample types	Blood, saliva, cell lines	Nanobind cell line, blood (RBC-lysis)
Library prep time	16hr for 96 samples	8hr for 24 samples
Samples / SMRT cell	96 on Revio 384 per run	8 – 48 on Revio and Vega

Coverage in automated 96-plex prep





Sample (N=96) Nanobind whole blood random donor

Figure 2. 96 samples of 1 µg of random donor blood extracted with Nanobind HT CBB kit from whole blood (protocol 102-573-Solution of the second carrier is 494-fold. Fragment size on RE panel is more uniform (3.3 – 5.6kb) than carrier panel (3.9 – 14kb) which results in more variable target coverage on carrier panel.

Carrier screening panel

F8 intron 1 Hemophilia Inversion	
FXN Friedrich's ataxia Repeat Expan	sion
FMR1 Fragile-X disease (FXS) Repeat Expan	sion
CYP21A2, Congenital adrenal hyperplasia, Small variants TNXB Classical-like Ehlers-Danlos syndrome copy number	and er
HBA1/2 Alpha thalassemia Deletion	
GBA Gaucher Small variants copy number	and er
SMN1/2 Spinal muscular atrophy Small variants copy number	and er
ARX Early-infantile epileptic Repeat Expan Partington syndrome (PRTS) (N=2)	sion
HBB Beta thalassemia Small variar	its
RPGR X-linked retinitis pigmentosa Small variar	its
AFF2 Fragile X syndrome, FRAXE type Repeat Expan	ision

CYP21A2 - CNV and fusions





Figure 3. A pair of gRNAs capture ~7.4 kb spanning the full CYP21A2 gene and 13 exons of TNXB on the 3' end. The gRNAs also cut the segmental duplication containing CYP21A1P and TNXA. HG03540 reference sample was sequenced on Revio + SPRQ and analyzed with paraphase 3.1 in the PureTarget carrier pipeline. Five haplotypes are captured showing 2 gene copies, 2 pseudogene copies, and a fusion allele

HBA1/2 - large deletion



Figure 4. A pair of gRNAs captures ~9.5kb spanning HBA1 and HBA2. The 3p7 deletion is neatly captured showing one copy of HBA2 and 2 copies of HBA1.

PureTarget carrier pipeline (PTCP)



Figure 5. Schematic of PureTarget carrier pipeline for variant calling. Pipeline will be available on DNA Nexus Marketplace and for flexible deployment on cloud or on-prem server via github repository.

Repeat expansion panel 2.0

Targets	Disease
ATN1, ATXN1, ATXN2, ATXN3, ATXN7, ATXN8, ATXN10, CACNA1A, PPP2R2B, TBP BEAN1, DAB1, FGF14, NOP56, ZFHX3	Spinocerebellar ataxia (SCA)
FMR1	Fragile-X disease (FXS)
AFF2	Fragile X syndrome, FRAXE type
AFF3	Intellectual disability associated with fragile site FRA2A
C9orf72	Frontotemporal dementia (FTD), amyotrophic lateral sclerosis (ALS)
FXN	Friedreich ataxia (FRDA)
RFC1	Cerebellar ataxia, neuropathy, and vestibular areflexia syndrome (CANVAS)
NOTCH2NLC	Neuronal intranuclear inclusion disease, Alzheimer disease and parkinsonism phenotype (NIID)
DMPK, CNBP	Myotonic dystrophy (DM)
HTT	Huntington disease (HD)
JPH3	Huntington's disease-like type2 (HDL2)
TCF4	Fuchs endothelial corneal dystrophy 3 (FECD3)
AR	Kennedy Disease, Spinal and bulbar muscular atrophy, (SBMA)
PABPN1	Oculopharyngeal muscular dystrophy (OPMD)
ABCD3, GIPC1, LRP12, RILPL1	Oculopharyngodistal myopathy (OPDM)
HOXD13	Syndactyly (SD5)
PHOX2B	Congenital central hypoventilation syndrome (CCHS)
PRNP	Creutzfeldt-Jakob disease (CJD)
CSTB	Progressive Myoclonic Epilepsy Type 1 (EPM1) Unverricht-Lundborg Disease (ULD)
	Construction of the second s

Heterozygous expansion in FGF14



Figure 6. 2 µg of DNA, extracted from blood with QIAsymphony, was manually prepared with Pure Target repeat expansion panel 2.0, sequenced on Revio + SPRQ in an 8plex, and analyzed with TRGT v3.0.0. Expected genotype from rpPCR.

Heterozygous expansion in FMR1



Figure 7. DNA from Coriel reference sample NA07537 was manually prepared with Pure Target repeat expansion panel 1.0, sequenced on Sequel IIe, and analyzed with TRGT v3.0.0. Expected genotype (28-29 and >200 repeats) from southern and PCR analysis. Observed genotype is 27 CGG, 2 AGG and 332 CGG and 1 AGG.

Conclusions

- PureTarget enables scalable targeted HiFi sequencing of native DNA
- Panels available for repeat expansion disease and carrier screening (rare disease)
- Manual and automated preps available
- PureTarget carrier pipeline available for genotyping with TRGT and Paraphase

Acknowledgements and References

The authors wish to thank Johannes Zschocke. Katharina Lackner and Katharina Herzog from Medical University Innsbruck and Dale Muzzey, Ben Obadia and Jack Roach from Myriad Women's Health for kindly sharing samples.

Dolzhenko et al., (2024). Nature Biotech., 1-9. https://doi.org/10.1038/s41587-023-02057-3. Chen X, Baker D, Dolzhenko E, et al., (2025). Nature Communications. doi:10.1038/s41467-025-57505-2