Preparing Kinnex[™] libraries using Kinnex single-cell RNA kit



Procedure & checklist

Before you begin

This procedure describes the workflow for constructing single-cell Kinnex libraries from 10x Chromium 3' or 5' cDNA using the *Kinnex single-cell RNA kit* (103-072-200) for library prep and sequencing on PacBio[®] Sequel[®] II, Sequel IIe, and Revio[™] systems.


This kit is intended for use with single-cell cDNA generated using the 10x Chromium Next GEM Single Cell 3' kit v3.1 or 10x Chromium Next GEM Single Cell 5' kit v2, standard throughput. It has not been tested for use on low throughput (LT) or high throughput (HT) kits which are currently unsupported.

Overview	
Samples per kit	12
Workflow time	3 days for up to 12 samples

cDNA input	
	>15 ng per library
Quantity	cDNA concentration should be >1ng/ μ L with up to 15 μ L in volume. See step 2.1 for 10x cDNA input requirement.
Average segment lengths	500-1,000 bp
Average 16-segment array lengths	10-15 kb



Workflow





Required materials and equipment

DNA sizing	
Femto Pulse system	Agilent Technologies, Inc. M5330AA
Femto Pulse gDNA 165 kb analysis kit	Agilent Technologies, Inc. FP-1002-0275
Agilent 2100 Bioanalyzer system	Agilent Technologies, Inc. G2939BA
Agilent High Sensitivity DNA Kit	Agilent Technologies, Inc. 5067-4627
DNA quantitation	
Qubit fluorometer	ThermoFisher Scientific Q33238
Qubit 1X dsDNA HS assay kit	ThermoFisher Scientific Q33230
SMRTbell [®] library preparation	
Kinnex™ single-cell concatenation kit	PacBio [®] 103-242-000*
Kinnex™ capture beads kit (12 rxn)	
(Includes Kinnex™ capture beads, Kinnex™ bead binding buffer, Kinnex™ bead washing buffer)	PacBio [®] 103-076-000*
SMRTbell® cleanup beads	PacBio [®] 102-158-300*
Elution buffer (50 mL)	PacBio [®] 101-633-500*
200 Proof ethanol, molecular biology or ACS grade	Any major lab supplier (MLS)
Nuclease-free water, molecular biology grade	Any MLS
8-channel pipettes	Any MLS
Wide orifice tips (200µL)	Rainin 30389241
0.2 mL 8-tube strips	USA Scientific TempAssure 1402-4708
Microcentrifuge	Any MLS
Magnetic separation rack compatible with 0.2 mL 8-tube strips	V&P Scientific VP 772F4-1
Magnetic bead rack	ThermoFisher Scientific 12321D
Thermocycler compatible with 0.2 mL 8-tube strips	Any MLS
Tube rotator	Any MLS
1.5 mL DNA LoBind tubes	Eppendorf 022431021
0.5 mL DNA LoBind tubes	Eppendorf 022431005

*Sold as part of bundle Kinnex[™] single-cell RNA kit (103-072-200)



Safety precautions

Refer to the Safety Data Sheet (SDS) for information on reagent hazards and protocols for safe handling, use, storage, and disposal.

cDNA Input

- Optimal range of 3,000–10,000 target cell recovery from the 10x Chromium 3' or 5' single cell workflow.
- Follow the best practices in the 10x Chromium user guide up until cDNA amplification, cleanup and QC (refer to 10x Chromium user guide for <u>5' kits</u> and <u>3' kits</u>). Input cDNA quality control is highly recommended before proceeding to the Kinnex workflow.

Reagent and sample handling

- Room temperature is defined as any temperature in the range of 18–23°C for this protocol.
- Thaw the repair buffer, nuclease buffer, Kinnex ligase buffer, Kinnex adapter mix, and primers at room temperature.
- Briefly vortex reagent buffers and Kinnex adapters prior to use. Enzyme mixes do not require vortexing.
- Quick spin all reagents in a microcentrifuge to collect liquid at tube bottom prior to use.
- Keep all temperature-sensitive reagents on ice.

Temperature-sensitive reagents			
Step used	Tube color	Reagent	
	Green	Kinnex single-cell PCR mix 103-244-500	
	Yellow	Kinnex 3' capture primer mix 103-182-400	
	Red	Kinnex 5' capture primer mix 103-182-200	
		Kinnex primers premix (A-PQ)	
		103-107-800 A	
	Orange	103-107-900 B	
		103-108-000 C	
TOO DOD and Kinney		103-108-100 D	
TSO PCR and Kinnex PCR		103-108-200 E	
		103-108-300 F	
		103-108-400 G	
		103-153-000 H	
		103-153-100 l	
		103-153-200 J	
		103-153-300 K	
		103-153-400 L	
		103-153-500 M	



		103-153-600 N
		103-153-700 0
		103-153-800 PQ
	Light green	Kinnex single-cell enzyme 103-243-800
	Yellow	Kinnex single-cell ligase 103-244-000
	White	Kinnex single-cell ligase buffer 103-244-100
	Red	Kinnex single-cell ligation additive 103-244- 400
Kinnex array formation		Kinnex adapter mix
		bc01 103-109-600
	Blue	bc02 103-109-700
		bc03 103-109-800
		bc04 103-109-900
	Green	DNA repair mix 103-110-000
DNA damage repair	Purple	Repair buffer 102-244-300
Nuclease treatment	Light green	Nuclease mix 103-110-100
	Light purple	Nuclease buffer 103-110-200

- Bring SMRTbell cleanup beads and Qubit 1X dsDNA HS reagents to room temperature for 30 minutes prior to use.
- Pipette-mix all bead binding and elution steps until beads are distributed evenly in solution.
- Wide-bore pipette tips help to minimize foaming specifically when resuspending Kinnex capture beads.
- Pipette-mix all library prep reactions by pipetting up and down 10 times or until fully resuspended.
- Samples can be stored at 4°C at all safe stopping points listed in the protocol.
- 1.5X SMRTbell cleanup is recommended before Kinnex array formation. If the cDNA contains smaller fragments

<200 bp, it is recommended to increase the SMRTbell cleanup ratio to 2X.

The Kinnex capture beads kit is not temperature-sensitive and contains the following contents:

Tube color	Reagent
Clear	Kinnex capture beads 102-144-900 Kinnex bead binding buffer 103-145-40 Kinnex bead washing buffer 103-145-500



Thermocycler programs

Program thermocycler(s) prior to beginning the protocol for the first time.

1. TSO PCR program (15 ng input)

Heated lid set at 105°C

Step	Time	Temperature	Cycles
1	3 min	98°C	1
2	20 sec	98°C	
3	30 sec	65°C	5
4	4 min	72°C	
5	5 min	72°C	1
6	Hold	4°C	1

or TSO PCR program (60-75 ng input)

Heated lid set at 105°C

Step	Time	Temperature	Cycles
1	3 min	98°C	1
2	20 sec	98°C	
3	30 sec	65°C	3
4	4 min	72°C	
5	5 min	72°C	1
6	Hold	4°C	1

2. TSO artifact removal program

Heated lid set at 47°C

Step	Time	Temperature
1	30 min	37°C
2	Hold	4°C



3. Kinnex PCR program

Heated lid set at 105°C

Step	Time	Temperature	Cycles
1	3 min	98°C	1
2	20 sec	98°C	
3	30 sec	68°C	9
4	4 min	72°C	
5	5 min	72°C	1
6	Hold	4°C	1

4. Kinnex primer digestion program

Heated lid set at 47°C

Step	Time	Temperature
1	30 min	37°C
2	Hold	4°C

5. Kinnex array ligation program

Heated lid set at 52°C

Step	Time	Temperature
1	60 min	42°C
2	Hold	4°C

6. DNA damage repair program

Heated lid set at 47°C

Step	Time	Temperature
1	30 min	37°C
2	Hold	4°C

7. Nuclease treatment program

Heated lid set at 47°C

Step	Time	Temperature
1	60 min	37°C
2	Hold	4°C



Workflow steps

1. Input cDNA quality control

This protocol requires at least 15 ng of 10x Chromium 3' single cell cDNA. Follow the 10x Chromium user guide up until cDNA amplification, cleanup and QC (refer to 10x Chromium user guide for <u>5' kits</u> and <u>3' kits</u>). Before you begin, evaluate the quantity and size distribution of input cDNA to determine whether it is suitable for the protocol (average size between 500–1500 bp).

~	Step	Instructions
	1.1	Bring the Qubit 1X dsDNA HS working solution and standards to room temperature.
	1.2	Pulse vortex or pipette mix each sample to homogenize the DNA in solution.
	1.3	Quick spin each sample to collect liquid.
	1.4	Take a 1 µL aliquot from each sample.
	1.5	Measure DNA concentration with a Qubit fluorometer using the 1X dsDNA HS kit.
	1.6	Dilute each sample to 1.0-1.5 ng/ μ L in elution buffer or water, based on the Qubit reading.
	1.7	Measure DNA size distribution with a Bioanalyzer system using the High Sensitivity DNA Kit.
	1.8	Proceed to the next step of the protocol if sample quality is acceptable.

2. TSO PCR

This PCR step generates biotinylated DNA-fragments to enable removal of TSO priming artifacts generated during cDNA synthesis.

cDNA amplification with Kinnex capture primers:

Red

✓	Step	Instructions				
		Normalize cDNA sample input to 15 ng if it is between 15 ng and 59 ng using elution buffer. Normalize cDNA sample input to 75 ng if it is higher than 75 ng using elution buffer. For cDNA amounts between 60–75ng, proceed without normalizing. Select either the Kinnex 3' or capture primer mix depending on the 10x Genomics kit used. Set up the following PCR reaction or ice (RM1).				
					•	
2.1		Reacti	on Mix 1 (RM	1):		
	2.1	~	Tube color	Component	Volume	
				Nuclease-free water	Up to 50 µL	
Croop Kinnex single-cell PCR mix 25 ul						

25 µL

10 µL

(103-244-500)

(103-182-200)

Kinnex 5' capture primer mix

10



Yellow	Kinnex 3' capture primer mix (103-182-400)	
	10x 5' or 3' cDNA library (1−5 ng/µL)	Up to 15 µL
	Total volume	50 µL

- 2.2 Pipette-mix **RM1**.
- 2.3 Quick spin **RM1** in a microcentrifuge to collect liquid.
- 2.4 Select the <u>TSO PCR program</u> based on cDNA input. Keep sample on ice until thermal cycler lid has heated to 105°C.

TSO PCR program (15 ng input)

Heated lid set at 105°C

Step	Time	Temperature	Cycles
1	3 min	98°C	1
2	20 sec	98°C	
3	30 sec	65°C	5
4	4 min	72°C	
5	5 min	72°C	1
6	Hold	4°C	1

or TSO PCR program (60-75 ng input)

Heated lid set at 105°C

Step	Time	Temperature	Cycles
1	3 min	98°C	1
2	20 sec	98°C	
3	30 sec	65°C	3
4	4 min	72°C	
5	5 min	72°C	1
б	Hold	4°C	1

Cleanup with 1.5X SMRTbell cleanup beads

2.5	Add 1.5X v/v (75 $\mu L)$ of resuspended, room-temperature SMRTbell cleanup beads to each tube of amplified cDNA.
2.6	Pipette mix the beads until evenly distributed.
2.7	Quick spin the tube strip in a microcentrifuge to collect liquid.
2.8	Leave at room temperature for 10 minutes to allow DNA to bind beads.

- 2.9 Place tube strip in a magnetic separation rack until beads separate fully from the solution.
- 2.10 Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.



2.11	Slowly dispense 200 μL , or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, pipette off the 80% ethanol and discard.
2.12	Repeat the previous step.
2.13	 Remove residual 80% ethanol: Remove tube strip from the magnetic separation rack. Quick spin tube strip in a microcentrifuge. Place tube strip back in a magnetic separation rack until beads separate fully from the solution. Pipette off residual 80% ethanol and discard.
2.14	Remove tube strip from the magnetic rack. Immediately add 42 µL of elution buffer to each tube and resuspend the beads by pipetting 10 times or until evenly distributed.
2.15	Quick spin the tube strip in a microcentrifuge to collect liquid.
2.16	Leave at room temperature for 5 minutes to elute DNA.
2.17	Place tube strip in a magnetic separation rack until beads separate fully from the solution.
2.18	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new tube strip. Discard old tube strip with beads.
2.19	Recommended: Evaluate sample concentration. Take a 1 µL aliquot from each tube. Measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit. Yield typically ranges from 100–600 ng.
2.20	Proceed to the next step of the protocol if sample quantity is acceptable (between 100 ng – 600 ng).

SAFE STOPPING POINT – Store at 4°C



3. TSO artifact removal

In this step, removal of DNA fragments containing TSO artifacts is performed using Kinnex capture beads. A <u>tutorial</u> <u>video</u> demonstrating this step is available.

✓	Step	Instructions
	3.1	Bring Kinnex capture beads kit to room temperature. Resuspend the beads by vortexing.
	3.2	Transfer 10 μ L resuspended Kinnex capture beads per sample to a PCR tube. Scale up the amount of beads if processing more than 4 samples (with 10% overage). If preparing more than 40 μ L of beads, use a 1.5 mL LoBind tube instead of PCR tube.
	3.3	Place the tube on the magnet until the beads separate fully from the solution.
	3.4	Carefully remove and discard the supernatant while the tube remains on the magnet. Avoid touching the bead pellet with the pipette tip.
	3.5	 Remove the tube from the magnet. Add 40 µL Kinnex bead binding buffer along the inside wall of the tube where the beads are collected and gently resuspend by pipetting using wide bore tips. DO NOT VORTEX. Note: the solution may be viscous. Highly recommend using wide bore tips to avoid foaming. When excess bubbles are present, lower cDNA recovery is expected. Quick-spin the tube in a microcentrifuge if needed. Note: Scale up the volume of Kinnex capture binding buffer accordingly if preparing more than 40 µL of beads.
	3.6	Place the tube on the magnet until the beads separate fully from the solution and remove the supernatant.
	3.7	 Resuspend the beads in 40 µL Kinnex bead binding buffer by pipetting slowly using wide bore tips. DO NOT VORTEX. Note: the solution may be viscous. Highly recommend using wide bore tips to avoid foaming. When excess bubbles are present, lower cDNA recovery is expected. Note: Scale up the volume of Kinnex capture binding buffer accordingly, if preparing more than 40 µL of beads. Distribute 40 µL of resuspended Kinnex capture beads into the appropriate number of PCR tubes before proceeding to Step 3.8.
	3.8	Add 40 µL of a solution containing the biotinylated DNA-fragments (from <u>Step 2.18</u>) to the resuspended beads. Mix carefully using wide bore tips to avoid foaming of the solution.
	3.9	Incubate the tube at room temperature for 15 minutes on a rotator to keep the beads in suspension. Quick-spin the tube in a microcentrifuge to collect liquid.
	3.10	Place the tube on the magnet until the beads separate fully from the solution and remove the supernatant.
	3.11	Resuspend the Kinnex capture beads/DNA-complex in 80 µL Kinnex bead washing buffer by rigorous pipette mixing until evenly distributed.
	3.12	Place the tube on the magnet until the beads separate fully from the solution and remove the



supernatant.

0.10	Remove the tube from the magnet. Resuspend the Kinnex capture beads/DNA-complex in 80 μ L
3.13	Kinnex bead washing buffer by rigorous pipette mixing until evenly distributed.

- 3.14 Place the tube on the magnet until the beads separate fully from the solution and remove the supernatant.
- 3.15 Remove the tube from the magnet. Resuspend the Kinnex capture beads/DNA complex in 80 μL nuclease free water by pipette mixing until evenly distributed.
- 3.16 Place the tube on the magnet until the beads separate fully from the solution and remove the supernatant.
- 3.17 Resuspend the capture beads/DNA-complex in 40 µL of elution buffer by pipette mixing until evenly distributed.
- 3.18 Add 2 µL Kinnex enzyme to the sample with capture beads to cleave the captured DNA products from Kinnex capture beads.
- 3.19 Pipette-mix each sample and a very quick spin in a microcentrifuge to collect liquid.

Run the TSO artifact removal program.

TSO artifact removal program

Heated lid set at 47°C

3.20

Step	Time	Temperature
1	30 min	37°C
2	Hold	4°C

3.21 Place the tube on the magnet for 1 minute and move the supernatant containing the library to a fresh tube.

- 3.22 Cleanup with 1.5X SMRTbell cleanup beads.
- 3.23 Add 1.5X v/v (63 μL) of resuspended, room-temperature SMRTbell cleanup beads to each sample.
- 3.24 Pipette-mix the beads until evenly distributed.
- 3.25 Quick-spin the tube strip in a microcentrifuge to collect liquid.
- 3.26 Leave at room temperature for 10 minutes to allow DNA to bind beads.
- 3.27 Place tube strip in a magnetic separation rack until the beads separate fully from the solution.
- 3.28 Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
 - 3.29 Slowly dispense 200 µL, or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, pipette off the 80% ethanol and discard.
 - 3.30 Repeat the previous step.

Remove residual 80% ethanol:

- 3.31 Remove the tube strip from the magnetic separation rack.
 - Quick-spin the tube strip in a microcentrifuge.



	• Place the tube strip back in a magnetic separation rack until the beads separate fully from the solution.
	Pipette off residual 80% ethanol and discard.
3.32	Remove the tube strip from the magnetic rack. Immediately add 46 μ L of elution buffer to each tube and resuspend the beads by pipetting 10 times or until evenly distributed.
3.33	Quick-spin the tube strip in a microcentrifuge to collect liquid.
3.34	Leave at room temperature for 5 minutes to elute DNA.
3.35	Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
3.36	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new 0.5 mL LoBind tube. Discard the old tube strip with beads.
	Recommended: Evaluate sample concentration.
3.37	• Take a 1 µL aliquot from each tube.
	Measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit.
3.38	Proceed to the next step of the protocol if sample quantity is acceptable (minimum 25 ng).
	SAFE STOPPING POINT – Store at 4°C

4. Kinnex PCR

Perform 16 parallel cDNA amplification reactions with Kinnex primers to generate DNA fragments containing orientationspecific Kinnex segmentation adapter sequences.

Note: All 16 PCR reactions are required. Missing/incorrect primer pairs will result in no/low SMRTbell yield.

✓	Step	Instructions						
		Set up the following PCR reaction Reaction Mix 2 (RM2):	mix per sample on ice.					
		 Master mix components 	Volume for 16X concatenation*					
		PCR-grade water	176-X µL					
	4.1	Kinnex single-cell PCR mix (103-244-500)	220 µL					
		25 ng of amplified cDNA from Step 3.38	ΧμL					
		Total volume	396 µL					
		X = 25 (ng) / purified pooled cDNA concentra *10% overage included	tions from Step 3.38					
	4.2	Quick-spin RM2 in a microcentrifuge to collect liquid.						
	4.3	Add 22.5 µL of RM2 to a new PCR	tube on ice. Repeat this					



sample (each containing 22.5 µL of RM2).

Add 2.5 µL of Kinnex primers premix into each of 16 PCR tubes on ice according to the table below.

PCR tube	Kinnex primers premix	PN
1	Kinnex primers premix A	103-107-800
2	Kinnex primers premix B	103-107-900
3	Kinnex primers premix C	103-108-000
4	Kinnex primers premix D	103-108-100
5	Kinnex primers premix E	103-108-200
6	Kinnex primers premix F	103-108-300
7	Kinnex primers premix G	103-108-400
8	Kinnex primers premix H	103-153-000
9	Kinnex primers premix I	103-153-100
10	Kinnex primers premix J	103-153-200
11	Kinnex primers premix K	103-153-300
12	Kinnex primers premix L	103-153-400
13	Kinnex primers premix M	103-153-500
14	Kinnex primers premix N	103-153-600
15	Kinnex primers premix O	103-153-700
16	Kinnex primers premix PQ	103-153-800

4.4

4.5 Pipette-mix each sample. The total volume of each tube should be 25.0 μL.

4.6 Quick-spin the strip tubes in a microcentrifuge to collect liquid.

Run the <u>Kinnex PCR program</u> with heated lid set to 105°C. Keep sample on ice until thermal cycler lid has heated to 105°C.

	Step	Time	Temperature	Cycles
	1	3 min	98°C	1
4.7	2	20 sec	98°C	
	3	30 sec	68°C	9
	4	4 min	72°C	
	5	5 min	72°C	1
	6	Hold	4°C	1

Clean up with 1.5X SMRTbell cleanup beads.

- 4.8 Pool entire volume of all 16 reactions into a single 1.5 mL LoBind tube.
- 4.9 Add 1.5X v/v (600 µL) of resuspended, room-temperature SMRTbell cleanup beads to the PCR pool.
- 4.10 Pipette-mix the beads until evenly distributed.



- 4.11 Quick-spin the tube strip in a microcentrifuge to collect all liquid from the sides of the tubes.
- 4.12 Leave at room temperature for 10 minutes to allow DNA to bind beads.
- 4.13 Place 1.5mL LoBind tube in a magnetic separation rack until beads separate fully from the solution.
- 4.14 Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
- 4.15 Slowly dispense 1 mL, or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, pipette off the 80% ethanol and discard.
- 4.16 Repeat the previous step.

4.17

Remove residual 80% ethanol:

- Remove the LoBind tube from the magnetic separation rack.
- Quick-spin the LoBind tube in a microcentrifuge.
- Place the LoBind tube back in a magnetic separation rack until beads separate fully from the solution.

Pipette off residual 80% ethanol and discard.

- $\begin{array}{c} \mbox{4.18} \\ \mbox{and resuspend the beads.} \end{array} \label{eq:4.18} Remove the LoBind tube from the magnetic rack. Immediately add 47 <math display="inline">\mu L$ of elution buffer to each tube and resuspend the beads. \end{array}
- 4.19 Quick-spin the LoBind tube in a microcentrifuge.
- 4.20 Incubate at room temperature for 5 minutes to elute DNA.
- 4.22 Place the LoBind tube in a magnetic separation rack until beads separate fully from the solution.
- 4.22 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new LoBind tube. Discard old tube with beads.

Recommended: Evaluate sample concentration.

- Take a 1 μ L aliquot from each tube, dilute with 9 μ L of elution buffer.
- 4.23 Using 1 µL of the dilution, measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit.
- 4.24 Proceed to the next step of the protocol if sample quantity is acceptable (maximum input: 5 µg). Do not proceed if less than 3 µg is available.

SAFE STOPPING POINT - Store at 4°C



5. Kinnex array formation

In this step, treat PCR amplified cDNA fragments with Kinnex enzyme to create single-stranded extensions to enable directional assembly of cDNA segments into a linear array.

	directional assembly of cDNA segments into a linear array.						
✓	Step	Instructions					
	5.1	In a 0.2 mL PCR tube, add 5 µg of sample from Step 4.22, in 47 µL of volume (106 ng/µl). Dilute with elution buffer going into this step if sample is too concentrated.					
	5.2	Add 10 µL of Kinnex enzyme to create single-stranded extensions on PCR-amplified cDNA fragments to enable subsequent directional assembly of 16 PCR products.					
	5.3	Pipette-mix each sample.					
	Run the Kinnex primer digestion program.						
		Kinnex p	rimer digesti	on program			
	5.4	Heated lid	l set at 47°C				
	0.1	Step	Time	Temperature			
		1	30 min	37°C			
		2	Hold	4°C			
	5.5			oter barcode 01–04 (sample for a total vo Component	· •	-	
			Blue	Kinnex barcode adap			
			Red	Kinnex ligation additiv		_	
				Total volume	23 µL	-	
	5.6	Pipette-m	ix each sampl	e.			
	Add the following components in the order and volume listed below to a new microcentrifuge Adjust component volumes for the number of samples being prepared, plus 10% overage. Fo individual preps, add components directly to each sample in the order and volume listed below Reaction Mix 3 (RM3):						
	5.7	~	Tube color	Component		Volume	
			White	Kinnex single-cell liga 103-244-100	se buffer	10 µL	
			Yellow	Kinnex single-cell ligas 103-244-000	se	10 µL	
				Total volume		20 µL	
	5.8	Pipette-m	ix RM3				
	5.9	Quick-spin RM3 in a microcentrifuge to collect liquid.					
	= 1.0	Add 20 yell of DN/2 to each communic					

- 5.10 Add 20 µL of **RM3** to each sample.
- 5.11 Pipette-mix each sample .



Run the Kinnex array ligation program.

Kinnex array ligation program

Heated lid set at 52°C

510

5.21

5.27

J. I Z	Step	Time	Temperature
	1	60 min	42°C
	2	Hold	4°C

Cleanup with 1.2X SMRTbell cleanup beads

- 5.13 Add 1.2X v/v (120 µL) of resuspended, room-temperature SMRTbell cleanup beads to each sample.
- 5.14 Pipette-mix the beads until evenly distributed.
- 5.15 Quick-spin the tube strip in a microcentrifuge to collect liquid.
- 5.16 Leave at room temperature for 10 minutes to allow DNA to bind beads.
- 5.17 Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
- 5.18 Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
- 5.19 Slowly dispense 200 μL, or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, pipette off the 80% ethanol and discard.
- 5.20 Repeat the previous step.

Remove residual 80% ethanol:

- Remove the tube strip from the magnetic separation rack.
- Quick-spin the tube strip in a microcentrifuge.
 - Place the tube strip back in a magnetic separation rack until beads separate fully from the solution.

Pipette off residual 80% ethanol and discard.

- 5.22 Remove the tube strip from the magnetic rack. Immediately add 43 µL of elution buffer to each tube and resuspend the beads by pipetting 10 times or until evenly distributed.
- 5.23 Quick-spin the tube strip in a microcentrifuge to collect liquid.
- 5.24 Leave at room temperature for 5 minutes to elute DNA.
- 5.25 Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
- 5.26 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new 0.5 mL LoBind tube. Discard old tube strip with beads.

Recommended: Evaluate sample concentration.

- Take a 1 μ L aliquot from each tube, dilute with 4 μ L of elution buffer.
- Measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit.

Proceed to DNA damage repair.

SAFE STOPPING POINT - Store at 4°C



6. DNA damage repair

✓	Step	Instruct	ions				
	6.1	ln a 0.2 r	mL PCR tube	e, add 42 µL of sample	e from Step	5.26	
		Add the following components in the order and volume listed below to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus 10% overage. For individual preps, add components directly to each sample in the order and volume listed below. Reaction Mix 4 (RM4):					
	6.2	V	Tube co	lor Component	Volume		
			Purple	Repair buffer 102-696-100	6 µL		
			Green	DNA repair mix 102-696-000	2 µL		
				Total volume	8 µL		
	6.3	Pipette-r	mix RM4 .				
	6.4	Quick-sp	pin RM4 in a	microcentrifuge to c	ollect liqu	id.	
	6.5	Add 8 µl	_ of RM4 to	each sample. Total v	olume sho	puld equal 50 μL.	
	6.6	Pipette-mix each sample. Quick-spin the strip tube in a microcentrifuge to collect liquid.					
	6.7						
	6.8 Run the <u>DNA damage repair program</u> . DNA damage repair program						
		Heated lid set at 47°C					
		Step	Time	Temperature			
		1	30 min	37°C			
		2	Hold	4°C			
				Clean up with 1.2X S	SMRTbell	cleanup beads	
	6.9	Add 1.22	X ν/ν (60 μL)	of resuspended, roo	m-temper	rature SMRTbell cleanup beads to each sample.	
	6.10			ds until evenly distrib			
	6.11			strip in a microcentri	0		
	6.12	Leave at	t room temp	erature for 10 minute	es to allow	DNA to bind beads.	
	6.13	Place th	e tube strip	in a magnetic separa	tion rack ı	until beads separate fully from the solution.	
	6.14	Slowly p	pipette off the	e cleared supernatan	t without	disturbing the beads. Discard the supernatant.	
	6.15) µL, or enough to cov ds, pipette off the 80		ads, of freshly prepared 80% ethanol into each and discard.	
	6.16	Repeat t	the previous	step.			



Remove residual 80% ethanol:

6.17

- Remove the tube strip from the magnetic separation rack.
- Quick-spin the tube strip in a microcentrifuge.
- Place the tube strip back in a magnetic separation rack until beads separate fully from the solution.
 - Pipette off residual 80% ethanol and discard.
- 6.18 Remove the tube strip from the magnetic rack. Immediately add 40 μL of elution buffer to each tube and resuspend the beads by pipetting 10 times or until evenly distributed.
- 6.19 Quick-spin the tube strip in a microcentrifuge to collect liquid.
- 6.20 Leave at room temperature for 5 minutes to elute DNA.
- 6.21 Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
- 6.22 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new tube strip. Discard old tube strip with beads.

SAFE STOPPING POINT - Store at 4°C

7. Nuclease treatment and cleanup

\checkmark	Step	Instruction	ons				
	71	Adjust co individua volume li	omponent volur	onents in the orde mes for the numbe mponents directly	er of samp		
	7.1	v	Tube color	Component	Volume		
			Light purple	Nuclease buffer 103-110-200	5 µL		
			Light green	Nuclease mix 103-110-100	5 µL		
				Total volume	10 µL		
	7.2	Pipette-n	nix RM5 .				
	7.3	Quick-sp	in RM5 in a mic	crocentrifuge to co	ollect liquid		
7.4 Add 10 μ L of RM5 to each sample. Total volume should equal 50 μ L.							
	7.5	Pipette-mix each sample.					
	7.6	Quick-sp	in the strip tube	e in a microcentrifu	uge to coll		



	Run the <u>nuclease treatment program</u> . Nuclease treatment program Heated lid set at 47°C						
7.7	Step ⁻	Time	Temperature				
	1	60 min	37°C				
	2	Hold	4°C				
7.8	Add 60 µL until evenly		cleanup beads to each sample from the previous step. Pipette-mix the beads ed.				
7.9	Quick-spin	the tube s	strip in a microcentrifuge to collect all liquid.				
7.10	Leave at ro	om tempe	erature for 10 minutes to allow DNA to bind beads.				
7.11	Place the t	ube strip ir	n a magnetic separation rack until beads separate fully from the solution.				
7.12	2 1 1		e cleared supernatant without disturbing the beads. It is recommended to save nother tube strip in case of poor DNA recovery.				
7.13	, ,		μL, or enough to cover the beads, of freshly prepared 80% ethanol into each ds, pipette off the 80% ethanol and discard.				
7.14	Repeat the	previous	step.				
	Remove re	sidual 80%	% ethanol:				
	Remove tube strip from the magnetic separation rack.						
7.15	• Quick	spin tube	e strip in a microcentrifuge.				
			back in a magnetic separation rack until beads separate fully from the solution.				
			0% ethanol and discard.				
7.16	and resusp	end the be	p from the magnetic rack. Immediately add 26 μL of elution buffer to each tube eads by pipetting 10 times or until evenly distributed.				
7.17	Quick-spin	the tube s	trip in a microcentrifuge to collect liquid.				
7.18	Leave at ro	om tempe	erature for 5 minutes to elute DNA.				
7.19	Place the t	ube strip ir	n a magnetic separation rack until beads separate fully from the solution.				
7.20	211		e cleared supernatant without disturbing the beads. Transfer supernatant to a ube. Discard old tube strip with beads.				
			rom each tube. Measure DNA concentration with a Qubit fluorometer using the 1: late the total mass.				
7.21	Recomme	nded: Furt	her dilute each aliquot to 250 pg/µL with Femto Pulse dilution buffer. Measure				
	final SMRT	bell library	y size distribution with a Femto Pulse system.				
7.22	concentra	tion must	nk Sample Setup to prepare the SMRTbell library for sequencing. DNA be less than 20 ng/µl to go into ABC. Using a concentration above 20 ng/µl bading during sequencing.				
7.23			ries at 4°C if sequencing within the week. Long-term storage should be at -20°C. v cycles when handling SMRTbell libraries.				
			PROTOCOL COMPLETE				

PacBi

Revision history (description)	Version	Date
Initial release	01	October 2023
Updated to incorporate PCR program into workflow steps, clarify sample quantity in TSO artifact removal, and add required DNA concentration before sequencing	02	March 2024

Research use only. Not for use in diagnostic procedures. © 2024 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at <u>pacb.com/license</u>. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, and PureTarget are trademarks of PacBio.

© 2024 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures. 103-254-300 REV 02 MAR2024

