# Preparing Kinnex<sup>™</sup> libraries using the Kinnex full-length RNA kit



Procedure & checklist


# Before you begin

This procedure describes the workflow for constructing Kinnex full-length RNA libraries from total RNA samples for sequencing on PacBio<sup>®</sup> Sequel<sup>®</sup> II, Sequel IIe, and Revio<sup>™</sup> systems.

Overview				
Samples	1-24			
Workflow time	1.5 days (for up to 24 samples)			
Number of SMRT <sup>®</sup> Cells per	>2 SMRT Cells for Revio system			
Kinnex library Prep	>4 SMRT Cells for Sequel II/IIe systems			
RNA input				
RNA input Quality/size distribution	RIN (RNA integrity number) ≥7.0			



# Workflow





# Required materials and equipment

RNA and DNA sizing	
2100 Bioanalyzer instrument	Agilent Technologies G2939BA
RNA 6000 Nano kit	Agilent Technologies 5067-1511
Agilent High Sensitivity DNA Kit	Agilent Technologies 5067-4626
Femto Pulse system	Agilent Technologies M5330AA
Genomic DNA 165 kb kit	Agilent Technologies FP-1002-0275
DNA quantitation	
Qubit Fluorometer	Thermo Fisher Scientific Q33238
Qubit 1X dsDNA HS Assay kit	Thermo Fisher Scientific Q33230
Qubit RNA HS Assay kit	Thermo Fisher Scientific Q32852
cDNA synthesis and amplification	
Iso-Seq <sup>®</sup> Express 2.0 Kit	PacBio <sup>®</sup> 103-071-500*
SMRTbell <sup>®</sup> cleanup beads	PacBio <sup>®</sup> 102-158-300*
Elution buffer (50 mL)	PacBio <sup>®</sup> 101-633-500*
Kinnex™ Library Prep	
Kinnex™ PCR 8-fold kit	PacBio <sup>®</sup> 103-071-600*
Kinnex™ concatenation kit	PacBio <sup>®</sup> 103-071-800*
	*Part of the Kinnex™ full-length RNA kit bundle (103-072-000)

#### Other Supplies

200 Proof ethanol, molecular biology or ACS grade	Any major lab supplier (MLS)
Nuclease-free water, molecular biology grade	Any MLS (e.g., Sigma-Aldrich W4502)
8-channel pipettes – P20 & P200	Any MLS
Single-channel pipette – P2, P10, P20, P100 or P200	Any MLS
0.2 mL 8-tube strips	USA Scientific TempAssure 1402-4708
Microcentrifuge	Any MLS
Magnetic separation rack compatible with 0.2 mL 8-tube strips	Any MLS (e.g., V&P Scientific, Inc. VP772F4-1)
Thermal cycler compatible with 0.2 mL 8-tube strips	Any MLS
DNA LoBind tubes	Eppendorf 022431021



# General best practices

Take care to accurately pipette SMRTbell cleanup beads because small changes in volume can significantly alter the size distribution of your sample.

Equilibrate the SMRTbell cleanup beads at room temperature for 30 mins prior to use.

In cDNA amplification and Kinnex PCR, keep sample(s) on ice until thermal cycler lid has reached 105°C to avoid digestion of primers by polymerase exonuclease activity.

This workflow takes 1.5 days to complete. If a stop is necessary, refer to the workflow for safe stopping points.

# Safety precautions

Refer to the Safety Data Sheet (SDS) for information on reagent hazards and protocols for safe handling, use, storage, and disposal.

# Multiplexing best practices

Multiplexing can be achieved with one of the three following methods:

- Barcoded cDNA primers using Iso-Seq primers bc01-12 in <u>step 3</u> of the protocol. To multiplex, use the Iso-Seq cDNA amplification primer in combination with Iso-Seq primers bc01-12 to amplify samples. After SMRTbell cleanup, Iso-Seq samples can be pooled and brought through a single Kinnex PCR reaction. Each barcoded primer is sufficient for 2 reactions, with the Iso-Seq kit supporting a total of 24 reactions.
- 2. Barcoded adapters using Kinnex adapters bc01–04. In this case, use barcoded adapters at <u>step 5</u> "Kinnex array formation" in the workflow.
- 3. A combination of the above 2 approaches to achieve 48-plex.

Note: if not performing multiplexing, the same Iso-Seq primer barcodes and Kinnex adapter barcodes are still used, but without pooling.



# **Reagents list**

Iso-Seq express 2.0 kit 103-071-500				
	Tube color	Reagent		
	Purple	Iso-Seq RT buffer 103-103-900		
	Orange	Iso-Seq RT primer mix 103-104-000		
	Yellow	Iso-Seq RT enzyme mix 103-104-100		
Kinnex array formation	Red	lso-Seq cDNA PCR mix 103-104-200		
	Blue	Iso-Seq Express TSO2.0 103-104-300		
	Green	Iso-Seq cDNA amplification primer 103-104-400		
		lso-Seq primer bc01 103-104-500		
		lso-Seq primer bc02 103-104-600		
		lso-Seq primer bc03 103-104-700		
	White	lso-Seq primer bc04 103-104-800		
		lso-Seq primer bc05 103-104-900		
		lso-Seq primer bc06 103-105-000		
		lso-Seq primer bc07 103-105-100		
		lso-Seq primer bc08 103-105-200		
		lso-Seq primer bc09 103-105-300		
		lso-Seq primer bc10 103-105-400		
		Iso-Seq primer bc11 103-105-500		
		lso-Seq primer bc12 103-105-600		



#### Kinnex PCR 8-fold kit 103-071-600

#### Tube color Reagent

Green	Kinnex PCR mix 103-107-700		
	Kinnex primer mix A 103-107-800		
	Kinnex primer mix B 103-107-900		
	Kinnex primer mix C 103-108-000		
Orongo	Kinnex primer mix D 103-108-100		
Orange	Kinnex primer mix E 103-108-200		
	Kinnex primer mix F 103-108-300		
	Kinnex primer mix G 103-108-400		
	Kinnex primer mix HQ 103-108-500		

#### Kinnex concatenation kit 103-071-800

Tube color	Reagent	
Red	Kinnex enzyme 103-110-400	
Yellow	Kinnex ligase 103-110-500	
White Kinnex array and repair buffer 103-110-300		
Green DNA repair mix 103-110-000		
Light Purple	e Nuclease buffer 103-110-200	
Light Green	Nuclease mix 103-110-100	
	Kinnex adapter bc01 mix 103-109-600	
Dhuo	Kinnex adapter bc02 mix 103-109-700	
Blue	Kinnex adapter bc03 mix 103-109-800	
	Kinnex adapter bc04 mix 103-109-900	



# Workflow steps

# 1. Input RNA quality control

This protocol requires high-quality RNA. Prior to library preparation, evaluate the size distribution of the input RNA to determine whether it is suitable for the protocol.

✓	Step	Instructions			
	1.1	Measure the RNA Integrity Number (RIN) with an Agilent 2100 Bioanalyzer Instrument using the RNA 6000 Nano kit.			
	Proceed to the next step of the protocol if sample quality is acceptable:				
	1 0	RIN Quality recommendations			
	1.2	≥7.0 Recommended. Proceed to next step of the protocol.			
		<7.0 Increased library failure rates or reduced data quality.			
		SAFE STOPPING POINT – Store at -70°C or below			

# 2. cDNA synthesis

### 2.1 cDNA synthesis

In this step, total RNA samples are converted to first-strand cDNA products.

✓	Step	Instructions
	2.1.1	Quick-spin the Iso-Seq RT enzyme mix in the microcentrifuge to collect liquid, then place on ice.
		Thaw the following components at room temperature, briefly vortex to mix, then quick-spin to

Thaw the following components at room temperature, briefly vortex to mix, then quick-spil collect liquid and place on ice.

	Tube color	Reagent
	Orange	Iso-Seq RT primer mix (103-104-000)
	Purple	Iso-Seq RT buffer (103-103-900)
2.1.2	Red	Iso-Seq cDNA PCR mix (103-104-200)
2.1.2	Green	Iso-Seq cDNA amplification primer (103-104-400)
	Blue	Iso-Seq Express TSO2.0 (103-104-300)
	White	Iso-Seq primer barcodes 01–12* (the number of primers thawed will depend on the number of samples processed) 103-104-500 through 103-105-600
	*If processing	only one sample, any of the 12 Iso-Seq barcoded primer



### 2.2 Primer annealing for first-strand synthesis

✓	Step	Instructions				
		For each RNA sample to be processed, prepare reagent mix 1 on ice by adding the following components to each tube in the PCR strip tube.				
		✓ Tube color	Components	Volume		
	2.2.1		Total RNA (300 ng)	<7 µL		
		Orange	Iso-Seq RT primer mix	2 µL		
			Nuclease-free water	Up to 9 µL		
			Total volume	9 µL		
	2.2.2	Thoroughly mix I	by pipetting up and dow	n 10 times.		
	2.2.3	Quick-spin the tube strip in a microcentrifuge to collect liquid.				
		Incubate in a the	rmal cycler with the foll	owing program.	Set the lid temperature to 80°C.	

2.2.4	Temperature	Time
2.2.4	70°C	5 min
	20°C	hold
	Proceed immediatel	y to the next step

### 2.3 Reverse transcription and template switching

$\checkmark$	Step	Instructions				
		For each RNA sample, prepare reagent mix 2 on ice by adding the following components in the order and volume listed below. Adjust component volumes for the number of samples being prepared, plus 10% overage.				
	2.3.1	✓ Tube colo	r Components		Volume	
	2.0.1	Purple	Iso-Seq RT buffer (vortex brief	ly before use)	5μL	
			Nuclease-free Water		3μL	
		Yellow	Iso-Seq RT enzyme mix		2 µL	
			Total volume added per reac	tion	10 µL	
	2.3.2	Pipette-mix and	d quick-spin in a microcentr	ll liquid.		
		Add 10 $\mu$ L of reaction mix 2 to the 9 $\mu$ L from reaction mix 1 ( <u>Section 2.2</u> ) for a total volume of 19 $\mu$ L.				
		✓ Tube	Reagent	Volume		
	2.3.3	Previous	Reagent mix 1 from step 2.2	9µL		
			Reagent mix 2	10 µL		
		Total volur	ne added per reaction			



**2.3.4** Thoroughly mix by pipetting up and down 10 times and then quick-spin to collect all liquid.

Incubate in a thermocycler with the following program. Set the lid temperature to 52°C.

	Temperature	Time
2.3.5	42°C	45 min
	20°C	Hold

Proceed immediately to the next step.

- **2.3.6** Remove the sample tube from the thermal cycler and add 2  $\mu$ L of Iso-Seq template switch oligo to the 19  $\mu$ L reaction at **room temperature** for a total volume of 21  $\mu$ L.
- 2.3.7 Mix by pipetting up and down 10 times and then quick-spin to collect all liquid from the sides of the tube.

Return sample tube to thermal cycler and incubate with the following program. Set the lid temperature to 52°C.

2.3.8	Temperature	Time
	42°C	15 min
	4°C	hold

### 2.4 1.3X SMRTbell bead cleanup

✓	Step	Instructions
	2.4.1	For each sample, add 29 $\mu$ L of elution buffer to the 21 $\mu$ L reverse transcription and template switching reaction (Section 2.3) for a total volume of 50 $\mu$ L.
	2.4.2	Add 65 µL of resuspended, room-temperature SMRTbell cleanup beads.
	2.4.3	Mix beads by pipetting 10 times or until evenly distributed.
	2.4.4.	Quick-spin strip tubes in a microcentrifuge to collect liquid.
	2.4.5	Leave at room temperature for 10 minutes to allow DNA to bind the beads.
	2.4.6	Place the strip tubes in a magnetic separation rack until the beads separate fully from the solution.
	2.4.7	Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
	2.4.8	Slowly dispense 200 $\mu$ L, or enough to cover the beads, of freshly prepared 80% ethanol into the strip tube. After 30 seconds, pipette off the 80% ethanol and discard.
	2.4.9	Repeat the previous step.



- 2.4.10 Remove residual 80% ethanol:
  - Remove the strip tube from the magnetic separation rack.
  - Quick-spin the strip tube in a microcentrifuge.
  - Place the strip tube back in a magnetic separation rack until beads separate fully from the solution.
  - Pipette off residual 80% ethanol and discard.
- 2.4.11 Remove the strip tube from the magnetic rack. Immediately add 21 µL of elution buffer to the strip tube and resuspend the beads by pipetting 10 times or until evenly distributed.
- 2.4.12 Quick-spin the strip tube in a microcentrifuge to collect liquid.
- 2.4.13 Leave at room temperature for 5 minutes to elute the DNA.
- 2.4.14 Place the strip tube in a magnetic separation rack until the beads separate fully from the solution.
- 2.4.15 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer 21 µL of the supernatant to a new strip tube. Discard the old strip tube with beads.
- 2.4.16 Proceed to the next step of the protocol.



# 3. cDNA amplification

First-strand cDNA products are PCR-amplified and barcoded using barcoded Iso-Seq primers at this step.

## 3.1 cDNA amplification

✓	Step	Instructions			
		and volume listed plus 10% overage	d below. Ad e. Pipette m	ljust component volumes for tl	he following components in the one number of samples being prep bc01–12 will be added to each mix.
	3.1.1	✓ Tube color	Componer	nts	Volume
		Red	lso-Seq cD	NA PCR mix	25 µL
		Green	Iso-Seq cDi	NA amplification primer	2 µL
			Total volum	ne	27 µL
	3.1.2	On ice, add 27 $\mu\text{L}$	of reaction	mix 3 to the 21 $\mu$ L of the eluted	cDNA (from previous <u>Section 2.</u>
	3.1.3	Add 2 $\mu$ L of the a	ppropriate	Iso-Seq primer barcode to eac	h sample for a total volume of 50
		✓ Tube color	Componer	nt	Volume
		White	lso-Seq pri	mer	2 µL
	3.1.4	Thoroughly mix b	y pipetting	up and down 10 times and the	en quick spin to collect all liquid.
	3.1.5			ram below with the lid temperanas heated to 105°C.	ature set to 105°C. Keep sample
		PCR program			
		45 seconds at 98°	C 1 cycle		
		10 seconds at 98°	C		
	15 second	15 seconds at 60°	C 10 cycles		
		3 minutes at 72°C	i o cycles		
		5 minutes at 72°C			
	Hold at 4°C				

SAFE STOPPING POINT -- Store at 4°C or -20°C for long-term storage



$\checkmark$	Step	Instructions
	3.2.1	Add 45 $\mu$ L (0.9x) of resuspended, room-temperature SMRTbell cleanup beads to the 50 $\mu$ L of cDNA amplified reaction from <u>Section 3.1</u> . The correct ratio of beads to sample is critical at this step. If targeting longer cDNA, 43 $\mu$ L of SMRTbell cleanup beads (0.86x) can be used.
	3.2.2	Mix beads by pipetting 10 times or until evenly distributed.
	3.2.3	Quick-spin strip tubes in a microcentrifuge to collect liquid.
	3.2.4	Leave at room temperature for 10 minutes to allow DNA to bind beads.
	3.2.5	Place the strip tubes in a magnetic separation rack until beads separate fully from the solution.
	3.2.6	Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
	3.2.7	Slowly dispense 200 $\mu$ L, or enough to cover the beads, of freshly prepared 80% ethanol into the strip tube. After 30 seconds, pipette off the 80% ethanol and discard.
	3.2.8	Repeat the previous step.
		<ul><li>Remove residual 80% ethanol:</li><li>Remove the strip tube from the magnetic separation rack.</li></ul>
	3.2.9	<ul> <li>Quick-spin the strip tube in a microcentrifuge.</li> <li>Place the strip tube back in a magnetic separation rack until beads separate fully from the solution.</li> <li>Pipette off residual 80% ethanol and discard.</li> </ul>
	3.2.10	Remove the strip tube from the magnetic rack. Immediately add 24 $\mu$ L of elution buffer to the strip tube and resuspend the beads by pipetting 10 times or until evenly distributed.
	3.2.11	Quick-spin the strip tube in a microcentrifuge to collect liquid.
	3.2.12	Leave at room temperature for 5 minutes to elute DNA.
	3.2.13	Place the strip tube in a magnetic separation rack until the beads separate fully from the solution.
	3.2.14	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer 24 $\mu$ L of the supernatant to a new strip tube. Discard the old strip tube with beads.
	3.2.15	<ul> <li>Recommended: Measure concentration and size distribution of each cDNA sample.</li> <li>Take a 1 µL aliquot from each strip tube. Dilute each aliquot with 4 µL of elution buffer.</li> <li>Measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit.</li> <li>Dilute 1:4 dilution further to 1.5 ng/µL based on the Qubit reading if needed.</li> <li>Run 1 µL on an Agilent Bioanalyzer using a High Sensitivity DNA kit.</li> </ul>

# 3.2 Cleanup of amplified cDNA using 0.9X SMRTbell Cleanup beads



3.2.16 The expected recovery after cDNA amplification SMRTbell clean-up is >80 ng. A minimum of 55 ng of total cDNA is recommended to proceed with Kinnex PCR (<u>Step 4</u>). If less than 55 ng but more than 25 ng is recovered, proceed with Kinnex PCR but expect lower yields. Do not proceed with less than 25 ng.

#### SAFE STOPPING POINT -- Store at 4°C or -20°C for long-term storage

#### 3.3 Pooling barcoded cDNA (skip if not multiplexing)

$\checkmark$	Step	Instructions
	3.3.1	Using the concentration reading from the Qubit fluorometer, pool an equal mass of each barcoded cDNA sample for a total mass of 55 ng. Store any remaining purified, amplified barcoded cDNA at 4°C for future use.
	3.3.2	Quick-spin the tube strip in a microcentrifuge to collect liquid.
	3.3.3	Proceed to next step of the protocol.

# 4. Kinnex PCR & pooling

#### 4.1 Kinnex PCR

Perform 8 parallel Kinnex PCR reactions with Kinnex primers to generate DNA fragments containing orientationspecific Kinnex segmentation sequences.

✓	Step	Instructions				
	4.1.1 Thaw primers. The entire volume of primers can be t use with a multi-channel pipette.					sferred to an 8-tube strip for ease of
			8X concatenation	Tube color	PN	
		1	Kinnex primer mix A		103-107-800	
		2	Kinnex primer mix B		103-107-900	
		3	Kinnex primer mix C		103-108-000	
		4	Kinnex primer mix D	Orongo	103-108-100	
		5	Kinnex primer mix E	Orange	103-108-200	
		6	Kinnex primer mix F		103-108-300	
		7	Kinnex primer mix G		103-108-400	
		8	Kinnex primer mix HQ		103-108-500	
	4.1.2	Briefly vortex to mix, then quick-spin to collect liquid and place the primer mixes on ice.				
	4.1.3		w the following comp ce on ice. Add the com			ien quick-spin to collect liquid and be.



Master mix components	Volume for 8X concatenation*	
PCR-grade water	88-X µL	X = 55 (ng)/purified pooled DN. concentration from <u>step 3.2.16</u>
Kinnex PCR mix (103-107-700)	110 µL	(single-plex) or <u>step 3.3.3</u> (multiplex)
55 ng of amplified cDNA from <u>Step 3.2.16</u>	XμL	*10% overage included
Total volume	198 µL	

- 4.1.4. Aliquot 22.5 µL of Master Mix 1 into each of the 8 PCR tubes (for 8X concatenation) on ice.
- 4.1.5 Add **2.5 μL** of Kinnex primer mix into each of 8 PCR tubes from step 4.1.4.
- 4.1.6 Thoroughly mix by pipetting up and down 10 times. Quick-spin to collect all liquid.

Set up the thermal cycler program listed below with the lid set to 105°C. Keep sample(s) on ice until the lid is heated to 105°C.

The duration of PCR is approximately 1 hour.

	Step	Temperature	Duration	Cycles
4.1.7	Initial denaturation	98°C	3 min	1
	Denaturation	98°C	20 s	
	Annealing	68°C	30 s	9
	Extension	72°C	4 min	9
	Final extension	72°C	5 min	1
	Hold	4°C	Hold	

### 4.2 Pooling of 8 Kinnex PCR products and 1.05X SMRTbell cleanup

✓	Step	Instructions
	4.2.1	Add exactly 23 $\mu$ L from each of the 8 PCR reactions into a 1.5 mL tube for a total volume of 184 $\mu$ L. An equal volume of each PCR product is necessary for efficient array assembly.
	4.2.2	Add 193 $\mu$ L ( <b>1.05X</b> v/v) of resuspended, room-temperature SMRTbell cleanup beads to a tube of pooled Kinnex PCR amplicon. The correct ratio of beads to pooled sample is critical at this step.
	4.2.3	Pipette-mix the beads until evenly distributed.
	4.2.4	Quick-spin the tube in a microcentrifuge to collect liquid.
	4.2.5	Leave at room temperature for 10 minutes to allow the DNA to bind beads
	4.2.6	Place the tube in a magnetic separation rack until the beads separate fully from the solution.
	4.2.7	Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
	4.2.8	Slowly dispense $200 \ \mu$ L, or enough to cover the beads, of <b>freshly prepared 80% ethanol</b> into the tube. After <b>30 seconds</b> , pipette off the 80% ethanol and discard.



4.2.9 Repeat the previous step.

Remove residual 80% ethanol:

- Remove the tube from the magnetic separation rack.
- Quick-spin the tube in a microcentrifuge.

# • Place the tube back in the magnetic separation rack until the beads separate fully from the solution.

• Pipette off residual 80% ethanol and discard.

4.2.11 Remove the tube from the magnetic rack. **Immediately** add  $40 \mu$ L of **elution buffer** to the tube and resuspend the beads by pipetting 10 times or until evenly distributed.

- 4.2.12 Quick-spin the tube in a microcentrifuge to collect liquid.
- 4.2.13 Leave at room temperature for 5 minutes to elute DNA.
- 4.2.14 Place tube in a magnetic separation rack until beads separate fully from the solution.
- 4.2.15 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a **new tube strip**. Discard old tube with beads.
- 4.2.16 Make a 1:10 dilution of the sample in elution buffer and measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit. Typical yield is  $6-12 \mu g$ .

SAFE STOPPING POINT -- Store at 4°C or -20°C for long-term storage

# 5. Kinnex array formation

### 5.1 Kinnex array formation

4.2.10

In this step, treat PCR-amplified cDNA fragments from <u>Step 4.2</u> with Kinnex enzyme, ligase, and barcoded Kinnex terminal adapters to assemble cDNA segments into a linear array.

✓	Step	Instructior	IS		
	5.1.1	Dilute with Add 2 µL c Note: if co	elution buffer going of Kinnex adapter bc(	into this st 01–04 (sele coded Kinr	from <u>Step 4.2.15</u> , in 39 µL of volume (128 ng/µL). ep if the sample is too concentrated. ect one barcode per library preparation). nex libraries for sequencing, make sure each library ed adapters.
		Tube			
		color	Components	Volume	
		Blue	Kinnex adapter bc	2 µL	
	5.1.2	Add the fo	llowing components	in the liste	d order.



If processing multiple samples, make a master mix with 10% overage. Pipette mix master mix.

Tube color	Components	Volume
White	Kinnex array and repair buffer (103-110-300)	7.0 µL
Red	Kinnex enzyme (103-110-400)	4.0 µL
Yellow	Kinnex ligase (103-110-500)	6.0 µL
	Total RM1 volume	17 µL

Add 17  $\mu$ L of master mix to the PCR tube containing sample and Kinnex adapter. Pipette-mix and run the Kinnex primer digestion/ligation program with the lid set to 55°C.

5.1.3	Step	Temperature	Duration
0.1.0	1	45°C	60 min
	2	4°C	Hold

After running the Kinnex primer digestion/ligation program, add 2  $\mu$ L of DNA repair mix directly to the Kinnex primer digestion/ligation sample.

5.1.4	Tube color	Components	Volume
	Green	DNA repair mix (103-110-000)	2 µL

5.1.5 Thoroughly mix by pipetting up and down 10 times and then quick spin to collect all liquid.

Run the DNA Damage Repair Program with the lid set to >55°C.

5.1.6	Step	Temperature	Duration
	1	45°C	30 min
	2	4°C	Hold

# 5.2 1X SMRTbell bead cleanup

Cleanup with 1X SMRTbell cleanup beads

<b>~</b>	Step	Instructions
	5.2.1	Add 1X v/v (60 $\mu L)$ of resuspended, room temperature SMRTbell cleanup beads to each sample.
	5.2.2	Pipette-mix the beads until evenly distributed and quick-spin in a microcentrifuge to collect liquid.
	5.2.3	Leave at room temperature for 10 minutes to allow the DNA to bind the beads.
	5.2.4	Place the tube strip in a magnetic separation rack until the beads separate fully from the solution.



- 5.2.5 Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
- 5.2.6 Slowly dispense 200 µL, or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, pipette off the 80% ethanol and discard.
- 5.2.7 Repeat the previous step.
- 5.2.8 Remove residual 80% ethanol:
  - Remove the tube strip from the magnetic separation rack.
  - Quick-spin the tube strip in a microcentrifuge.
  - Place the tube strip back in a magnetic separation rack until beads separate fully from the solution.
  - Pipette off residual 80% ethanol and discard.
- 5.2.9 Remove the tube strip from the magnetic rack. Immediately add 40 µL of elution buffer to each tube and resuspend the beads by pipetting 10 times or until evenly distributed. Quick-spin the tube strip in a microcentrifuge to collect liquid.
- 5.2.10 Leave at room temperature for 5 minutes to elute DNA.
- 5.2.11 Place the tube strip in a magnetic separation rack until the beads separate fully from the solution.
- 5.2.12 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a **new PCR strip tube**. Discard old tube with beads.

#### SAFE STOPPING POINT -- Store at 4°C or -20°C for long-term storage

# 6. Nuclease treatment & cleanup

#### 6.1 Nuclease treatment

✓	Step	Instructions			
		volumes for t For individual	he number of sam	ples being p	to a new microcentrifuge tube. Adjust component prepared, plus 10% overage. Pipette mix master mix. tly to each sample from the previous step in the
	611	Nuclease master mix			
	6.1.1	Tube	Component	Volume	
		Light purple	Nuclease buffer (103-110-200)	5 µL	
		Light green	Nuclease mix (103-110-100)	5 µL	
		Total volume		10 µL	
	612	Add 10 ul of	Nucleopa Maatar r	miv to ooob	cample. The total volume should be 50 ul

6.1.2 Add 10  $\mu$ L of Nuclease Master mix to each sample. The total volume should be 50  $\mu$ L.



- 6.1.3 Thoroughly mix by pipetting up and down 10 times and then quick spin to collect all liquid.
- 6.1.4 Run the nuclease treatment program with the lid set to >47°C.

Step	Temperature	Duration
1	37°C	15 min
2	4°C	Hold

### 6.2 Final cleanup with SMRTbell cleanup beads

🖌 Step	Instructions
6.2.1	Add 50 $\mu$ L SMRTbell cleanup beads to each sample from the previous step. Pipette-mix the beads until evenly distributed.
6.2.2	Quick-spin the tube strip in a microcentrifuge to collect all liquid.
6.2.3	Leave at room temperature for 10 minutes to allow DNA to bind the beads.
6.2.4	Place the tube strip in a magnetic separation rack until the beads separate fully from the solution.
6.2.5	Slowly pipette off the cleared supernatant without disturbing the beads and discard the supernatant.
6.2.6	Slowly dispense 200 $\mu$ L, or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, pipette off the 80% ethanol and discard.
6.2.7	Repeat the previous step.
6.2.8	<ul> <li>Remove residual 80% ethanol:</li> <li>Remove the tube strip from the magnetic separation rack.</li> <li>Quick-spin the tube strip in a microcentrifuge.</li> <li>Place the tube strip back in a magnetic separation rack until beads separate fully from the solution.</li> <li>Pipette off residual 80% ethanol and discard.</li> </ul>
6.2.9	Remove the tube strip from the magnetic rack. Immediately add 20 $\mu$ L of elution buffer to each tube and resuspend the beads by pipetting 10 times or until evenly distributed.
6.2.1	0 Quick-spin the tube strip in a microcentrifuge to collect liquid.
6.2.1	1 Leave at room temperature for 5 minutes to elute DNA.
6.2.1	Place the tube strip in a magnetic separation rack until the beads separate fully from the solution.



6.2.13 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new PCR tube strip. Discard old tube strip with beads.

Take a 1 µL aliquot from each tube. Make a 1:5 dilution of the sample in elution buffer and measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit. Calculate the total mass. Expect 10-25% recovery of the starting Kinnex-PCR product.

6.2.14 total mass. Expect 10-25% recovery of the starting Kinnex-PCR product.
 Recommended: Further dilute each aliquot to 250 pg/µL with the Femto Pulse dilution buffer.
 Measure the final SMRTbell library size distribution with a Femto Pulse system.

Proceed to SMRT Link Sample Setup to prepare the SMRTbell library for sequencing. DNA

- 6.2.15 **concentration must be less than 60 ng/µl to go into ABC.** Using a concentration above 60 ng/µl will result in lower loading during sequencing.
- 6.2.16 Store SMRTbell libraries at 4°C if sequencing within the week. Long-term storage should be at -20°C. Minimize freeze-thaw cycles when handling SMRTbell libraries.

#### PROTOCOL COMPLETE

Revision history (description)	Version	Date
Initial release	01	October 2023
Minor updates throughout	02	December 2023
Updated to incorporate PCR program and tube color information into workflow steps	03	March 2024

Research use only. Not for use in diagnostic procedures. © 2024 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at <u>pacb.com/license</u>. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, and PureTarget are trademarks of PacBio.

© 2024 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures. 103-238-700 REV 03 MAR2024

