

Before you begin	Step 1. Repair and A-tailing	Step 2. Adapter ligation and cleanup																																																		
<ul style="list-style-type: none"> <li>Use <b>300–1000 ng</b> of DNA per SMRT® Cell 8M.</li> <li>Divide <b>1000 ng</b> by the number of samples when multiplexing with SMRTbell barcoded adapters.</li> <li>Iso-Seq® samples require <math>\geq 160</math> ng of cDNA per SMRT Cell 8M.</li> </ul> <table border="1" data-bbox="111 321 684 448"> <thead> <tr> <th>DNA shearing</th> <th>Human, plant, animal</th> <th>Microbial or metagenomic</th> </tr> </thead> <tbody> <tr> <td>Megarupter 3</td> <td>Speed 31</td> <td>Speed 40</td> </tr> <tr> <td>Modal size</td> <td>15–18 kb</td> <td>7–12 kb</td> </tr> </tbody> </table> <ul style="list-style-type: none"> <li>Resuspend samples in <b>47 <math>\mu</math>L</b> of low TE buffer using the SMRTbell cleanup bead steps prior to SMRTbell prep kit 3.0 steps.</li> <li>Adjust component volumes for the number of samples being prepared, plus 10% overage when preparing the reagent mixes.</li> <li>Add reagent components directly to sample(s) if not preparing reagent mixes. Pipette mix and spin down prior to incubation.</li> </ul>	DNA shearing	Human, plant, animal	Microbial or metagenomic	Megarupter 3	Speed 31	Speed 40	Modal size	15–18 kb	7–12 kb	<p><b>1.1</b> Make reaction mix 1 (RM1) in new tube.</p> <table border="1" data-bbox="737 175 1308 378"> <thead> <tr> <th colspan="3">Reaction mix 1 (RM1) per sample volumes</th> </tr> <tr> <th>Tube</th> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>Purple</td> <td>Repair buffer</td> <td>8 <math>\mu</math>L</td> </tr> <tr> <td>Blue</td> <td>End repair mix</td> <td>4 <math>\mu</math>L</td> </tr> <tr> <td>Green</td> <td>DNA repair mix</td> <td>2 <math>\mu</math>L</td> </tr> </tbody> </table> <p><b>1.2</b> Pipette mix and quick spin <b>RM1</b> to collect liquid.</p> <p><b>1.3</b> Add <b>14 <math>\mu</math>L</b> of the <b>RM1</b> to each sample.</p> <p><b>1.4</b> Incubate samples with the following program:</p> <table border="1" data-bbox="793 524 1140 683"> <thead> <tr> <th>Temperature</th> <th>Time</th> </tr> </thead> <tbody> <tr> <td>37°C</td> <td>30 minutes</td> </tr> <tr> <td>65°C</td> <td>5 minutes</td> </tr> <tr> <td>4°C</td> <td>Hold</td> </tr> </tbody> </table> <p><b>1.5</b> Proceed to next step.</p>	Reaction mix 1 (RM1) per sample volumes			Tube	Component	Volume	Purple	Repair buffer	8 $\mu$ L	Blue	End repair mix	4 $\mu$ L	Green	DNA repair mix	2 $\mu$ L	Temperature	Time	37°C	30 minutes	65°C	5 minutes	4°C	Hold	<p><b>2.1</b> Add <b>4 <math>\mu</math>L</b> of SMRTbell adapter (barcoded or non-barcoded) to each sample.</p> <p><b>2.2</b> Make reaction mix 2 (RM2) in new tube.</p> <table border="1" data-bbox="1377 264 1934 431"> <thead> <tr> <th colspan="3">Reaction mix 2 (RM2) per sample volumes</th> </tr> <tr> <th>Tube</th> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>Yellow</td> <td>Ligation mix</td> <td>30 <math>\mu</math>L</td> </tr> <tr> <td>Red</td> <td>Ligation enhancer</td> <td>1 <math>\mu</math>L</td> </tr> </tbody> </table> <p><b>2.3</b> Pipette mix and quick spin <b>RM2</b> to collect liquid.</p> <p><b>2.4</b> Add <b>31 <math>\mu</math>L</b> of <b>RM2</b> to each sample.</p> <p><b>2.5</b> Incubate samples using the following steps:</p> <table border="1" data-bbox="1415 553 1764 672"> <thead> <tr> <th>Temperature</th> <th>Time</th> </tr> </thead> <tbody> <tr> <td>20°C</td> <td>30 minutes</td> </tr> <tr> <td>4°C</td> <td>Hold</td> </tr> </tbody> </table> <p><b>2.6</b> Clean up with SMRTbell cleanup beads and resuspend in <b>40 <math>\mu</math>L</b> of elution buffer.</p>	Reaction mix 2 (RM2) per sample volumes			Tube	Component	Volume	Yellow	Ligation mix	30 $\mu$ L	Red	Ligation enhancer	1 $\mu$ L	Temperature	Time	20°C	30 minutes	4°C	Hold
DNA shearing	Human, plant, animal	Microbial or metagenomic																																																		
Megarupter 3	Speed 31	Speed 40																																																		
Modal size	15–18 kb	7–12 kb																																																		
Reaction mix 1 (RM1) per sample volumes																																																				
Tube	Component	Volume																																																		
Purple	Repair buffer	8 $\mu$ L																																																		
Blue	End repair mix	4 $\mu$ L																																																		
Green	DNA repair mix	2 $\mu$ L																																																		
Temperature	Time																																																			
37°C	30 minutes																																																			
65°C	5 minutes																																																			
4°C	Hold																																																			
Reaction mix 2 (RM2) per sample volumes																																																				
Tube	Component	Volume																																																		
Yellow	Ligation mix	30 $\mu$ L																																																		
Red	Ligation enhancer	1 $\mu$ L																																																		
Temperature	Time																																																			
20°C	30 minutes																																																			
4°C	Hold																																																			
Step 3. Nuclease treatment and cleanup	SMRTbell cleanup bead steps	AMPure PB bead size selection steps																																																		
<p><b>3.1</b> Make reaction mix 3 (RM3) in new tube.</p> <table border="1" data-bbox="111 865 657 1027"> <thead> <tr> <th colspan="3">Reaction mix 3 (RM3) per sample volumes</th> </tr> <tr> <th>Tube</th> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>Light purple</td> <td>Nuclease buffer</td> <td>5 <math>\mu</math>L</td> </tr> <tr> <td>Light green</td> <td>Nuclease mix</td> <td>5 <math>\mu</math>L</td> </tr> </tbody> </table> <p><b>3.2</b> Pipette mix <b>RM3</b> and quick spin to collect liquid.</p> <p><b>3.3</b> Add <b>10 <math>\mu</math>L</b> of <b>RM3</b> to each sample.</p> <p><b>3.4</b> Incubate samples using the following steps:</p> <table border="1" data-bbox="170 1149 516 1268"> <thead> <tr> <th>Temperature</th> <th>Time</th> </tr> </thead> <tbody> <tr> <td>37°C</td> <td>15 minutes</td> </tr> <tr> <td>4°C</td> <td>Hold</td> </tr> </tbody> </table> <p><b>3.5</b> Clean up samples with SMRTbell cleanup or AMPure® PB beads depending on chosen size selection workflow.</p>	Reaction mix 3 (RM3) per sample volumes			Tube	Component	Volume	Light purple	Nuclease buffer	5 $\mu$ L	Light green	Nuclease mix	5 $\mu$ L	Temperature	Time	37°C	15 minutes	4°C	Hold	<ol style="list-style-type: none"> <li>Add the appropriate volume of beads to each sample (e.g., 1X for WGS and 1.3X for amplicons).</li> <li>Mix and bind for <b>10 minutes</b> at room temp.</li> <li>Place samples on magnetic stand, allow beads to pellet, and remove supernatant.</li> <li>Wash bead pellet with freshly prepared 80% EtOH. Repeat for a total of 2 washes.</li> <li>Remove all residual EtOH.</li> <li>Remove samples from magnet and immediately resuspend beads in the appropriate volume for the next step of the library prep.</li> <li>Mix and leave at room temp for <b>5 minutes</b>.</li> <li>Place samples back on magnet, allow beads to pellet, and transfer supernatant to new tube.</li> </ol>	<ol style="list-style-type: none"> <li>Prepare a <b>35% (v/v)</b> dilution of AMPure PB beads or use previously prepared dilution.</li> <li>Add <b>3.1X (v/v)</b> of diluted beads to each sample.</li> <li>Mix and leave at room temp for <b>20 minutes</b>.</li> <li>Place samples on magnetic stand, allow beads to pellet, and remove supernatant.</li> <li>Wash bead pellet with freshly prepared 80% EtOH. Repeat for a total of 2 washes.</li> <li>Remove all residual EtOH.</li> <li>Remove sample from magnet and immediately resuspend beads in <b>15 <math>\mu</math>L</b> of elution buffer.</li> <li>Mix and leave at room temp. for <b>5 minutes</b>.</li> <li>Place samples back on magnet, allow beads to pellet, and transfer supernatant to an Eppendorf DNA LoBind or PCR strip tube.</li> </ol>																																
Reaction mix 3 (RM3) per sample volumes																																																				
Tube	Component	Volume																																																		
Light purple	Nuclease buffer	5 $\mu$ L																																																		
Light green	Nuclease mix	5 $\mu$ L																																																		
Temperature	Time																																																			
37°C	15 minutes																																																			
4°C	Hold																																																			

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable [license terms](#). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omnione, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, and SBB are trademarks of Pacific Biosciences of California, Inc. (PacBio). All other trademarks are the sole property of their respective owners.