




SMRT[®] Link Python API reference (v13.0)



Research use only. Not for use in diagnostic procedures.

P/N 103-285-100 Version 01 (December 2023)

© 2023 Pacific Biosciences of California, Inc. ("PacBio")

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document.

PACBIO DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESS, STATUTORY, IMPLIED OR OTHERWISE, INCLUDING, BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY, SATISFACTORY QUALITY, NONINFRINGEMENT OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL PACBIO BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, PURSUANT TO ANY STATUTE, OR ON ANY OTHER BASIS FOR SPECIAL, CONSEQUENTIAL, INCIDENTAL, EXEMPLARY OR INDIRECT DAMAGES IN CONNECTION WITH (OR ARISING FROM) THIS DOCUMENT, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT PACBIO IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license.

Trademarks:

Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omnione, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, and Kinnex are trademarks of Pacific Biosciences of California Inc. (PacBio).

See <https://github.com/broadinstitute/cromwell/blob/develop/LICENSE.txt> for Cromwell redistribution information.

PacBio
1305 O'Brien Drive
Menlo Park, CA 94025
www.pacb.com



SMRT[®] Link Python API reference (v13.0)

Introduction	3
Viewing API information.....	3
Python reference client	3
Python reference client setup instructions.....	4
How to query one or more runs and their details.....	4
How to retrieve all dataset IDs and movie names produced by a run.....	6
How to query all the dataset details for a top-level dataset and all of its child datasets	7
How to query P0, P1, and P2 metrics.....	9
How to export demultiplexed datasets	9
How to query datasets.....	12
How to terminate a SMRT Link job.....	12
How to get all created, running, successful, failed, or terminated analysis jobs	13
How to poll for a job to successfully complete within a specified timeout	14
How to query a Run XML.....	14
How to query a run design.....	14
How to create a run design by importing a run design CSV	15
How to find the Run QC reports associated with an analysis job	16
How to find the contents of specific Run QC reports	19
How to query instrument status, such as Running or Complete.....	19
How to query a job.....	19
How to start a SMRT Link job for a specific workflow	20
How to combine a sample split across multiple cells.....	21
How to poll every 10 minutes until a collection is complete, then launch a HiFi Mapping job using the official PacBio hg38 reference, and poll until it completes successfully	21
How to find information about the Run QC reports associated with an analysis job	22
How to export demultiplexed dataset metrics as a CSV file	22

Introduction

The SMRT Link web services API, provided by PacBio, allows integration of SMRT Link with third-party software. It is also used for accessing features such as designing and performing quality control on instrument runs, querying new data from the instrument, and starting analyses on the sequence data.

This document includes Python examples for using the API to perform the tasks listed in the Table of Contents.

Viewing API information

For detailed information on the SMRT Link web services API calls, including definitions and examples, see Swagger-generated information here:

```
https://<server name>:8243/sl/docs/services/#/default
```

where `<server name>:8243` is the name and port number of your local SMRT Link server.

Python reference client

SMRT Link v13.0 includes a new Python reference client covering the most common API methods. The Python reference client provides concise and self-documenting examples of how to use the PacBio API, and is easily portable to other languages. As the Python reference client is **not** a full SDK, it deals only with simple Python types (`int`, `str`, `float`, `list`, `dict`) rather than more complex objects.

The code for the Python reference client can be downloaded on github [here](#).

The Python reference client is redistributable as a standalone library with no non-standard dependencies (only the `requests` library available on PyPI).

The Python reference client is **required** to use the Python examples in this document. The Python reference client includes methods for all of the following API calls:

Runs:

- Get run(s)
- Get run xml
- Get run collection(s)
- Get run from collection
- Get run collection reports
- Get run collection barcodes
- Get run collection [hifi](#) reads
- Get run collection [hifi](#) reads barcoded datasets
- Get run reports
- Get run design
- Import run design CSV
- Delete run
- Import run xml
- Update run xml

Chem bundle:

- Get active bundle metadata
- Get chemistry bundle metadata
- Get active bundle file
- Get chemistry bundle file

Jobs:

- Get job
- Get job reports
- Download job report resources
- Get job datastore
- Get job entry points
- Get job datasets
- Get job options
- Download job datastore file
- Get analysis jobs
- Get analysis jobs by state
- Get analysis jobs by parent
- Get SMRT Analysis nested jobs
- Create analysis job
- Terminate analysis job
- Get import dataset jobs
- Create import dataset [zip] job
- Create import collection job
- Create merge datasets job
- Get pipeline(s)
- Poll for successful job

Datasets:

- Get consensus read sets
- Get consensus read sets by movie
- Get barcoded child datasets
- Get subread sets
- Get references
- Get barcode sets
- Get consensus read set
- Get subread set
- Get reference set
- Get barcode set
- Get consensus read set reports
- Get barcode set contents
- Get barcode set record names
- Get dataset metadata
- Get dataset jobs
- Get dataset search

Python reference client setup instructions

If you downloaded the Python reference client to your local directory:

1. Make sure the `requests` module is installed for Python3. This will vary depending on your system, but `pip3 install requests` is one option.
2. Start the Python3 interpreter. (Version 3.9 or later is required.)
3. Import the Python reference client:


```
from smrtlink_client import SmrtLinkClient
```
4. Set up the Python reference client with your SMRT Link server:


```
client = SmrtLinkClient.connect("<your-SMRT-Link-server-name>", "smrt-link-username", "password")
```

If you have SMRT Link v13.0 installed, do the following:

1. Start the Python3 interpreter bundled with SMRT Link:


```
$SMRT_ROOT/smrtcnds/bin/python3
```
2. Import the Python reference client:


```
from pbcommand.services.smrtlink_client import SmrtLinkClient
```
3. Set up the Python reference client with your SMRT Link server:


```
client = SmrtLinkClient.connect("<your-SMRT-Link-server-name>", "smrt-link-username", "password")
```

Security note: We recommend that you use a dedicated API client user that is **separate** from your usual network-wide login. SMRT Link administrators can define local users that only have access to the GUI and API **without** being able to log in to any other system.

How to query one or more runs and their details

You can find the run ID by using the `runs` endpoint and specifying a query to search for the desired run:

```
GET /smrt-link/runs
```

Following is an example of how to do this using Python:

```
runs = client.get_runs(name=None, reserved=None, instrumentType=None,
                       chipType=None, collectionUuid=None, movieName=None)
```

Example request:

```
runs = client.get_runs(instrumentType = "Revio")
```

In this example, you can query a run based on the following parameters:

- `name (str)`: Filter by run name; partial matches are supported.
- `reserved (bool)`: Filter by reservation status; set to `true` for runs selected on instrument.
- `instrumentType (str)`: Filter by instrument type (Revio, Sequel2e, or Sequel2).
- `chipType (str)`: Filter by chip type (8mChip or 25mChip).
- `collectionUuid (str)`: Retrieve the run for a specific collection UUID.
- `movieName (str)`: Filter by movie name associated with one of the runs.

Following is an example of what this returns:

```
{
  "reserved": true, "numLRCells": 0,
  "name": "20230414_Revio",
  "completedAt": "2023-04-17T04:19:50.718Z",
  "chemistrySwVersion": "1.2.3.11111",
  "instrumentType": "Revio",
  "chipType": "25mChip",
  "instrumentName": "12345e",
  "context": "r1234_20230414_212018",
  "instrumentSwVersion": "1.1.1.11111",
  "numCellsCompleted": 1,
  "totalCells": 1,
  "primaryAnalysisSwVersion": "12.0.0.0",
  "status": "Complete",
  "numStandardCells": 1,
  "createdAt": "2023-04-14T18:17:23.849Z",
  "startedAt": "2023-04-14T21:20:44.057Z",
  "createdBy": "bobsmith",
  "totalSamples": 97,
  "numCellsFailed": 0,
  "plate2": "1234567800301580073920231018",
  "instrumentSerialNumber": "12345",
  "transfersCompletedAt": "2023-04-18T00:13:42.032Z",
  "uniqueId": "0a12bcd2-5916-4142-9b82-220a7bb04d13",
  "ccsExecutionMode": "OnInstrument"
}
```

- reserved means that the run design was used on an instrument. If false, that means the run design was created in SMRT Link, but has not yet been used in a run.
- createdAt is the time the run design was created.
- startedAt is the time the user pressed the button to close the Revio instrument door and begin a run.
- completedAt is the time when all of the instrument's movies have finished acquiring and after the cleanup; that is, instrument state = Complete, Aborted, or Terminated. The cleanup is skipped if another run has been pre-loaded.
- uniqueID is the Run UUID you see in the SMRT Link URL on the Run Details page:

<https://my-smrtlink-server:8243/sl/runs/12a3bc45-1299-49c6-9f84-9ceac118ce9e>

- context is the Run ID you see in the SMRT Link GUI on the Run Details page:

The screenshot shows the PacBio SMRT Link Run Details page for a run named "My_Run_Name". The run status is "COMPLETE". The "Overview" section contains the following details:

Run Created: 2023-10-07, 10:17:27 PM	Completed Cells: 4	Run ID: r12345e_20231008_191807
Run Start: 2021-10-08, 02:23:02 PM	Failed Cells: 0	Instrument SN: 12345e
Run Complete: 2023-10-10, 02:28:14 AM	Time remaining for PostProcessing: -	Instrument Control SW Version: 11.1.1.11111
Created By: bsmith	Transfer Complete: 2023-10-10, 05:01:48 AM	Instrument Chemistry Bundle: 11.1.1.11111
Instrument Name: my-instrument		Primary SW Version: 11.1.1.11111

How to retrieve all dataset IDs and movie names produced by a run

You can find all the datasets produced by a run by using the `runs` endpoint and specifying a query to search for collections with a matching run ID:

```
GET /smrt-link/runs/{run_id}/collections
```

Following is an example of how to do this using Python:

```
collections = client.get_run_collections(run_id = "<run_id>")
```

Following is an example of what this returns:

```
[
  {
    "name": "Sequel2e_hifi_20231010_cell14",
    "completedAt": "2023-10-10T12:01:48.810Z",
    "instrumentName": "12345e",
    "context": "m64263e_211009_220021",
    "well": "D01",
    "projectId": 1,
    "sequencingKit": "123456101826100062822",
    "status": "Complete",
    "instrumentId": "12345e",
    "startedAt": "2021-10-09T22:00:22.030Z",
    "createdBy": "bobsmith",
    "cellType": "Standard",
    "uniqueId": "702d5232-9c2d-4f92-a35a-784949f9db0e",
    "collectionPathUri": "/data/hifiviral/r64263e_20211008_191807/4_D01",
    "ccsExecutionMode": "OnInstrument",
    "runId": "186f258a-dac7-479d-a122-aa0e6bbcab57",
    "ccsId": "9dbcaaf6-298a-465a-ad56-3175d1151c57",
    "movieMinutes": 480.0
  },
  {
    "name": "Sequel2e_hifiviral_20211010_cell14",
    "completedAt": "2021-10-10T04:39:54.955Z",
    "instrumentName": "12345e",
    "context": "m64263e_211009_134807",
    "well": "C01",
    "projectId": 1,
    "sequencingKit": "122894101826100062822",
    "status": "Complete",
    "importedAt": "2023-08-25T08:18:42.900Z",
    "instrumentId": "64263e",
    "startedAt": "2021-10-09T13:48:08.677Z",
    "createdBy": "bobsmith",
    "cellType": "Standard",
    "uniqueId": "c2690435-3336-4cb2-9bc1-b509d6901bbc",
    "collectionPathUri": "/data/hifiviral/r64263e_20211008_191807/3_C01",
    "ccsExecutionMode": "OnInstrument",
    "runId": "186f258a-dac7-479d-a122-aa0e6bbcab57",
    "ccsId": "ca24d2be-dae-480f-b3e6-c352104aff4d",
    "movieMinutes": 480.0
  }
]
```

```

},
...
]

```

- `ccsId` is the Dataset UUID that is visible at the top level of Data Management.
- `runId` is the Run ID you see in the URL in SMRT Link on the Runs/Run Details page:



- `context` is the movie name that is part of the output dataset you see on the Run Details page in SMRT Link.

How to query all the dataset details for a top-level dataset and all of its child datasets

Following is an example of how to do this using Python:

```

parent_dataset = client.get_consensusreadset(dataset_id = <dataset_id>)
child_datasets = client.get_barcode_child_datasets(parent_dataset_id = <ccsID>,
                                                    barcode_name=None,
                                                    biosample_name=None)

```

You can use this, for example, to find all the UUID's of your child datasets.

For Revio or Sequel IIe runs demultiplexed on-instruments, the parent dataset ID is the same as the `ccsId` referenced in the `collections` endpoint.

Following is an example of what this returns:

```

{
  "instrumentControlVersion": "1.1.1.111111",
  "tags": "ccs,testdata,barcoded",
  "instrumentName": "my-instrument",
  "uuid": "d717f468-e8f2-44cc-9a25-a50d699319b0",
  "dnaBarcodeName": "bc2082--bc2082",
  "totalLength": 1003747,
  "projectId": 1,
  "numRecords": 160,
  "wellSampleName": "My-sample_A01",
  "bioSampleName": "EC_40_82",
  "version": "3.0.1",
  "cellId": "DA120589",
  "id": 286,
  "md5": "c93c21ab91289d707e7adf31df57a949",
  "parentName": "My-sample_A01-Cell11 (CCS) (all samples)",
  "importedAt": "2023-08-25T12:16:50.898Z",
  "jobId": 164,
  "createdAt": "2023-08-25T12:13:48.489Z",
  "isActive": true,
  "createdBy": "bsmith",
  "wellName": "A01",
  "cellIndex": 0,
  "parentUuid": "f6dba83f-36b5-4443-b6fd-c8730ae350ad",
  "metadataContextId": "m64004_210910_192345",

```



```

    "numChildren": 0,
    "numResources": 1,
    "runName": "20231110_my_run",
    "datasetType": "PacBio.DataSet.ConsensusReadSet",
    "comments": "ccs dataset converted"
  },
  {
    "name": "My-sample_A01-Cell11 (CCS) (BS_43_85)",
    "updatedAt": "2023-08-25T12:13:48.560Z",
    "path": "<path>",
    "instrumentControlVersion": "1.1.1.111111",
    "tags": "ccs,testdata,barcoded",
    "instrumentName": "my-instrument",
    "uuid": "65ba42ac-a5f1-4f43-be2a-a26d32118dc4",
    "dnaBarcodeName": "bc2085--bc2085",
    "totalLength": 1176220,
    "projectId": 1,
    "numRecords": 146,
    "wellSampleName": "My-sample_A01",
    "bioSampleName": "BS_43_85",
    "version": "3.0.1",
    "cellId": "DA120589",
    "id": 285,
    "md5": "90607244e2d16c7035d6b652aa6be68a",
    "parentName": "My-sample_A01-Cell11 (CCS) (all samples)",
    "importedAt": "2023-08-25T12:16:50.851Z",
    "jobId": 164,
    "createdAt": "2023-08-25T12:13:48.560Z",
    "isActive": true,
    "createdBy": "bsmith",
    "wellName": "A01",
    "cellIndex": 0,
    "parentUuid": "f6dba83f-36b5-4443-b6fd-c8730ae350ad",
    "metadataContextId": "m64004_210910_192345",
    "numChildren": 0,
    "numResources": 1,
    "runName": "20231110_my_run",
    "datasetType": "PacBio.DataSet.ConsensusReadSet",
    "comments": "ccs dataset converted"
  },
  {
    "name": "My-sample_A01-Cell11 (CCS) (BS_40_81)",
    "updatedAt": "2023-08-25T12:13:48.456Z",
    "path": " <path>.consensusreadset.xml",
    "instrumentControlVersion": "1.1.1.111111",
    "tags": "ccs,testdata,barcoded",
    "instrumentName": "my-instrument",
    "uuid": "2b94dbb5-0be6-42b9-aa6c-a630c48a0e2d",
    "dnaBarcodeName": "bc2081--bc2081",
    "totalLength": 780298,
    "projectId": 1,
    "numRecords": 87,
    "wellSampleName": "My-sample_A01",
    "bioSampleName": "BS_40_81",
    "version": "3.0.1",

```

```

"cellId": "DA120589",
"id": 284,
"md5": "2e34c910f9e9856d8215382024e34109",
"parentName": "My-sample_A01-Cell11 (CCS) (all samples)",
"importedAt": "2023-08-25T12:16:50.804Z",
"jobId": 164,
"createdAt": "2023-08-25T12:13:48.456Z",
"isActive": true,
"createdBy": "bsmith",
"wellName": "A01",
"cellIndex": 0,
"parentUuid": "f6dba83f-36b5-4443-b6fd-c8730ae350ad",
"metadataContextId": "m64004_210910_192345",
"numChildren": 0,
"numResources": 1,
"runName": "20231110_my_run ",
"datasetType": "PacBio.DataSet.ConsensusReadSet",
"comments": "ccs dataset converted"
},
...
]

```

How to query P0, P1, and P2 metrics

You can query P0, P1, and P2 loading metrics using the SMRT Link API. The following example searches run collection reports for the Loading Report and extracts the P0, P1, and P2 values from it. You need the Run UUID and Collection UUID.

Following is an example of how to do this using Python:

```

reports = client.get_run_collection_reports(<run_uuid>, <collection_uuid>)
for r in reports:
    if "loading" in r["reportTypeId"]:
        report_uuid = r["dataStoreFile"]["uuid"]
        report = client.load_datastore_report_file(report_uuid)
        print({attr["id"]:attr["value"] for attr in report["attributes"]})

```

Following is an example of what this returns:

```

{'loading_xml_report.productive_zmws': 8014671,
'loading_xml_report.productivity_0_n': 808345,
'loading_xml_report.productivity_1_n': 6161738,
'loading_xml_report.productivity_2_n': 1044588}

```

How to export demultiplexed datasets

You can retrieve the demultiplexed "child" datasets for a PacBio instrument run, optionally filtering by barcode name (such as bc2001--bc2001) or biosample name.

Note: Biosample name = well name = well sample name = collection name.

Following is an example of how to do this using Python:

```

#first get the collection_id
collection_id = client.get_run_collections(<run_id>)

```

This returns:

```
{
  "name": "20230414_84026_16kbHG002_97barcodes ",
  "completedAt": "2023-04-18T00:13:42.032Z",
  "instrumentName": "my-instrument",
  "context": "m84026_230415_224020_s3",
  "well": "C01",
  "projectId": 1,
  "sequencingKit": "030158102118800101823",
  "labelNumber": "00739",
  "status": "Complete",
  "importedAt": "2023-04-17T23:55:12.727Z",
  "instrumentId": "84026",
  "startedAt": "2023-04-15T22:40:24.006Z",
  "createdBy": "admin",
  "cellType": "Standard",
  "uniqueId": "69b09865-5c23-4717-92ad-75968a43f443",
  "collectionPathUri": "/collections/349/r84026_20230414_212018/2_C01",
  "ccsExecutionMode": "OnInstrument",
  "runId": "0f99fea6-5916-4142-9b82-220a7bb04d13",
  "ccsId": "3acce2c3-d904-4d08-aba0-2628d0dcccbf",
  "movieMinutes": 1440.0
}
```

#uniqueId is the collection_id

```
demuxed_datasets = client.get_run_collection_hifi_reads_barcode_datasets(<run_id>,
<collection_id>,
barcode_name=None,
biosample_name=None)
```

Following is an example of what this returns:

```
[
  {
    "name": "20230414_84026_225pM_97Barcodes_bsmith",
    "updatedAt": "2023-04-17T23:55:12.721Z",
    "path": "<path>.consensusreadset.xml",
    "instrumentControlVersion": "1.1.1.11111",
    "tags": "ccs,barcoded",
    "instrumentName": "my-instrument",
    "uuid": "ab1234cd-5efg-7h99-0000-0a1209170a1b",
    "dnaBarcodeName": "bc2096--bc2096",
    "totalLength": 891005473,
    "projectId": 1,
    "numRecords": 56787,
    "wellSampleName": "my_wellsamplename",
    "bioSampleName": "my_biosamplename",
    "version": "3.0.1",
    "cellId": "1000000475102202200071423",
    "id": 132,
    "md5": "33c8697d37a3feb9fbb844c913ea916c",
    "parentName": "20230414_84026_225pM_97Barcodes",
    "importedAt": "2023-04-17T23:55:12.721Z",
    "jobId": 27,
```

```

    "createdAt": "2023-04-17T23:55:12.721Z",
    "isActive": true,
    "createdBy": "bsmith",
    "wellName": "C01",
    "cellIndex": 6,
    "parentUuid": "3acce2c3-d904-4d08-aba0-2628d0dcccbf",
    "metadataContextId": "m12345_230415_224020_s3",
    "numChildren": 0,
    "numResources": 1,
    "runName": "20230414_84026_16kbHG002_97Barcodes",
    "datasetType": "PacBio.DataSet.ConsensusReadSet",
    "comments": "Record generated by runqc-reports"
  },
  {
    "name": "20230414_84026_16kbHG002_97Barcodes-Cell17 (225pM_16kb_HG002)",
    "updatedAt": "2023-04-17T23:55:12.647Z",
    "path": "<path>.consensusreadset.xml",
    "instrumentControlVersion": "1.1.1.11111",
    "tags": "ccs,barcoded",
    "instrumentName": "my-instrument",
    "uuid": "7960507a-437a-455b-a1b0-9cc6a9c2386d",
    "dnaBarcodeName": "bc2095--bc2095",
    "totalLength": 553068822,
    "projectId": 1,
    "numRecords": 37200,
    "wellSampleName": "20230414_84026_16kbHG002_97Barcodes",
    "bioSampleName": "225pM_16kb_HG002",
    "version": "3.0.1",
    "cellId": "1000000475102202200071423",
    "id": 131,
    "md5": "f21fcf132a20e009b13836cfc403ce15",
    "parentName": "20230414_84026_16kbHG002_97Barcodes-Cell17 (all samples)",
    "importedAt": "2023-04-17T23:55:12.647Z",
    "jobId": 27,
    "createdAt": "2023-04-17T23:55:12.647Z",
    "isActive": true,
    "createdBy": "admin",
    "wellName": "C01",
    "cellIndex": 6,
    "parentUuid": "3acce2c3-d904-4d08-aba0-2628d0dcccbf",
    "metadataContextId": "m84026_230415_224020_s3",
    "numChildren": 0,
    "numResources": 1,
    "runName": "20230414_84026_16kbHG002_97Barcodes",
    "datasetType": "PacBio.DataSet.ConsensusReadSet",
    "comments": "Record generated by runqc-reports"
  },
  ...
]

```

How to query datasets

You can retrieve a list of HiFi datasets, with optional search parameters.

Following is an example of how to do this using Python:

```
client.get_consensusreadsets(name=None,
                             bioSampleName=None,
                             wellSampleName=None,
                             metadataContextId=None)
```

Following is a partial list of supported search terms:

- name (dataset name)
- bioSampleName
- wellSampleName
- metadataContextId (movie name)

You can find the name (dataset name), wellSampleName, and bioSampleName in the SMRT Link GUI here:



You can find the Movie Name by clicking the dataset name and viewing the Dataset Details page.

Search term notes:

- String searches are always **case-insensitive**.
- Most of the non-timestamp string fields in the data model are searchable using partial strings by adding the prefix `like:` to the search term, such as `client.get_consensusreadsets(bioSampleName="like:HG002")`
- The prefixes `not:` (inequality), `unlike:`, `start:` and `end:` are also supported.
- For numerical fields, `not:`, `lt:`, `lte:`, `gt:`, and `gte:` are supported, as well as `range:{start},{end}`.

How to terminate a SMRT Link job

You can immediately terminate an analysis job.

Following is an example of how to do this using Python:

```
client.terminate_analysis_job(job_id)
```

This returns:

```
{'message': 'Cromwell workflow <uuid> TERMINATED for <job_id>'}
```

How to get all created, running, successful, failed, or terminated analysis jobs

You can search for all analysis jobs of a specific state: CREATED, SUBMITTED, RUNNING, SUCCESSFUL, FAILED, TERMINATED, or ABORTED.

Following is an example of how to do this using Python:

```
terminated_jobs = client.get_analysis_jobs_by_state(state)
```

Following is an example that returns all jobs that were terminated:

```
terminated_jobs = client.get_analysis_jobs_by_state(state = 'TERMINATED')
```

Following is an example of what this returns:

```
[{
  "subJobTypeId": "cromwell.workflows.pb_assembly_hifi",
  "name": "test",
  "updatedAt": "2023-09-20T07:51:05.637Z",
  "workflow": "{}",
  "path": "<path>",
  "state": "TERMINATED",
  "tags": "",
  "uuid": "1b717f91-6651-45d3-b8cf-e1b4b4427f6a",
  "externalJobId": "547075f8-674d-4fbb-a0ba-69f55c627892",
  "jobStartedAt": "2023-09-20T07:50:10.373Z",
  "applicationName": "Genome Assembly",
  "projectId": 1,
  "childJobsCount": 0,
  "jobCompletedAt": "2023-09-20T07:51:05.637Z",
  "jobTypeId": "analysis",
  "id": 627,
  "smrtlinkVersion": "1.1.1.11111",
  "comment": "Description for job Run Analysis Application",
  "isNested": false,
  "createdAt": "2023-09-20T07:50:05.567Z",
  "isActive": true,
  "createdBy": "bsmith",
  "createdByEmail": "bsmith@company.com",
  "isMultiJob": false,
  "jsonSettings": "",
  "jobUpdatedAt": "2023-09-20T07:51:05.637Z"
}]
```

How to poll for a job to successfully complete within a specified timeout

You can poll a submitted job of any type until it completes successfully within the specified timeout, or raise an exception.

Following is an example of how to do this using Python:

```
client.poll_for_successful_job(job_id, sleep_time, max_time):
```

Following is an example of what this returns:

```
finished_job = client.poll_for_successful_job(job_id=<job_id>,
                                             sleep_time=120,
                                             max_time=28800)
```

How to query a Run XML

You can retrieve the XML data model for a PacBio instrument run.

Following is an example of how to do this using Python:

```
run_xml = client.get_run_xml(run_id)
```

This returns the XML as a string; this is the primary format that SMRT Link uses to send instructions to the instrument, and much of it will end up in the output dataset XMLs (with further modifications from the instrument software).

How to query a run design

Following is an example of how to do this using Python:

```
run_design = client.get_run_design(run_id)
```

Following is an example of what this returns for an AAV run:

```
{
  "chipType": "25mChip",
  "createdBy": "bsmith",
  "experimentDescription": "",
  "experimentId": "",
  "experimentName": "",
  "instrumentType": "Revio",
  "runDescription": "",
  "runName": "Bsmith Run 09.28.2023 18:27",
  "samples": [
    {
      "application": "AAVAnalysis",
      "includeCpG": true,
      "isBarcoded": true,
      "sampleDescription": "",
      "readSegmentation": false,
      "demuxMode": "OnInstrument",
      "uuid": "9aaee069-8ecf-4b4f-8f37-d0bd6a3abb5a",
      "ccsUuid": "0fab0fc2-4c13-496c-aa06-352a30539e6f",
      "includeLowQuality": true,
    }
  ]
}
```

```

"barcodeCsvFileName": "",
"sampleName": "test-well",
"ccsMode": "OnInstrument",
"includeKinetics": false,
"useDynamicLoading": true,
"loadingConcentration": 0.0,
"projectId": 1,
"adapter": "mas16",
"barcodedSamples": [
  {
    "bioSampleName": "Bio Sample 1",
    "dnaBarcodeName": "bc2001--bc2001",
    "uuid": "560148e0-c4a4-4512-9e2c-2c3a89a80f6d"
  }
],
"scIsoSeq": false,
"controlKit": "Lxxxxx102798000123199",
"templatePrepKit": "Lxxxxx999999001123199",
"movieTimeHours": 24.0,
"barcodesFasta":
">bc2001\nATCGTGCGACGAGTAT\n>bc2002\nTGCATGTTCATGAGTAT\n>bc2003\nACGAGTGCTCGAGTAT\n>bc200
4\nTGCAGTGTCTCGAGTAT\n>bc2005\nTGACTCGATCGAGTAT\n>bc2006\nCATGCGATCTGAGTAT\n>...\n",
"libraryType": "AAV",
"automationParameters": [
  {
    "name": "MovieLength",
    "valueDataType": "Double",
    "simpleValue": "1440"
  }
],
"sequencingKit": "000000102118800110723",
"labelNumber": "12345",
"primaryAutomationName": "",
"emitSubreadsPercent": 0,
"copyFiles": [],
"barcodeUuid": "43f950a9-8bde-3855-6b25-c13368069745",
"insertSize": 500,
"wellName": "A01",
"consensusMode": "strand",
"symmetricBarcodes": true,
"cellType": "human",
"minBarcodeScore": 80,
"heteroduplexDetection": false,
"plateNumber": 1,
"bindingKit": "Lxxxxx102739100123123"
}
],
"uuid": "ed271e1b-efa4-4999-ab7f-a945350348c3"
}

```

How to create a run design by importing a run design CSV

Following is an example of how to do this using Python:

```
client.import_run_design_csv(csv_file)
```


How to find the Run QC reports associated with an analysis job

You can obtain **all** collection-level reports associated with a run. (**Note:** This was introduced in SMRT Link v13.0.)

Following is an example of how to do this using Python:

```
run_reports = client.get_run_reports(run_id)
```

This returns information about the following reports, including when the report was created and the path to the report. (**Note:** To get the **content** of the reports, see **How to find the contents of specific Run QC reports.**)

- import_dataset.report_detect_cpg_methyl
- collection.barcode_preview_report
- import_dataset.report_adapters
- import_dataset.report_loading
- import_dataset.report_raw_data
- import_dataset.report_barcode
- import_dataset.report_control
- import_dataset.report_ccs2
- import_dataset.report_ccs_processing

Following is an example of what this returns:

```
[
  {
    "dataStoreFile": {
      "createdAt": "2023-09-28T08:21:56.426Z",
      "description": "detect_cpg_methyl",
      "fileSize": 1420,
      "fileTypeId": "PacBio.FileTypes.JsonReport",
      "importedAt": "2023-09-28T08:22:03.997Z",
      "isActive": true,
      "modifiedAt": "2023-09-28T08:21:56.426Z",
      "name": "Report detect_cpg_methyl",
      "path": "<directories>/reports/detect_cpg_methyl.report.json",
      "sourceId": "import_dataset.report_detect_cpg_methyl",
      "uuid": "4229c152-4e49-4199-adaa-d6d91ffd83b3"
    },
    "reportTypeId": "import_dataset.report_detect_cpg_methyl"
  },
  {
    "dataStoreFile": {
      "createdAt": "2023-09-28T08:16:51.670Z",
      "description": "PacBio Report barcode_preview_report (6ad9e5df-5aaf-44d3-9499-5ad7a0cc163d)",
      "fileSize": 5381,
      "fileTypeId": "PacBio.FileTypes.JsonReport",
      "importedAt": "2023-09-28T08:16:51.683Z",
      "isActive": true,
      "modifiedAt": "2023-09-28T08:16:51.670Z",
      "name": "Barcode Preview Report",
      "path": "<path>/barcode_preview_report.report.json",
      "sourceId": "collection.barcode_preview_report",
      "uuid": "6ad9e5df-5aaf-44d3-9499-5ad7a0cc163d"
    }
  }
]
```

```

    },
    "reportTypeId": "collection.barcode_preview_report"
  },
  {
    "dataStoreFile": {
      "createdAt": "2023-09-28T08:21:56.448Z",
      "description": "PacBio Report adapter_xml_report (d340b519-6fdc-401a-b3b0-731245ec6c26)",
      "fileSize": 836,
      "fileTypeId": "PacBio.FileTypes.JsonReport",
      "importedAt": "2023-09-28T08:22:03.998Z",
      "isActive": true,
      "modifiedAt": "2023-09-28T08:21:56.448Z",
      "name": "Adapter Report",
      "path": "<path> /adapter.report.json",
      "sourceId": "import_dataset.report_adapters",
      "uuid": "d340b519-6fdc-401a-b3b0-731245ec6c26"
    },
    "reportTypeId": "import_dataset.report_adapters"
  },
  {
    "dataStoreFile": {
      "createdAt": "2023-09-28T08:21:56.451Z",
      "description": "PacBio Report loading_xml_report (a33dfa8c-e591-453d-a9e2-d4f6aab379b6)",
      "fileSize": 4106,
      "fileTypeId": "PacBio.FileTypes.JsonReport",
      "importedAt": "2023-09-28T08:22:03.998Z",
      "isActive": true,
      "modifiedAt": "2023-09-28 <path> /loading.report.json",
      "sourceId": "import_dataset.report_loading",
      "uuid": "a33dfa8c-e591-453d-a9e2-d4f6aab379b6"
    },
    "reportTypeId": "import_dataset.report_loading"
  },
  {
    "dataStoreFile": {
      "createdAt": "2023-09-28T08:21:56.462Z",
      "description": "PacBio Report raw_data_report (924e472c-caa0-4d3d-bf10-df09b850915b)",
      "fileSize": 3390,
      "fileTypeId": "PacBio.FileTypes.JsonReport",
      "importedAt": "2023-09-28T08:22:03.999Z",
      "isActive": true,
      "modifiedAt": "2023-09-28T08:21:56.462Z",
      "name": "Raw Data Report",
      "path": "<path>/raw_data.report.json",
      "sourceId": "import_dataset.report_raw_data",
      "uuid": "924e472c-caa0-4d3d-bf10-df09b850915b"
    },
    "reportTypeId": "import_dataset.report_raw_data"
  },
  {
    "dataStoreFile": {
      "createdAt": "2023-09-28T08:21:56.619Z",
      "description": "PacBio Report barcode (63facc29-49f2-4df8-84df-75711813efac)",

```

```

    "fileSize": 33184,
    "fileTypeId": "PacBio.FileTypes.JsonReport",
    "importedAt": "2023-09-28T08:22:04.000Z",
    "isActive": true,
    "modifiedAt": "2023-09-28T08:21:56.619Z",
    "name": "Report barcode",
    "path": "<path> /barcodes.report.json",
    "sourceId": "import_dataset.report_barcode",
    "uuid": "63facc29-49f2-4df8-84df-75711813efac"
  },
  "reportTypeId": "import_dataset.report_barcode"
},
{
  "dataStoreFile": {
    "createdAt": "2023-09-28T08:21:56.600Z",
    "description": "PacBio Report control (2155f2a7-dfa9-4cc2-8e14-b7677406a8ae)",
    "fileSize": 2155,
    "fileTypeId": "PacBio.FileTypes.JsonReport",
    "importedAt": "2023-09-28T08:22:03.999Z",
    "isActive": true,
    "modifiedAt": "2023-09-28T08:21:56.600Z",
    "name": "Control Report",
    "path": "<path>/control.report.json",
    "sourceId": "import_dataset.report_control",
    "uuid": "2155f2a7-dfa9-4cc2-8e14-b7677406a8ae"
  },
  "reportTypeId": "import_dataset.report_control"
},
{
  "dataStoreFile": {
    "createdAt": "2023-10-25T08:23:00.230Z",
    "description": "PacBio Report ccs2 (e4363794-b72b-40bf-b881-180bf733416a)",
    "fileSize": 9777,
    "fileTypeId": "PacBio.FileTypes.JsonReport",
    "importedAt": "2023-10-25T08:23:16.895Z",
    "isActive": true,
    "modifiedAt": "2023-10-25T08:23:00.230Z",
    "name": "CCS Analysis Report",
    "path": "<path>/reports/ccs.report.json",
    "sourceId": "import_dataset.report_ccs2",
    "uuid": "e4363794-b72b-40bf-b881-180bf733416a"
  },
  "reportTypeId": "import_dataset.report_ccs2"
}
]

```

How to find the contents of specific Run QC reports

You can obtain the contents of the following Run QC reports:

- **Detect CpG Methylation Report**
- **Barcode Preview Report**
- **Adapter Report**
- **Loading Report**
- **Raw Data Report**
- **Report barcode**
- **Control Report**
- **CCS Analysis Report**

Following is an example of how to do this using Python:

```
reports = client.get_run_reports(run_id)
```

For `r` in reports:

```
    report_uuid = r["dataStoreFile"]["uuid"]
    report = client.load_datastore_report_file(report_uuid)
```

How to query instrument status, such as Running or Complete

You can obtain the instrument state of an instrument using the instrument's serial number.

The possible instrument states are:

Starting, WarmUp, SelfTest, Ready, Running, ShuttingDown, Problem. (Ready indicates that the prior run was complete and the instrument is ready to begin sequencing again.)

Following is an example of how to do this using Python:

```
instrument_state = client.get_instrument_state(serial)
```

To obtain the instrument states of **all** your instruments:

```
instrument_states = client.get_instrument_states()
```

How to query a job

You can retrieve a job of any type by integer ID or UUID.

Following is an example of how to do this using Python:

```
job = client.get_job(job_id)
```

Following is an example of what this returns:

```
{
  "subJobTypeId": "cromwell.workflows.pb_segment_reads",
  "name": "Read Segmentation with multiple cells",
  "updatedAt": "2023-09-28T08:28:49.099Z",
  "workflow": "{}",
  "path": "<path>",
  "state": "SUCCESSFUL",
```

```

"tags": "testkit",
"uuid": "d6e7b447-4af2-4fb2-9b8b-4188a930aa28",
"externalJobId": "d698e9f9-c521-449d-961d-b7946931ea5f",
"jobStartedAt": "2023-09-28T08:24:36.986Z",
"applicationName": "Read Segmentation",
"projectId": 1,
"childJobsCount": 0,
"jobCompletedAt": "2023-09-28T08:28:49.099Z",
"jobTypeId": "analysis",
"id": 131,
"smrtlinkVersion": "1.1.1.11111",
"comment": "Description for job Run Analysis Application",
"isNested": false,
"createdAt": "2023-09-28T08:24:35.732Z",
"isActive": true,
"createdBy": "bsmith",
"createdByEmail": "bsmith@comp[any.com]",
"isMultiJob": false,
"jsonSettings": "{\"name\": \"Read Segmentation with multiple cells\"...}",
"jobUpdatedAt": "2023-09-28T08:28:49.099Z"
}

```

How to start a SMRT Link job for a specific workflow

Following is an example of how to do this using Python:

```
job = client.create_analysis_job(options)
```

You can submit a SMRT Analysis job to be run as soon as possible. (This requires that **all** input datasets have already been imported.)

Options:

Job options schema:

- `pipelineId`: String such as 'cromwell.workflows.pb_align_ccs'.
- `name`: Job name string.
- `entryPoints`: List of dataset entry points.
- `taskOptions`: List of workflow task options.
- `projectId`: int or null.
- `presetId`: string or null.

Entry point model:

- `entryId`: Pre-set identifier, can be any of `eid_ccs`, `eid_barcode`, `eid_ref_dataset`, `eid_barcode_2`, or `eid_subread`.
- `fileTypeId`: Dataset MetaType, from the top-level XML tag.
- `datasetId`: Dataset UniqueID (UUID).

Task/workflow option model:

- `optionId`: String ID such as `mapping_min_length`.
- `value`: string, float, int, bool, or occasionally null.
- `optionTypeId`: Type of value field.

How to combine a sample split across multiple cells

Following is an example of how to do this using Python:

```
DS_TYPE = "PacBio.DataSet.ConsensusReadSet"
datasets = client.get_consensusreadsets(bioSampleName="MySample1234")
job = client.create_merge_datasets_job([d["uuid"] for d in datasets])
job = client.poll_for_successful_job(job["id"])
datastore = client.get_job_datastore(job["id"])
merged_datasets = [f for f in datastore if f["fileTypeId"] == DS_TYPE]
merged_datasets = client.get_consensusreadsets(jobId=job["id"])
```

How to poll every 10 minutes until a collection is complete, then launch a HiFi Mapping job using the official PacBio hg38 reference, and poll until it completes successfully

Following is an example of how to do this using Python:

```
import time
collection = client.get_run_collection(run_id, collection_id)
while True:
    dataset = client.get_dataset_search(collection["ccsId"])
    if dataset:
        break
    else:
        time.sleep(600)
job = client.create_analysis_job({
    "name": "My Mapping Job",
    "pipelineId": "cromwell.workflows.pb_align_ccs",
    "entryPoints": [
        {
            "entryId": "my_id",
            "datasetId": collection["ccsId"],
            "fileTypeId": "PacBio.DataSet.ConsensusReadSet"
        },
        {
            "entryId": "my_id_ref_dataset",
            "datasetId": "ab1234cd-5efg-7h99-0000-0a1209170a1b",
            "fileTypeId": "PacBio.DataSet.ReferenceSet"
        }
    ],
    "taskOptions": []
})
job = client.poll_for_successful_job(job["id"])
```

How to find information about the Run QC reports associated with an analysis job

You can obtain **all** collection-level reports associated with a run. (**Note:** This was Introduced in SMRT Link 13.0.)

Following is an example of how to do this using Python:

```
entry_points = client.get_job_entry_points(job_id)
movie_names = set([])
for entry_point in entry_points:
    if entry_point["datasetType"] == "PacBio.DataSet.ConsensusReadSet":
        dataset = client.get_consensusreadset(entry_point["datasetUUID"])
        movie_names.append(dataset["metadataContextId"])
qc_reports = []
for movie_name in movie_names:
    runs = client.get_runs(movieName=movie_name)
    if len(runs) == 1:
        collections = client.get_run_collection(runs[0]["unique_id"])
        for collection in collections:
            if collection["context"] == movie_name:
                reports = client.get_collection_reports(run_id,
                    collection["unique_id"])
                qc_reports.append(reports)
```

How to export demultiplexed dataset metrics as a CSV file

Following is an example of how to do this using Python. **Note:** This **also** requires the `pbcommand` library distributed with SMRT Link, but it is also possible to extract these data data from the raw report JSON.

```
from pbcommand.pb_io import load_report_from_json
from pbcommand.services.smrtlink_client import SmrtLinkClient

def export_barcode_report_csv(server, user, password, dataset_uuid, output_file):
    client = SmrtLinkClient.connect (server, user, password)
    reports = client.get_consensusreadset_reports(dataset_uuid)
    for report_file in reports:
        if "barcode" in report_file["reportTypeId"]:
            rpt_uuid = report_file["dataStoreFile"]["uuid"]
            report_json = client.load_datastore_report_file(rpt_uuid)
            report = load_report_from_json(report_json)
            report.tables[0].to_csv(output_file)
            break
    else:
        raise RuntimeError(f"Can't find barcodes report for {dataset_uuid}")
```