

Abstract


Single Molecule, Real-Time (SMRT®) Sequencing holds promise for addressing new frontiers to understand molecular mechanisms in evolution and gain insight into adaptive strategies. With read lengths exceeding 10 kb, we are able to sequence high-quality, closed microbial genomes with associated plasmids, and investigate large genome complexities, such as long, highly repetitive, low-complexity regions and multiple tandem-duplication events. Improved genome quality, observed at 99.9999% (QV60) consensus accuracy, and significant reduction of gap regions in reference genomes (up to and beyond 50%) allow researchers to better understand coding sequences with high confidence, investigate potential regulatory mechanisms in noncoding regions, and make inferences about evolutionary strategies that are otherwise missed by the coverage biases associated with short-read sequencing technologies.

Additional benefits afforded by SMRT Sequencing include the simultaneous capability to detect epigenomic modifications and obtain full-length cDNA transcripts that obsolete the need for assembly. With direct sequencing of DNA in real-time, this has resulted in the identification of numerous base modifications and motifs, which genome-wide profiles have linked to specific methyltransferase activities. Our new offering, the Iso-Seq™ Application, allows for the accurate differentiation between transcript isoforms that are difficult to resolve with short-read technologies. PacBio® reads easily span transcripts such that both 5'/3' primers for cDNA library generation and the poly-A tail are observed. As such, exon configuration and intron retention events can be analyzed without ambiguity. This technological advance is useful for characterizing transcript diversity and improving gene structure annotations in reference genomes.

We review solutions available with SMRT Sequencing, from targeted sequencing efforts to obtaining reference genomes (>100 Mb). This includes strategies for identifying microsatellites and conducting phylogenetic comparisons with targeted gene families. We highlight how to best leverage our long reads that have exceeded 20 kb in length for research investigations, as well as currently available bioinformatics strategies for analysis. Benefits for these applications are further realized with consistent use of size selection of input sample using the BluePippin™ device from Sage Science as demonstrated in our genome improvement projects. Using the latest P5-C3 chemistry on model organisms, these efforts have yielded an observed contig N50 of ~6 Mb, with the longest contig exceeding 12.5 Mb and an average base quality of QV50.


Single Molecule Real Time (SMRT®) Sequencing

PacBio® RS II Sequencing Chemistries Provide Long Read Lengths >20 Kb



Example read length distribution from a SMRT Sequencing run with 20 kb size-selected *E. coli* library using a 180-min movie. Average throughput of 350 Mb per SMRT Cell with ~50,000 reads.


Universal SMRTbell™ Template



Schematic of SMRTbell Template Sequencing


Top 10 Applications for Evolutionary Biology

(1) PacBio-only de novo Genome Assemblies




Public datasets, SMRT Analysis and compatible third party software are available from PacBio DevNet: <http://pacbiodevnet.com/>

(2) Hierarchical Genome Assembly Process (HGAP)




A bioinformatics solution to make best use of all PacBio reads for improved consensus accuracy and genome assemblies.


(3) Obtain Complete Microbial Genomes with Epigenetic Modifications and Associated Plasmids



(4) Conduct Whole-Genome Phylogeny




(5) Compare Genome Organization between Strains




Collaboration with A. Zeddeman, H. van der Heide, M. Bart & F. Mooi National Institute for Public Health and the Environment (RIVM), Netherlands


(6) Discover de novo Microsatellites



(7) Detect Complex Duplication & Triplication Events




(8) Direct Full-length Isoform Sequencing



(9) Survey Targeted Gene Families



(10) Whole Genome Survey of Duplication Events



Recent Publications

- PacBio Blog: *A. thaliana* genome assembly effort
- PacBio Blog: *D. melanogaster* genome assembly effort
- PacBio Blog: PAC PacBio Workshop Showcases User's Large Genome Efforts
- Chin CS, et al. "Nonhybrid, finished microbial genome assemblies from long-read SMRT sequencing data." *Nat Methods*. Jun;10(6):563-9 (2013).
- Case Study: Beyond Four Bases: Epigenetic Modifications Prove Critical to Understanding *E. Coli* Outbreak (2012)
- Cooper KK, et al. "Comparative genomics of enterohemorrhagic *Escherichia coli* O145:H28 demonstrates a common evolutionary lineage with *Escherichia coli* O156:H7 (2014)
- Grohme M.A., et al. "Microsatellite marker discovery using single molecule real-time circular consensus sequencing on the Pacific Biosciences RS" *Biotechniques* 55(5), 253-256 (2013).
- Maron, L.G. et al. (2013) A rare gene copy-number variant that contributes to maize aluminum tolerance and adaptation to acid soils. *PNAS* 110(13):5241-6
- Poster: Tseng, E. et al. "Isoform Sequencing: Unveiling the Complex Landscape of the Eukaryotic Transcriptome on the PacBio® RS II" (PAG conference 2014)
- Zhang W et al. (2013) "PacBio sequencing of gene families - A case study with wheat gluten genes." *Gene* doi:10.1016/j.gene.2013.10.009
- Rogers RL, et al. (2014) "Landscape of Standing Variation for Tandem Duplications in *Drosophila yakuba* and *Drosophila simulans*." *Mol Biol Evol* doi: 10.1093/molbev/msu124

