



## New Features

### Sequel<sup>®</sup> Instrument Control Software (ICS)

- Support for the new 3.0 versions of the Sequel Binding Kit, Sequencing Plates, and DNA Internal Control.
- Support for the new SMRT<sup>®</sup> Cell 1M v3 and SMRT<sup>®</sup> Cell 1M v3 LR products.
- ICS now accesses SMRT Link via authenticated routes.

## System Improvements

### Sequel<sup>®</sup> Instrument Control Software (ICS)

- Numerous signal processing and base caller improvements resulting in longer reads, and improved accuracy, in low-complexity sequence contexts. Improved spike-in control read filtering resulting in a ten-fold lower leak rate.
- Reduced adapter false positive calls.
- Robotics and workflow improvements.
- Increased chiller reliability.
- Improved IUI touch screen response and system status LED behavior.
- The IUI Quick Run module now selects diffusion loading by default.

## Fixed Issues

### Sequel<sup>®</sup> Instrument Control Software (ICS)

- Fixed an intermittent issue where data from a failed transfer could be deleted by system housekeeping after seven days.
- Fixed an issue where the system could become unresponsive if it attempted to send a notification email while the user's network was down.
- Fixed an issue with logging the reagent plate used in a run.

## Known Issues

### Sequel<sup>®</sup> Instrument Control Software (ICS)

- The instrument self-test intermittently becomes unresponsive. This requires a reboot of the Instrument Control Software.
- The virtual keyboard can cover input fields and selections when using the **Admin** dialog or when using the Quick Run functionality.
- When using a barcode scanner, it must be unplugged **before** restarting the instrument. If the barcode scanner is **not** unplugged, the USB controller will fail to correctly allocate the USB ports, causing a failure to properly start up the instrument.
- The on-instrument expiration time for Sequencing Plate 3.0, 8-reaction version is incorrectly set to 72 hours. The correct expiration time is 100 hours.

---

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2015 - 2018, Pacific Biosciences of California, Inc. All rights reserved. Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences products and/or third party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and the applicable license terms at <http://www.pacb.com/legal-and-trademarks/product-license-and-use-restrictions/>.

Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq and Sequel are trademarks of Pacific Biosciences. BluePippin and SageELF are trademarks of Sage Science, Inc. NGS-go and NGSengine are trademarks of GenDx. FEMTO Pulse and Fragment Analyzer are trademarks of Advanced Analytical Technologies. All other trademarks are the sole property of their respective owners.

P/N 100-715-100 Version 14 (October 2018)