

The background of the slide is a blurred image of a laboratory microplate. A pipette tip is visible at the top right, dispensing a drop of bright pink liquid into one of the wells. The wells contain varying amounts of this pink liquid. The overall scene is brightly lit, with a soft focus on the background.

PacBio

Technical overview – Kinnex library preparation for full-length 16S rRNA gene sequencing

Sequel II and IIe systems ICS v11.0

Revio system ICS v13.0+

SMRT Link v13.0+


PN 103-344-800 Rev 01 | March 2024

Technical overview – Kinnex library preparation for full-length 16S rRNA gene sequencing

Technical Overview

1. Kinnex 16S rRNA method overview
2. Kinnex 16S rRNA library preparation workflow details
3. Kinnex 16S rRNA sequencing preparation workflow details
4. Kinnex 16S rRNA example sequencing performance data
5. Kinnex 16S rRNA data analysis workflow overview
6. Technical documentation & applications support resources

Kinnex library preparation for full-length 16S rRNA gene sequencing: Getting started



PacBio

Application note

Kinnex 16S rRNA kit for full-length 16S sequencing

Introduction

Microbes are an essential part of the ecosystem for human, plant, and animal species and the environments they live in. Microbes perform metabolic activities, produce and degrade compounds, and play a role in health, fitness, phenotype, and ecology. The human microbiome has been shown to be pivotal for human health, with dysbiosis in the gut microbiome having been linked to conditions such as inflammatory bowel disease (IBD), diabetes, cardiovascular disease, colon cancer, and neurological disease. Similarly, both soil and marine microbes play an active role in organism health. These relationships demonstrate a growing need and appreciation for more comprehensively characterizing the species within microbiomes and associating them with biological outcomes.

All bacteria have a 16S rRNA gene, making targeted 16S sequencing a reliable and cost-effective approach

for assessing the composition of metagenomic communities. This is especially true for low bacterial biomass samples where amplicon sequencing is the best approach. However, the high similarity in the 16S rRNA genes among related bacteria mean that sequencing the entirety of the 16S gene (~1.5 kb) with high accuracy is essential for characterizing at the species or strain level.

Recent comparative studies have shown that PacBio® full-length 16S sequencing outperforms other sequencing methods (Notario et al., 2023, Figure 1). The Kinnex™ 16S rRNA kit takes amplified 16S gDNA amplicons as input and outputs a sequencing-ready library that results in a 12-fold throughput increase compared to other 16S libraries. Allowing up to 1,536-plex per library, Kinnex technology enables highly accurate, cost-effective 16S sequencing for microbiome studies.

© 2023 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
102-326-601 REV01 DEC2023

PacBio

Application note – Kinnex 16S rRNA kit for full-length 16S sequencing (102-326-601)

Summary overview of application-specific library preparation and data analysis workflow recommendations.

PacBio

Preparing Kinnex™ libraries from 16S rRNA amplicons

Procedure & checklist

Overview

This procedure provides instructions for generating Kinnex libraries from full-length 16S amplicons for sequencing on PacBio® Sequel II, Sequel IIe, and Revio™ systems.

1. Amplification of full-length 16S genes (V1–V9 regions) from metagenomic samples using barcoded Forward and Reverse 16S primers
2. Concatenation of 16S amplicons to ~19 kb
3. Multiplexed sequencing on the Sequel II/IIe and Revio systems

Barcoded 16S-specific primers (12 forward and 32 reverse) can be used in different combinations allowing for the multiplexing of up to 384 samples on one SMRT™ Cell. If combined with barcoded Kinnex adapters (4-plex), a total of 1,536 samples can be sequenced.

© 2024 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
103-238-800 REV03 FEB2024

PacBio

Procedure & checklist – Preparing Kinnex libraries from 16S rRNA amplicons (103-238-800)

Technical documentation containing application-specific library preparation protocol details.

PacBio

Technical overview – Kinnex library preparation for full-length 16S rRNA gene sequencing

Sequel II and IIe systems ICS v11.0
Revio system ICS v13.0
SMRT Link v13.0

PN 103-344-800 Rev 01 | December 2023


Example sequencing performance for Kinnex 16S rRNA libraries prepared from mock microbial community genomic DNA

Revio system example data* (1,536-plex data set)

PacBio

Technical overview – Kinnex library preparation for full-length 16S rRNA gene sequencing (103-344-800)

Technical overview presentations describe sample preparation details for constructing Kinnex HiFi libraries for specific applications. Example sequencing performance data for a given application are also summarized.






Kinnex 16S rRNA method overview

Kinnex 16S rRNA method overview

Use Kinnex full-length RNA kit to perform high-accuracy, full-length isoform sequencing with PacBio long-read systems



- Official protocol to generate barcoded 16S amplicons compatible with Kinnex 16S rRNA kit
- Protocol supports up to 384-plex multiplexing
- 2-day Kinnex library preparation using **Kinnex 16S rRNA kit (103-072-100)**
- SMRT Link Run Design support for 'Kinnex 16S rRNA' application type option with auto-analysis (read segmentation only)¹
- Demultiplex 16S amplicon barcodes in SMRT Link to generate per-sample read BAM files
- Analyze per-sample BAM files using GitHub tools or other custom 16S analysis pipeline

Kinnex 16S rRNA method overview (cont.)

Full-length 16S gene amplification




Bacterial gDNA fragments isolated from metagenomic samples

Recommended gDNA input for full-length 16S PCR amplification = 1-2 ng

PCR amplification of full-length 16S genes

Can barcode up to 384 samples during 16S gene amplification step using PacBio barcoded Fwd and Rev primers¹



Sample plate layout for 96-plex PCR design using 12 different 16S Barcoded Forward Primers and 8 different 16S Barcoded Reverse Primers.

To multiplex 384 samples, use 12 barcoded Fwd primers + 32 barcoded Rev primers and set up four 96-well PCR reaction plates.

Forward (F) PCR primer (Kinnex16S_Fwd_XX¹)

Reverse (R) PCR primer (Kinnex16S_Rev_XX¹)



Dual-barcoded full-length 16S amplicon product

Contains Kinnex PCR handle 1

Equal-volume pooling & cleanup



Pooled barcoded full-length 16S amplicons for input into Kinnex library construction

Contains Kinnex PCR handle 2

Kinnex library prep, sequencing & analysis




Amplified 16S products containing a different PacBio barcode at each end (Fwd BC ≠ Rev BC)

Kinnex PCR (12 parallel PCR Rx to append Kinnex segmentation adapters)



Pooling & Kinnex array formation

Barcoded Kinnex terminal adapter (bcM0001-bcM0004)²




Can also barcode samples during Kinnex array formation using up to 4 different barcoded Kinnex terminal SMRTbell adapters²

DNA damage repair & nuclease Tx



Anneal Kinnex sequencing primer / Bind polymerase / Cleanup (ABC)



Purified polymerase-bound 12-segment Kinnex library

Long-read sequencing & data analysis

Sequel II, Sequel IIe or Revio system



Third-party metagenomic secondary analysis tools

¹ 12 different 16S barcoded Forward PCR primers + 32 different 16S barcoded Reverse PCR primers are available for 16S gene amplification step to multiplex up to 384 samples.

² Kinnex adapter barcode sequences can be downloaded from [SMRT Link](https://www.pacb.com/kinnex-data-management/) Data Management module.

Kinnex 16S rRNA library preparation procedure description

Procedure & checklist – Preparing Kinnex libraries from 16S rRNA amplicons ([103-238-800](https://www.pacb.com/support/103-238-800)) describes the workflow for generating Kinnex libraries from full-length 16S amplicons for sequencing on PacBio Sequel II, Sequel IIe, and Revio systems

Overview


This procedure provides instructions for generating Kinnex libraries from full-length 16S amplicons for sequencing on PacBio® Sequel® II, Sequel IIe, and Revio™ systems.

1. Amplification of full-length 16S genes (V1–V9 regions) from metagenomic samples using barcoded Forward and Reverse 16S primers
2. Concatenation of 16S amplicons to ~19 kb
3. Multiplexed sequencing on the Sequel II/IIe and Revio systems

Barcoded 16S-specific primers (12 forward and 32 reverse) can be used in different combinations allowing for the multiplexing of up to 384 samples on one SMRT® Cell. If combined with barcoded Kinnex adapters (4-plex), a total of 1536 samples can be sequenced.



Kinnex 16S rRNA kit
103-072-100 (12 rxn)



Kinnex 16S rRNA library template (~18 kb)
Contains 12 concatenated full-length 16S amplicon segments

Preparing Kinnex™ libraries from 16S rRNA amplicons



Procedure & checklist

Overview

This procedure provides instructions for generating Kinnex libraries from full-length 16S amplicons for sequencing on PacBio® Sequel® II, Sequel IIe, and Revio™ systems.

1. Amplification of full-length 16S genes (V1–V9 regions) from metagenomic samples using barcoded Forward and Reverse 16S primers
2. Concatenation of 16S amplicons to ~19 kb
3. Multiplexed sequencing on the Sequel II/IIe and Revio systems

Barcoded 16S-specific primers (12 forward and 32 reverse) can be used in different combinations allowing for the multiplexing of up to 384 samples on one SMRT® Cell. If combined with barcoded Kinnex adapters (4-plex), a total of 1536 samples can be sequenced.

© 2024 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
103-238-800 REV02 MAR2024



PacBio [Documentation](https://www.pacb.com/support/103-238-800) (103-238-800)

- Kinnex 16S rRNA library prep protocol uses **Kinnex 16S rRNA kit**
→ **Do not use** SMRTbell prep kit 3.0 with this protocol

Kinnex 16S rRNA kit components

Kinnex 16S rRNA kit provides full support for Kinnex library prep workflow starting with 16S amplicon DNA as input

Kinnex 16S rRNA kit bundle (103-072-100)

Includes Kinnex PCR kit, Kinnex concatenation and ancillary DNA cleanup reagents needed for incorporation of Kinnex segmentation adapters and Kinnex array formation for generating Kinnex 16S rRNA libraries.

Kinnex 16S rRNA kit components


Component	Description
1	 Kinnex PCR 12-fold kit (12 rxn) <ul style="list-style-type: none">Contains reagents for Kinnex PCR to incorporate segmentation adapters
2	 Kinnex concatenation kit (12 rxn) <ul style="list-style-type: none">Contains reagents for Kinnex array formation and SMRTbell template constructionIncludes barcoded Kinnex adapter mixes (bcM0001 – bcM0004)
3	 SMRTbell cleanup beads <ul style="list-style-type: none">For DNA cleanup
4	 Elution buffer <ul style="list-style-type: none">For DNA cleanup

Kinnex 16S rRNA library barcoding options for sample multiplexing

Kinnex 16S rRNA library preparation procedure supports up to 1,536-plex sample multiplexing

Kinnex 16S rRNA library preparation procedure supports **up to 1,536-plex** sample multiplexing through combined use of:

- 12 different 16S barcoded Forward PCR primers¹ (Kinnex16S_Fwd_01 – Kinnex16S_Fwd_12)
- 32 different 16S barcoded Reverse PCR primers¹ (Kinnex16S_Rev_13 – Kinnex16S_Rev_44)
- 4 different barcoded Kinnex terminal SMRTbell adapters (bc01 – bc04)



Kinnex 16S rRNA experimental design considerations

Kinnex 16S rRNA application use case recommendations for PacBio systems


	Sequel II and Ile systems		Revio system	
Experimental goal	Determine the microbial diversity (phylogeny and taxonomy) of bacteria in a metagenomic sample			
Sample multiplexing ¹	Up to 384 samples per SMRT Cell 8M (384-plex)		Up to 1,536 samples per Revio SMRT Cell (1536-plex)	
Expected coverage per sample ²	96-plex	260 K	96-plex	625 K
	192-plex	130 K	192-plex	313 K
	384-plex	65 K	384-plex	156 K
	768-plex	33 K	768-plex	78 K
	1,536-plex	16 K	1,536-plex	39 K
Kinnex library prep protocol	Procedure & checklist – Preparing Kinnex libraries from 16S rRNA amplicons (103-238-800)			
Metagenomic DNA input amount input into 16S gene amplification	1-2 ng of input gDNA per metagenomic sample			
16S amplicon DNA input into Kinnex library prep workflow	35 ng of purified pooled 16S amplicon DNA			
SMRT Link data analysis workflows	Read Segmentation			
Community data analysis tools	pb-16S-nf			

¹ Kinnex concatenation kit (103-071-800) can support up to 1,536-plex sample multiplexing through the combined use of 12 different 16S barcoded Forward PCR primers + 32 different 16S barcoded Reverse PCR primers and 4 different barcoded Kinnex terminal SMRTbell adapters during Kinnex 16s rRNA library construction.

² With proper full array formation and adequate sequencing, one SMRT Cell on the Sequel II, Ile, and Revio systems are expected to achieve 20–25 million and 50–60 million 16S sequences, respectively. For most 16S analysis applications, typically aim for ~30-50 K reads/sample.

Kinnex 16S rRNA experimental design considerations (cont.)

PacBio full-length 16S protocol and PacBio long-read sequencing platform recommendations for different sample types



Stool samples typically need 10–20k average reads/sample; soil samples need 30–50k

¹ Refer to *Procedure & checklist – Preparing Kinnex libraries from 16S rRNA amplicons* (103-238-800).

² Refer to *Procedure & checklist – Amplification of bacterial full-length 16S rRNA gene with barcoded primers* (101-599-700) and *Procedure & checklist – Preparing multiplexed amplicon libraries using SMRTbell prep kit 3.0* (102-359-000).



Kinnex 16S rRNA library preparation workflow details

Procedure & checklist – Preparing Kinnex libraries from 16S amplicons (103-238-800)

Procedure & checklist [103-238-800](#) describes the workflow for generating Kinnex libraries from full-length 16S amplicons for sequencing on PacBio Sequel II, Sequel IIe, and Revio systems¹

Procedure & checklist contents

1. **Barcoded 16S gene-specific forward and reverse primer** oligo synthesis and storage recommendations.
2. **Bacterial input genomic DNA QC recommendations** and **general best practices** for reagent & sample handling.
3. Enzymatic workflow steps for **PCR amplification of 16S gene with barcoded primers**.
4. Enzymatic workflow steps for **construction of 12-segment Kinnex arrays** from amplified 16S DNA products.
5. Enzymatic workflow steps for **DNA damage repair & nuclease treatment** of Kinnex libraries.
6. Workflow steps for **final cleanup of Kinnex SMRTbell libraries** using SMRTbell cleanup beads.

Preparing Kinnex™ libraries from 16S rRNA amplicons PacBio

Procedure & checklist

Overview

This procedure provides instructions for generating Kinnex libraries from full-length 16S amplicons for sequencing on PacBio® Sequel® II, Sequel IIe, and Revio™ systems.

1. Amplification of full-length 16S genes (V1–V9 regions) from metagenomic samples using barcoded Forward and Reverse 16S primers
2. Concatenation of 16S amplicons to ~19 kb
3. Multiplexed sequencing on the Sequel II/IIe and Revio systems

Barcoded 16S-specific primers (12 forward and 32 reverse) can be used in different combinations allowing for the multiplexing of up to 384 samples on one SMRT® Cell. If combined with barcoded Kinnex adapters (4-plex), a total of 1536 samples can be sequenced.


© 2024 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
103-238-800 REV02 MAR2024

PacBio

PacBio [Documentation](#) ([103-238-800](#))

Kinnex 16S rRNA library construction workflow overview

Procedure & checklist – Preparing Kinnex libraries using the Kinnex 16S rRNA kit (103-238-800)



General best practices recommendations for preparing Kinnex 16S rRNA libraries

Preparation of barcoded 16S gene-specific forward and reverse primers

- We recommend resuspending stock oligos with a target concentration of 100 μM in 10 mM Tris-HCl pH 8.0–8.5 (elution buffer) or low TE (10 mM Tris-HCl with 0.1 mM EDTA)
- To prepare oligo working solutions, dilute each primer individually to 2.5 μM in 10 mM
- Tris-HCl pH 8.0–8.5 (elution buffer) or low TE. For example, add 5 μL of 100 μM primer stock to 195 μL of 10 mM Tris-HCl pH 8.0–8.5 buffer. This volume of diluted oligo is sufficient for running more than 50 PCR reactions
- Always mix primer stocks well before preparing dilutions. Prior to use, verify that the concentration of each diluted oligo solution is 2.5 μM by directly measuring the OD260 value using a Nanodrop system
- Aliquot the diluted oligos in 96-well plates in the format shown in the plate maps below

Plate map for preparing barcoded 16S gene-specific forward primers

	1	2	3	4	5	6	7	8	9	10	11	12
A	Fwd_01	Fwd_02	Fwd_03	Fwd_04	Fwd_05	Fwd_06	Fwd_07	Fwd_08	Fwd_09	Fwd_10	Fwd_11	Fwd_12
B												
C	Fwd_01	Fwd_02	Fwd_03	Fwd_04	Fwd_05	Fwd_06	Fwd_07	Fwd_08	Fwd_09	Fwd_10	Fwd_11	Fwd_12
D												
E	Fwd_01	Fwd_02	Fwd_03	Fwd_04	Fwd_05	Fwd_06	Fwd_07	Fwd_08	Fwd_09	Fwd_10	Fwd_11	Fwd_12
F												
G	Fwd_01	Fwd_02	Fwd_03	Fwd_04	Fwd_05	Fwd_06	Fwd_07	Fwd_08	Fwd_09	Fwd_10	Fwd_11	Fwd_12
H												

Plate map for preparing barcoded 16S gene-specific reverse primers

	1	2	3	4	5	6	7	8	9	10	11	12
A	Rev_13		Rev_21		Rev_29		Rev_37					
B	Rev_14		Rev_22		Rev_30		Rev_38					
C	Rev_15		Rev_23		Rev_31		Rev_39					
D	Rev_16		Rev_24		Rev_32		Rev_40					
E	Rev_17		Rev_25		Rev_33		Rev_41					
F	Rev_18		Rev_26		Rev_34		Rev_42					
G	Rev_19		Rev_27		Rev_35		Rev_43					
H	Rev_20		Rev_28		Rev_36		Rev_44					

General best practices recommendations for preparing Kinnex 16S rRNA libraries (cont.)

DNA input and reagent & sample handling

Required DNA input amount

- Recommended input bacterial genomic DNA amount per sample is **1–2 ng**
- Expected 16S amplicon size is approximately **1,500 bases** and typical 16S PCR product yields are **50–300 ng** per sample

DNA QC & handling recommendations for bacterial gDNA isolated from metagenomic samples


- For best results, characterize bacterial gDNA samples thoroughly and **normalize gDNA concentrations** before use
- Bring gDNA samples to room temperature and mix well by pipetting to ensure sample homogeneity, then measure gDNA concentration using Qubit dsDNA HS assay reagents
- Assess sample purity using a Nanodrop system → OD260/280 should be between **1.8 and 2.0** for purified double-stranded DNA
- To ensure pipetting accuracy, plan to **deliver 1–2 ng of gDNA to each individual PCR reaction in a constant 5 µL volume**
 - Normalize sample gDNA concentration to 0.2–0.4 ng/µL in 10 mM Tris-HCl pH 8.0–8.5 (elution buffer) prior to setting up PCR reactions.
 - Note: Nuclease-free water and Elution buffer (EB) from PacBio can be used in place of 10mM Tris-HCl pH 8.0–8.5 for gDNA normalization.
- Based on prior PacBio experience, QIAGEN Powerfecal Pro kit extracts DNA of sufficient quality for this workflow.

Reagent and sample handling

- **Thaw PCR Ready Mix on ice** and mix well before use
- **Note: All PCR reactions described in this procedure must be set up and kept on ice** until read to load onto thermal cycler instrument
 - High proofreading activity of the enzyme in the PCR Ready Mix will rapidly degrade primers at room temperature

16S full-length sequence capture

Perform PCR amplification of 16S gene with barcoded forward and reverse primers



1. PCR amplification of 16S gene with barcoded primers

✓ Step Instructions

Thaw the PCR Ready Mix, briefly vortex to mix, and place on ice. Note that all PCR reactions must be set up on ice; the high proofreading activity of the enzyme will result in rapid primer degradation at room temperature.

Thaw plates containing the diluted forward and reverse primers. Briefly spin the plate to ensure that the entire volume is at the bottom of each well.

1.1

16S PCR Master Mix 1 components	1 sample	N	For 96-plex	For 192-plex	For 384-plex
PCR-grade Water	1.5 µL	1.5 x N x 1.1	158 µL	317 µL	634 µL
2X KAPA HiFi HotStart ReadyMix	12.5 µL	12.5 x N x 1.1	1320 µL	2640 µL	5280 µL
Total	14 µL	14 x N x 1.1	1478 µL	2957 µL	5914 µL

Transfer 14 µL of the prepared 16S PCR Master Mix 1 into a 96-well PCR plate for each 96-plex. For a 96-plex experiment design, use one 96-well plate. For a 384-plex experiment design, use four 96-well plates. Add 5 µL (1ng) of each diluted gDNA sample to each well containing 16S PCR Master Mix 1 on ice.

1.2

The figures below illustrate an example plate layout for setting up a 96-plex PCR design using twelve different 16S Barcoded Forward Primers and eight different 16S Barcoded Reverse Primers. Please refer to Appendix - 384 barcodes layout for all of the plate designs.


1.3

- Prepare PCR Master Mix¹ on ice for processing up to 384 metagenomic DNA samples (include 10% overage)
- Ensure all reagents are thawed and mixed prior to use

- **IMPORTANT:** All PCR reactions described in this procedure must be set up and kept on ice until read to load onto thermal cycler instrument
 - High proofreading activity of the enzyme in the PCR Ready Mix will rapidly degrade primers at room temperature

16S full-length sequence capture (cont.)

Perform PCR amplification of 16S gene with barcoded forward and reverse primers



1. PCR amplification of 16S gene with barcoded primers

1.4 On ice, add 3 μL of the Barcoded Forward Primers (2.5 μM) to wells containing 19 μL of gDNA and 16S PCR Master Mix followed by 3 μL of the Barcoded Reverse Primers (2.5 μM). The final concentration of the barcoded forward and reverse primers in each well is 0.3 μM . The final reaction volume in each well is 25 μL . Mix well by pipetting. Seal the plates to prevent evaporation during PCR. Briefly spin the plate in a refrigerated centrifuge (4°C) to ensure that the entire sample volume is at the bottom of each well.

1.5 Set up a thermal cycler with the program shown below. Set the lid temperature to 105°C and pre-heat the thermal cycler until the lid temperature reaches 105°C and before adding the 96-well PCR plate. Keep the 96-well PCR plates on ice until the lid is pre-heated. The duration of PCR is around 1 hour.

Step	Temperature	Duration	Cycle
Initial Denaturation	95 °C	3 min	1
Denaturation	98 °C	20 s	20
Annealing	57 °C	30 s	
Extension	72 °C	75 s	
Final Extension	72 °C	5 min	1
Hold	4 °C	Hold	

1.6 Spot-check amplification results by directly loading 1 μL of one or more PCR products onto an Agilent Bioanalyzer Chip. The expected target amplicon size is ~1500 bp, and the amount of amplicon material generated from each sample should be comparable as assessed by analyzing the relative intensity of the ~1500 bp PCR product. (Figure 1)

1.7 Proceed to pooling and SMRTbell cleanup in the next step.


SAFE STOPPING POINT -- Store at 4°C or -20°C for long-term storage

- **IMPORTANT:** Pre-heat thermal cycler until lid temperature reaches 105°C before loading the 96-well PCR plate.
 - Keep the 96-well PCR plates on ice until the lid is pre-heated.

- Expected 16S amplicon size is approximately 1,500 bases and typical 16S PCR product yields are 50–300 ng per sample

16S PCR amplicon pooling & SMRTbell cleanup

Pool barcoded 16S PCR amplicons and perform cleanup using 1.1X SMRTbell cleanup beads



• Minimum amount of pooled 16S amplicon DNA needed to proceed with Kinnex PCR = 35 ng

2.1 Pooling of barcoded 16S PCR amplicons

Step	Instructions
2.1.1	<p>If PCR products are of the expected size and comparable quantity, pool equal volumes of each PCR reaction in a clean DNA LoBind microcentrifuge tube according to the recommendations below:</p> <ul style="list-style-type: none"> For a 96-plex experiment design, we recommend pooling 10 μL from each PCR reaction. For a 192-plex or higher-plex experiment design, we recommend pooling 5 μL from each PCR reaction. <p>Typical total yield from each 25 μL PCR reaction is ~50-300 ng. If doing less than 96-plex, pool 20 μL from each PCR reaction into subsequent steps, but sure there is at least 35ng into the Kinnex PCR step.</p> <p>Store unused PCR reactions at -20°C for future use if desired.</p>
2.1.2	Proceed to SMRTbell cleanup in the next step.

- For 96-plex: Pool 10 μ L from each PCR reaction
- For 192-plex or higher: Pool 5 μ L from each PCR reaction
- For <96-plex: Pool 20 μ L from each PCR reaction into subsequent steps
- Typical total yield from each 25 μ L PCR reaction is ~50-300 ng
- Ensure there is at least 35 ng of pooled amplicon DNA to proceed to the Kinnex PCR step


2.2 Cleanup of pooled 16S PCR amplicon using 1.1X SMRTbell cleanup beads

Step	Instructions
2.2.1	Add 1.1X v/v (volume over volume) of resuspended, room-temperature SMRTbell cleanup beads to the tube of pooled 16S amplicon. Note: Please use a 5 mL LoBind tube if the volume is more than 2mL.
2.2.2	Pipette-mix the beads until evenly distributed.
2.2.3	Quick-spin the tube in a microcentrifuge to collect liquid.
2.2.14	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer the supernatant to a single new 1.5mL LoBind tube or tube strip. Discard the old tube with beads.
2.2.15	Measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit using 1 μ L aliquot from the LoBind tube. Typical total yield from each 25 μ L PCR reaction is ~50-300 ng.



SAFE STOPPING POINT -- Store at 4°C or -20°C for long-term storage

- Perform 1.1X SMRTbell bead cleanup at room temp.





Example Bioanalyzer DNA sizing QC analysis results for pooled 16S amplicon DNA samples generated from mock microbial community genomic DNA (ATCC MSA-1003 20 Strain Staggered Mix).

Kinnex PCR

In this step, incorporate programmable Kinnex segmentation adapter sequences into amplified DNA products

- 1 16S full-length sequence capture
- ↓ ○ Safe stop
- 2 Pooling & SMRTbell cleanup
- ↓ ○ Safe stop
- 3 **Kinnex PCR & SMRTbell cleanup**
- ↓ ○ Safe stop
- 4 Kinnex array formation & DNA repair
- ↓
- 5 SMRTbell cleanup
- ↓ ○ Safe stop
- 6 Nuclease treatment
- ↓
- 7 Final cleanup with SMRTbell cleanup beads

Set up 12 parallel Kinnex PCR reactions per 16S sample with premixed Kinnex primers to generate amplified DNA products containing programmable sequences at both ends




Example: Amplified DNA products from Kinnex PCR 3 contain flanking segmentation adapter sequences that are **complementary** to the ends of DNA products from Kinnex PCR 2 & Kinnex PCR 4



Kinnex PCR (cont.)

Procedural notes



3.1 Prepare Kinnex primers premix

Step	Instructions		
3.1.1	Thaw the following components. The entire volume of primers can be transferred to an 8-strip tube for ease of use with a multi-channel pipette.		
12X concatenation	Tube color	P/N	
1	Kinnex primer mix A	Orange	103-107-800
2	Kinnex primer mix B	Orange	103-107-900
3	Kinnex primer mix C	Orange	103-108-000
4	Kinnex primer mix D	Orange	103-108-100
5	Kinnex primer mix E	Orange	103-108-200
6	Kinnex primer mix F	Orange	103-108-300
7	Kinnex primer mix G	Orange	103-108-400
8	Kinnex primer mix H	Orange	103-153-000
9	Kinnex primer mix I	Orange	103-153-100
10	Kinnex primer mix J	Orange	103-153-200
11	Kinnex primer mix K	Orange	103-153-300
12	Kinnex primer mix LQ	Orange	103-144-000


3.1.2 Briefly vortex to mix, and quick-spin to collect liquid. Place the primer mixes on ice and proceed to the preparation of the Kinnex PCR master mix.

- Can transfer entire volume of primers to PCR tubes for ease of use with multi-channel pipettes (12 primer mix tubes)



Kinnex PCR (cont.)

Procedural notes



3.2 Kinnex PCR

Step	Instructions																										
3.2.1	<p>Thaw the following components on ice, briefly vortex to mix, and quick-spin to collect liquid. Place reagents on ice. Add the below components to a LoBind tube, pulse vortex to mix and quick-spin. Place master mix on ice.</p> <table border="1"><thead><tr><th>Master mix components</th><th>Volume for 12X concatenation*</th></tr></thead><tbody><tr><td>PCR grade water</td><td>132-X μL</td></tr><tr><td>Kinnex PCR Mix (103-107-700)</td><td>165 μL</td></tr><tr><td>35 ng of purified amplicons from Step 2.2.14*</td><td>X μL</td></tr><tr><td>Total volume</td><td>297 μL</td></tr></tbody></table> <p>X= 35 (ng)/ purified pooled DNA concentration from step 2.2.14</p> <p>*10% overage included.</p>	Master mix components	Volume for 12X concatenation*	PCR grade water	132-X μ L	Kinnex PCR Mix (103-107-700)	165 μ L	35 ng of purified amplicons from Step 2.2.14 *	X μ L	Total volume	297 μ L																
Master mix components	Volume for 12X concatenation*																										
PCR grade water	132-X μ L																										
Kinnex PCR Mix (103-107-700)	165 μ L																										
35 ng of purified amplicons from Step 2.2.14 *	X μ L																										
Total volume	297 μ L																										
3.2.2	Distribute 22.5 μ L of Master Mix 2 into each 12 PCR tubes (for 12X concatenation) on ice.																										
3.2.3	Add 2.5 μ L of Kinnex primers premix into each of 12 PCR tubes of Step 3.2.2 on ice.																										
3.2.4	<p>Set up the thermal cycler as shown below with lid temperature set to 105°C. Keep samples on ice and do not add samples to thermal cycler until the lid has reached 105°C. The duration of PCR is approximately 35 minutes.</p> <table border="1"><thead><tr><th>Step</th><th>Temperature</th><th>Duration</th><th>Cycle</th></tr></thead><tbody><tr><td>Initial Denaturation</td><td>98°C</td><td>3 min</td><td>1</td></tr><tr><td>Denaturation</td><td>98°C</td><td>20 s</td><td rowspan="3">9</td></tr><tr><td>Annealing</td><td>68°C</td><td>30 s</td></tr><tr><td>Extension</td><td>72°C</td><td>90 s</td></tr><tr><td>Final Extension</td><td>72°C</td><td>5min</td><td>1</td></tr><tr><td>Hold</td><td>4°C</td><td>Hold</td><td></td></tr></tbody></table>	Step	Temperature	Duration	Cycle	Initial Denaturation	98°C	3 min	1	Denaturation	98°C	20 s	9	Annealing	68°C	30 s	Extension	72°C	90 s	Final Extension	72°C	5min	1	Hold	4°C	Hold	
Step	Temperature	Duration	Cycle																								
Initial Denaturation	98°C	3 min	1																								
Denaturation	98°C	20 s	9																								
Annealing	68°C	30 s																									
Extension	72°C	90 s																									
Final Extension	72°C	5min	1																								
Hold	4°C	Hold																									

SAFE STOPPING POINT -- Store at 4°C or -20°C for long-term storage


- Set up Kinnex PCR reactions **ON ICE**
- PCR polymerase 3'→5' exonuclease activity negatively impacts amplification yield if prepared at room temp.

- **Critical step! Correct setup of all 12 Kinnex PCR reactions is required** – any missing/incorrect MAS primer pairs will result in no/low SMRTbell yield

- Set up on ice and add PCR reaction to thermal cycler after lid has preheated to 105°C to avoid digestion of primers by polymerase exonuclease activity

Kinnex PCR (cont.)

Procedural notes



3.3 Pooling of 12 Kinnex PCR products and 1.1x SMRTbell cleanup

✓ Step	Instructions
3.3.1	Add 23 μL from each of the 12 PCR reactions into a 1.5 mL tube for a total volume of 276 μL .
3.3.2	Add 1.1X v/v (volume over volume, 304 μL) of resuspended, room-temperature SMRTbell cleanup beads to the tube of pooled Kinnex PCR amplicon.
3.3.3	Pipette-mix the beads or invert the tube until evenly distributed.
3.3.4	Quick-spin the tube in a microcentrifuge to collect liquid.
3.3.5	Incubate at room temperature for 10 minutes to allow DNA to bind the beads.
3.3.6	Place the tube in a magnetic separation rack until the beads separate fully from the solution.
3.3.7	Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
3.3.8	Slowly dispense 1000 μL , or enough to cover the beads, of freshly prepared 80% ethanol into the tube. After 30 seconds, remove the 80% ethanol and discard.
3.3.9	Repeat the previous step.
	Remove residual 80% ethanol:
	<ul style="list-style-type: none">Remove the tube from the magnetic separation rack.Quick-spin the tube in a microcentrifuge.Place the tube back in a magnetic separation rack until the beads separate fully from the solution.Pipette off residual 80% ethanol and discard.
3.3.10	
3.3.11	Remove the tube from the magnetic rack. Immediately add 40 μL of Elution buffer to the tube and resuspend the beads by pipetting 10 times or until evenly distributed.
3.3.12	Quick-spin the tube in a microcentrifuge to collect liquid.
3.3.13	Incubate at room temperature for 5 minutes to elute DNA.
3.3.14	Place the tube in a magnetic separation rack until the beads separate fully from the solution.
3.3.15	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer the supernatant to a PCR tube strip. Discard the old tube with beads.
3.3.16	Make a 1:10 dilution of the sample and measure the DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit. Typical yield is 5–9 μg .

- Pool exactly 23 μL from each Kinnex PCR reaction in a clean 1.5 mL DNA LoBind tube¹ for a total combined volume of 184 μL

- Add exactly 304 μL of SMRTbell cleanup beads (1.1X)
- Kinnex PCR mix significantly increases stringency of SMRTbell clean up beads, so accurate pipetting is critical


- Perform DNA concentration QC using Qubit dsDNA HS kit
- Typical yield of Kinnex PCR products is 5–9 μg
- Verify there is sufficient yield of Kinnex PCR products (min. 4 μg) to proceed to Kinnex array formation step

Kinnex array formation

In this step, assemble 16S amplicons (“segments”) containing programmable ends into a linear array

- 1 16S full-length sequence capture
- 2 Pooling & SMRTbell cleanup
- 3 Kinnex PCR & SMRTbell cleanup
- 4 **Kinnex array formation & DNA repair**
- 5 SMRTbell cleanup
- 6 Nuclease treatment
- 7 Final cleanup with SMRTbell cleanup beads

Treat pooled Kinnex PCR products with Kinnex enzyme to create single-stranded overhangs to enable subsequent directional assembly of 16S amplicons into a linear, ordered array




Kinnex library prep workflow supports incorporation of PacBio barcodes at the array formation step to enable up to **4-plex sample multiplexing at the SMRTbell library level**

Kinnex array formation (cont.)

Kinnex terminal adapters incorporate barcode sequences to enable up to 4-plex sample multiplexing at the library level


- Kinnex adapters contain **barcode sequences¹** to enable (optional) sample multiplexing at the SMRTbell library level (**up to 4-plex**)
 - Forward and reverse Kinnex adapter pairs are pre-mixed in Kinnex concatenation kits
 - Kinnex concatenation kits contain a total of **4 barcoded Kinnex adapter mixes (bc01-bc04)** to enable multiplexing of up to 4 samples per SMRT Cell



Kinnex array formation (cont.)

Kinnex terminal adapters use a new design that enables improved SMRT sequencing performance

- Kinnex adapters enable:
 - Longer polymerase read length → Improved HiFi conversion rate (HiFi reads/Total *P1* reads)
 - Improved *P1* loading efficiency




- New Kinnex adapter design requires a **different sequencing primer (Kinnex sequencing primer 103-179-000)**

The diagram compares two sequencing setups. On the left, a Sequel II binding kit 3.2 and cleanup beads (102-333-300) is shown next to a Sequel II sequencer. Two vials are shown: a blue one labeled 'Sequencing primer 3.2' and a pink one labeled 'Kinnex sequencing primer (103-179-000)'. On the right, a Revio polymerase kit² (102-817-600) is shown next to a Revio sequencer. Two vials are shown: a blue one labeled 'Standard sequencing primer' and a pink one labeled 'Kinnex sequencing primer (103-179-000)'.

Kinnex array formation & DNA damage repair (cont.)

Procedural notes



4. Kinnex array formation

Step	Instructions															
4.1.1	In a 0.2 mL PCR tube, add 5 µg of sample from Step 3.3.15 , in 39 µL of volume (128 ng/µL). Proceed with Elution buffer going into this step if the sample is too concentrated. Add 2 µL of Kinnex adapter barcode 01–04 mix (select a single barcode per sample). Note: if not barcoding, select any Kinnex adapter barcode for use.															
	<table border="1"> <thead> <tr> <th>Tube color</th> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>Blue</td> <td>Kinnex adapter</td> <td>2.0 µL</td> </tr> </tbody> </table>	Tube color	Component	Volume	Blue	Kinnex adapter	2.0 µL									
Tube color	Component	Volume														
Blue	Kinnex adapter	2.0 µL														
4.1.2	Add the following components in the listed order. If processing multiple samples, make a master mix with 10% overage. Pipette to mix.															
	<table border="1"> <thead> <tr> <th>Tube color</th> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>White</td> <td>Kinnex array and repair buffer (103-110-300)</td> <td>7.0 µL</td> </tr> <tr> <td>Red</td> <td>Kinnex enzyme (103-110-400)</td> <td>4.0 µL</td> </tr> <tr> <td>Yellow</td> <td>Kinnex ligase (103-110-500)</td> <td>6.0 µL</td> </tr> <tr> <td colspan="2">Total RM1 volume</td> <td>17 µL</td> </tr> </tbody> </table>	Tube color	Component	Volume	White	Kinnex array and repair buffer (103-110-300)	7.0 µL	Red	Kinnex enzyme (103-110-400)	4.0 µL	Yellow	Kinnex ligase (103-110-500)	6.0 µL	Total RM1 volume		17 µL
Tube color	Component	Volume														
White	Kinnex array and repair buffer (103-110-300)	7.0 µL														
Red	Kinnex enzyme (103-110-400)	4.0 µL														
Yellow	Kinnex ligase (103-110-500)	6.0 µL														
Total RM1 volume		17 µL														
	Add 17 µL of master mix to the PCR tube containing sample and Kinnex barcode adapter. Pipette-mix and run the Kinnex primer digestion/ligation program with the lid set to 55°C.															
4.1.3	<table border="1"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Duration</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>45°C</td> <td>60 min</td> </tr> <tr> <td>2</td> <td>4°C</td> <td>Hold</td> </tr> </tbody> </table>	Step	Temperature	Duration	1	45°C	60 min	2	4°C	Hold						
Step	Temperature	Duration														
1	45°C	60 min														
2	4°C	Hold														
	After running the Kinnex primer digestion/ligation program, add 2 µL of DNA repair mix directly to the Kinnex primer digestion/ligation sample.															
4.1.4	<table border="1"> <thead> <tr> <th>Tube color</th> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>Green</td> <td>DNA repair mix (103-110-000)</td> <td>2 µL</td> </tr> </tbody> </table>	Tube color	Component	Volume	Green	DNA repair mix (103-110-000)	2 µL									
Tube color	Component	Volume														
Green	DNA repair mix (103-110-000)	2 µL														
4.1.5	Thoroughly mix by pipetting up and down 10 times and then quick-spin to collect all liquid.															
4.1.6	Run the DNA Damage Repair Program with the lid set to >55°C.															
	<table border="1"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Duration</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>45°C</td> <td>30 min</td> </tr> <tr> <td>2</td> <td>4°C</td> <td>Hold</td> </tr> </tbody> </table>	Step	Temperature	Duration	1	45°C	30 min	2	4°C	Hold						
Step	Temperature	Duration														
1	45°C	30 min														
2	4°C	Hold														


- Recommended input amount to proceed with Kinnex array formation is **5 µg** of Kinnex PCR amplicons (from Step 3)
- Proceeding with <3 µg is **not recommended** since lower input amounts may lead to insufficient final library yields to enable optimal sequencing results

- IMPORTANT:** If combining multiple barcoded Kinnex libraries for sequencing, make sure each library uses one of the 4 different Kinnex barcoded adapters

- Perform **DNA Damage Repair** step to repair nicked / damaged DNA sites within newly formed Kinnex array products

SMRTbell bead cleanup

Procedural notes





5. 1X SMRTbell bead cleanup

✓	Step	Instructions
	5.1	Add 60 μ L (1X v/v) of resuspended, room temperature SMRTbell cleanup beads to each sample.
	5.2	Pipette-mix the beads until evenly distributed and quick-spin in a microcentrifuge to collect liquid.
	5.3	Incubate at room temperature for 10 minutes to allow the DNA to bind the beads.
	5.4	Place the tube strip in a magnetic separation rack until the beads separate fully from the solution.
	5.5	Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
	5.6	Slowly dispense 200 μ L, or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, remove the 80% ethanol and discard.
	5.7	Repeat the previous step.
	5.8	Remove residual 80% ethanol: <ul style="list-style-type: none">Remove the tube strip from the magnetic separation rack.Quick-spin the tube strip in a microcentrifuge.Place the tube strip back in a magnetic separation rack until the beads separate fully from the solution.Pipette off residual 80% ethanol and discard.
	5.9	Remove the tube strip from the magnetic rack. Immediately add 40 μ L of Elution buffer to each tube and resuspend the beads by pipetting 10 times or until evenly distributed. Quick-spin the tube strip in a microcentrifuge to collect liquid.
	5.10	Leave at room temperature for 5 minutes to elute DNA.
	5.11	Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
		SAFE STOPPING POINT -- Store at 4°C or -20°C for long-term storage

- Perform **1.0X** SMRTbell bead cleanup at room temp.

Nuclease treatment & final SMRTbell bead cleanup


Perform nuclease treatment and final SMRTbell bead cleanup to remove incomplete SMRTbell template arrays



- After nuclease treatment, most remaining SMRTbell templates are complete (full-length) 12-segment transcript array molecules capped with Kinnex adapters

Nuclease treatment

Procedural notes



Final Kinnex library yield is typically sufficient to load **≥2 SMRT Cells**


6. Nuclease treatment

Step	Instructions												
6.1	Add the following components to a new microcentrifuge tube. Adjust the component volumes for the number of samples being prepared, plus 10% overage. Pipette-mix the master mix. For individual preps, add components directly to each sample from the previous step in the order and volume listed below. Nuclease master mix <table border="1"><thead><tr><th>Tube</th><th>Component</th><th>Volume</th></tr></thead><tbody><tr><td>Light purple</td><td>Nuclease buffer (103-110-200)</td><td>5 µL</td></tr><tr><td>Light green</td><td>Nuclease mix (103-110-100)</td><td>5 µL</td></tr><tr><td colspan="2">Total volume</td><td>10 µL</td></tr></tbody></table>	Tube	Component	Volume	Light purple	Nuclease buffer (103-110-200)	5 µL	Light green	Nuclease mix (103-110-100)	5 µL	Total volume		10 µL
Tube	Component	Volume											
Light purple	Nuclease buffer (103-110-200)	5 µL											
Light green	Nuclease mix (103-110-100)	5 µL											
Total volume		10 µL											
6.2	Add 10 µL of Nuclease Master mix to each sample. The total volume should be 50 µL.												
6.3	Thoroughly mix by pipetting up and down 10 times and then quick-spin to collect all liquid.												
6.4	Run the nuclease treatment program with the lid set to >47°C. <table border="1"><thead><tr><th>Step</th><th>Temperature</th><th>Duration</th></tr></thead><tbody><tr><td>1</td><td>37°C</td><td>15 min</td></tr><tr><td>2</td><td>4°C</td><td>Hold</td></tr></tbody></table>	Step	Temperature	Duration	1	37°C	15 min	2	4°C	Hold			
Step	Temperature	Duration											
1	37°C	15 min											
2	4°C	Hold											

• Perform nuclease treatment for 15 min

Final cleanup with SMRTbell cleanup beads

Procedural notes



Final Kinnex library yield is typically sufficient to load **≥2 SMRT Cells**

7. Final cleanup with 1.0X SMRTbell cleanup beads

Step	Instructions
7.1	Add 50 μL (1X v/v) of resuspended, room temperature SMRTbell cleanup beads to each sample from the previous step. Pipette-mix the beads until evenly distributed.
7.2	Quick-spin the tube strip in a microcentrifuge to collect all liquid.
7.3	Incubate at room temperature for 10 minutes to allow DNA to bind to the beads.
7.4	Place the tube strip in a magnetic separation rack until the beads separate fully from the solution.
7.5	Slowly pipette off the cleared supernatant without disturbing the beads.
7.6	Slowly dispense 200 μL , or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, remove the 80% ethanol and discard.
7.7	Repeat the previous step.
7.13	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer the supernatant to a new 0.5 mL LoBind tube or a PCR tube strip. Discard the old tube strip with beads.
7.14	Take a 1 μL aliquot from each tube. Make a 1:5 dilution of the sample in Elution buffer and measure the DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit. Calculate the total mass. Expect 10–20% recovery of the starting Kinnex-PCR product.
7.14	Recommended: Further dilute each aliquot to 250 pg/ μL with EB. Measure the final SMRTbell library size distribution with a Femto Pulse system (see Figure 1).
7.15	Proceed to the SMRT® Link Sample Setup to prepare the SMRTbell library for sequencing. DNA concentration must be less than 60 ng/μL to proceed to ABC. Using a concentration >60 ng/ μL will result in low loading on the sequencer.
7.16	Store SMRTbell libraries at 4°C if sequencing within the week. Long-term storage should be at -20°C. Minimize freeze-thaw cycles when handling SMRTbell libraries.

Note: Due to diverse sources of bacterial genomic DNA, there might be contaminants that affect the sequencing performance. An additional clean-up of final SMRTbell library using 3.1X diluted Ampure PB (35% v/v, part number 100-265-900) or Monarch Genomic DNA Purification Kit (#T3010S) has been shown to remove contaminants effectively.

PROTOCOL COMPLETE



- Perform **DNA concentration QC** on final purified Kinnex RNA library using a Qubit dsDNA HS assay
 - Typical final SMRTbell library yield from 5 μg of input DNA into Kinnex array formation is **~10 – 20%** – a much higher observed yield might suggest incomplete digestion of partial SMRTbell templates
 - **Troubleshooting tip:** If SMRTbell library yield is higher than expected and *P1* loading is lower than expected, consider repeating the nuclease treatment step





- Perform **DNA sizing QC** on final purified Kinnex 16S rRNA library using a Femto Pulse system (expected final library insert size is **~17 – 18 kb**)

- Kinnex 16S rRNA final SMRTbell library concentration must be **≤60 ng/ μL** to proceed with SMRT Link sample setup (ABC¹)
 - Using a concentration above 60 ng/ μL will result in lower loading during sequencing

Kinnex 16S rRNA library prep inputs & expected step yields

Final Kinnex library yield is typically sufficient to load ≥ 2 SMRT Cells



Example Femto Pulse DNA sizing QC analysis results for final Kinnex 16S rRNA library prepared with pooled 16S amplicon DNA samples generated from mock microbial community genomic DNA (ATCC MSA-1003 20 Strain Staggered Mix).

Example Kinnex 16S rRNA library prep yields

gDNA input for 16S PCR	1.1 ng
16S amplicon DNA input for Kinnex array formation	6000 ng
Post-nuclease treatment & final library cleanup yield (%) ¹	1080 ng (18%)

¹ Post-nuclease treatment & final cleanup yields typically ranged from ~10% to ~20% when using 16S amplicon DNA samples generated from mock microbial community genomic DNA (ATCC MSA-1003 20 Strain Staggered Mix).


Final Kinnex library yield is typically sufficient to load ≥ 2 SMRT Cells



Kinnex 16S rRNA sequencing preparation workflow details

Sample Setup & Run Design recommendations for Kinnex 16S rRNA libraries

SMRT Link supports Kinnex 16S rRNA sequencing preparation & analysis workflow for PacBio systems¹



SMRT Link module	Key setup parameters For Kinnex libraries	Sequel II/III system recommended settings for Kinnex libraries	Revo system recommended settings for Kinnex libraries
Sample setup	Library type	Kinnex	
	Primer	Kinnex sequencing primer	
	Binding/Polymerase kit ¹	Sequel II binding kit 3.2 (includes Kinnex sequencing primer)	Revo polymerase kit (includes Kinnex sequencing primer)
	Concentration on plate	40 – 60 pM	100 – 150 pM
Runs → Run design	Adapter / Library type	SMRTbell Adapter Design = SMRTbell Kinnex Prep Kit	Library type = Kinnex
	Movie collection time	30 hrs	24 hrs
	Use adaptive loading	YES	
	On-instrument CCS	CCS Analysis Output - Include Low Quality Reads = NO CCS Analysis Output - Include Kinetics Information = NO	Consensus Mode = MOLECULE

SMRT Link Sample Setup and Run Design for Kinnex kits video demonstration


Video demonstration of SMRT Link Sample Setup and Run Design setup procedure for Kinnex kits supporting full-length RNA sequencing, single-cell RNA sequencing and full-length 16S rRNA sequencing

Sample Setup / Sample Calculation
Sequel II binding kit 3.1/3.2, Revio polymerase kit

Conversion Calculator
Autosaved at 2023-11-20, 09:23:31 AM

+ Add Sample Group

< Sample group >	
	Copy Remove Lock Download CSV
Name	My Batch of Samples
Application	Kinnex full-length RNA
Library type	Kinnex
Polymerase / Binding kit	Revio polymerase kit
Number of samples	1 samples
SMRT Cells per sample	1 cells
Available volume per sample	20 uL
Insert size	16000 bp
Sample concentration	40 ng/uL
Cleanup anticipated yield	75 %
Concentration on plate	130 pM Recommended: 100-150 pM
Minimum pipetting volume	1 uL
Comment	

 YouTube

[Demo video for Sample Setup and Run Design for Kinnex kits \(SMRT Link v13.0+\)](#)

- Demo video for Sample Setup and Run Design for Kinnex kits in SMRT Link v13.0+
- Kinnex kits support full-length RNA sequencing (Kinnex full-length RNA kit), full-length 16S rRNA sequencing (Kinnex 16S rRNA kit) and full-length single-cell RNA sequencing (Kinnex single-cell RNA kit)

SMRT Link Sample Setup procedure for Kinnex 16S rRNA libraries



Revio system



Sequel II and Ile systems

< Sample group >		< Sample group >	
<input type="button" value="Copy"/> <input type="button" value="Remove"/> <input type="button" value="Lock"/> <input type="button" value="Download CSV"/>		<input type="button" value="Copy"/> <input type="button" value="Remove"/> <input type="button" value="Lock"/> <input type="button" value="Download CSV"/>	
Name	Kinnex full-length 16S rRNA library demo	Name	Kinnex full-length 16S rRNA library demo
Application	Kinnex 16S rRNA	Application	Kinnex full-length RNA
Library type	Kinnex	Library type	Kinnex
Polymerase / Binding kit	Revio polymerase kit	Polymerase / Binding kit	Sequel II Binding Kit 3.2
Number of samples	1 samples	Number of samples	1 samples
SMRT Cells per sample	1 cells	SMRT Cells per sample	1 cells
Available volume per sample	20 uL	Available volume per sample	20 uL
Insert size	18000 bp	Insert size	18000 bp
Sample concentration	40 ng/uL	Sample concentration	40 ng/uL
Cleanup anticipated yield	75 %	Cleanup anticipated yield	75 %
Concentration on plate	130 pM Recommended: 100-150 pM	Concentration on plate	50 pM Recommended: 40-60 pM
Minimum pipetting volume	1 uL	Minimum pipetting volume	1 uL
Comment	Kinnex library containing array of 12 FL 16S amplicon segments	Comment	Kinnex library containing array of 12 FL 16S amplicon segments

• Select **application type** to autofill fields in green

IMPORTANT: Specify **Library type = Kinnex**

• Library type field determines sequencing primer type to use for annealing step
→ Kinnex libraries require use of **Kinnex sequencing primer¹**

• Select **Revio polymerase kit** for Revio system and **Sequel II Binding Kit 3.2** for Sequel II/Ile systems

• Recommended Kinnex 16S rRNA library input concentration for sample setup is **20 – 60 ng/μL**

• Recommended OPLC range is **100 – 150 pM** for Revio system and **40 – 60 pM** for Sequel II/Ile systems

• **Recommended target P1 loading range**

• Revio system: **~50 – 70%**

• Sequel II and Ile systems: **~60 – 80%**

SMRT Link Run Design procedure for **Revio system**

Sample and run information



Kinnex 16S rRNA

▼ **Plate 1, Well A01: Kinnex 16S rRNA library demo**

Application <small>Required</small>	Kinnex 16S rRNA
Plate Well <small>Required</small>	Plate 1, Well A01
Well Name <small>Required</small>	Kinnex 16s rRNA library demo
Well Comment	
Library Type <small>Required</small>	Kinnex
Insert Size (bp) <small>Required</small>	18000
Polymerase Kit <small>Required</small>	Revio polymerase kit
Movie Acquisition Time (hours)	24

Use Adaptive Loading
 YES NO

Specify Insert Size

Recommend 24 hrs movie collection for Revio Kinnex samples

Specify Use Adaptive Loading = YES

Select desired **application type** to autofill Library Type, Polymerase Kit & Movie Acquisition Time recommended settings

Specify **Kinnex** library type (instead of Standard or AAV)¹

Specify **Revio** polymerase kit



Forward and reverse standard terminal adapters have the same structure



Forward and reverse Kinnex terminal adapters have different structures

Library Type field determines which adapter finding algorithm is used during post-primary analysis¹

¹ **Note:** When sequencing a Kinnex library sample, if 'Standard' library type is mistakenly selected instead of 'Kinnex' then a higher missing adapter rate (> 95%) and a slight degradation in barcode demultiplexing performance (~93-96% barcoded HiFi read yield) will be observed.

SMRT Link Run Design procedure for **Revio system** (cont.)

Sample indexing (barcoding) information



Kinnex 16S rRNA

Default = YES for Sample is indexed

Samples

Sample is indexed YES NO

Indexes Required: MAS SMRTbell barcoded adapters (v2)


Same Barcodes on Both Ends of Sequence YES NO

Biosample names Required: Interactively | From a File

Specify Indexes FASTA = MAS SMRTbell barcoded adapters (v2)

Specify YES for Same barcodes on both ends of sequences

Example complete Kinnex library molecule containing barcoded Kinnex terminal MAS adapters¹ at both ends



Example interactive biosample name specification for a multiplexed Kinnex library sample

Barcode Selector and Sample Name Editor

Barcode ID	Bio Sample ID
<input type="checkbox"/> bcM0001--bcM0001	
<input type="checkbox"/> bcM0002--bcM0002	
<input type="checkbox"/> bcM0003--bcM0003	
<input type="checkbox"/> bcM0004--bcM0004	

Barcode Selector and Sample Name Editor

Barcode ID	Bio Sample ID
<input type="checkbox"/> bcM0003--bcM0003	
<input type="checkbox"/> bcM0004--bcM0004	
<input type="checkbox"/> bcM0001--bcM0001	Kinnex adapter-barcoded library 1
<input type="checkbox"/> bcM0002--bcM0002	Kinnex adapter-barcoded library 2

SMRT Link



MAS SMRTbell barcoded adapter indexes

```
>bcM0001
ACAGTC
>bcM0002
ATGACG
>bcM0003
CACGTG
>bcM0004
CATCGC
```

SMRT Link Run Design procedure for **Revio system** (cont.)

Run options and data options



Kinnex 16S rRNA

Run Options

Library Concentration (pM)
Required

On-plate loading concentration is required for Revio samples

Data Options

Include Base Kinetics YES NO

Consensus Mode MOLECULE STRAND

Assign Data To Project [?](#)

Default = NO for Include Base Kinetics

Default Consensus Mode = MOLECULE¹

Can leave Include Base Kinetics and Consensus Mode fields at their default settings for Kinnex library samples

SMRT Link Run Design procedure for Sequel II/Ile systems

Sample information and run information



Kinnex 16S rRNA

- Select desired **Kinnex application** from the **Application** field drop-down menu
- The following fields are **auto-populated** with default recommended values and high-lighted in **green**:

SMRTbell Adapter Design

→ SMRTbell Kinnex Prep Kit

Binding Kit

→ Sequel II Binding Kit 3.2

Sequencing Kit

→ Sequel II Sequencing Plate 2.0 (4 rxn or 1 rxn)

DNA Control Complex

→ Sequel II DNA Internal Control Complex 3.2

Movie Time per SMRT Cell

→ 30 hrs

Pre-Extension Time

→ 2 hrs

SMRTbell Adapter Design field determines which adapter finding algorithm is used during post-primary analysis¹

Default SMRTbell adapter design for Kinnex samples is **SMRTbell Kinnex Prep Kit**

Recommended OPLC for Sequel II/Ile Kinnex library samples is **40 – 60 pM**

Recommended movie time = **30 hrs**

Select desired Kinnex application type from drop-down menu

SAMPLE 1: Kinnex full-length RNA library demo , A01, 30 hour movie, 16000 bp insert

Copy Delete

Import from Sample Setup Select Sample

Application Required Kinnex 16S rRNA

Well Sample Name Required Kinnex 16S rRNA library demo

Bio Sample Name

Sample Comment

Sample Well A01

SMRTbell Adapter Design Required SMRTbell® Kinnex Prep Kit

Binding Kit Required Sequel® II Binding Kit 3.2

Sequencing Kit Required Sequel® II Sequencing Plate 2.0 (4 rxn)

DNA Control Complex Sequel® II DNA Internal Control Complex 3.2

Insert Size (bp) Required 18000

Recommended Concentration on Plate (pM) 40 – 60 pM

On-Plate Loading Concentration (pM) Required 45

Movie Time per SMRT Cell (hours) 30

Use Pre-Extension YES NO

Pre-Extension Time (hours) 2

CCS Analysis will be performed on-instrument to produce HiFi .bam files.

Example sample information entered into a Sequel IIe system run design worksheet for a Kinnex 16S rRNA library sample.

SMRT Link Run Design procedure for Sequel II/Ie systems (cont.)

Advanced options



Kinnex 16S rRNA

- For all Kinnex library samples, leave the following **Advanced Options** fields at their **default settings**
 - Use Adaptive Loading**
→ YES
 - Loading Target (P1 + P2)**
→ 0.85
 - Maximum Loading Time**
→ 2 hours
 - CCS Analysis Output - Include Low Quality Reads**
→ NO
 - CCS Analysis Output - Include Kinetics Information**
→ NO
 - Pre-Extension Time**
→ 2 hrs
- If desired, specify to use an alternative project folder for the **Add Data to Project** field

Advanced Options

Use Adaptive Loading YES NO

Loading Target (P1 + P2)

Maximum Loading Time (hours)

CCS Analysis Output - Include Low Quality Reads YES NO

CCS Analysis Output - Include Kinetics Information YES NO

Add Data to Project

Example default Advanced Options settings entered into a Sequel IIe system run design worksheet for a Kinnex 16S rRNA library sample.

SMRT Link Run Design procedure for Sequel II/Ile systems (cont.)

Barcoded sample options



Kinnex 16S rRNA

- For Kinnex library samples, can leave most **Barcoded Sample Options** fields at their **default settings**

Specify Bio Sample Names, either interactively or by downloading a CSV file (**Interactively** or **From a file**)

If desired, specify to perform barcode demultiplexing on-instrument or in SMRT Link (default = **On-instrument** for Sequel IIe system)

Can leave most of these fields at their default values

Barcoded Sample Options

Sample Is Barcoded YES NO

Barcode Set Required MAS SMRTbell barcoded adapters (v2)

Same Barcodes on Both Ends of Sequence YES NO

Assign Bio Sample Names to Barcodes Required

Demultiplex Barcodes ON INSTRUMENT IN SMRT LINK DO NOT GENERATE

Example default Barcoded Sample Options settings entered into a Sequel IIe system run design worksheet for a Kinnex full-length RNA library sample.

SMRT Link Run Design analysis options for **Revio system** and **Sequel II/IIe systems**

Kinnex 16S rRNA

Analysis Options Default = YES for Add Analysis

Add Analysis YES NO


Analysis Name Required Kinnex_16S_rRNA_Demo_Analysis_Job_Name

Select Analysis Workflow Required Read Segmentation

Segmentation Adapter Set MAS-Seq Adapter v2 (MAS12)

Advanced Parameters

Read Segmentation



Read Segmentation data utility splits arrayed HiFi reads at segmentation adapter positions, generating segmented reads (S-reads) which are the comprising fragments.

SMRT Link Run Design analysis options for **Revio system** and **Sequel II/IIe systems** (cont.)



Analysis Options

Add Analysis YES NO


Analysis Name
Required Kinnex_16S_rRNA_Demo_Analysis_Job_Name

Select Analysis Workflow
Required Read Segmentation

Segmentation Adapter Set
MAS-Seq Adapter v2 (MAS12)

Advanced Parameters

Specify **Segmentation Adapter Set** that corresponds to the Kinnex library concatenation method used
→ For Kinnex 16S rRNA samples, specify **MAS-Seq Adapter v2 (MAS12)**






Kinnex 16S rRNA example sequencing performance data

Example Kinnex 16S rRNA library preparation QC results

Kinnex full-length 16S RNA library prepared from mock microbial community genomic DNA


Final Kinnex library yield is typically sufficient to load ≥ 2 SMRT Cells

16S amplicon DNA QC



Example Bioanalyzer DNA sizing QC analysis results for pooled 16S amplicon DNA samples generated from mock microbial community genomic DNA (ATCC MSA-1003 20 Strain Staggered Mix).

Final Kinnex 16S rRNA library QC



Example Femto Pulse DNA sizing QC analysis results for final Kinnex 16S rRNA library.

Sample multiplexing design for Kinnex 16S rRNA library

384-plex 16S PCR x 4-plex Kinnex adapter barcoding \rightarrow Total sample multiplex level = **1,536-plex**

 ZymoBIOMICS Fecal Reference 96 barcoded samples	 ZymoBIOMICS Gut Microbiome Standard 96 barcoded samples	 ATCC MSA-1002 20 Strain Even Mix 96 barcoded samples	 ATCC MSA-1003 20 Strain Staggered Mix 96 barcoded samples	X 4 Kinnex barcoded adapters
--	--	---	---	-------------------------------------


gDNA input for 16S PCR	1.1 ng
16S amplicon DNA input Kinnex PCR products for Kinnex array formation	6000 ng
Post-nuclease treatment & final library cleanup yield (%) ¹	1080 ng (18%)

¹ Post-nuclease treatment & final cleanup yields typically ranged from ~10% to ~20% when using mock microbial community genomic DNA for Kinnex full-length 16S rRNA library construction.


Example sequencing performance for Kinnex 16S rRNA libraries prepared from mock microbial community genomic DNA

Revio system example data¹ (1,536-plex data set)


Raw Data Report



HiFi Read Length



Read Segmentation Metrics



For Revio system Kinnex 16S rRNA samples

- Yield of HiFi reads is ideally >4.5 Million
- Yield of segmented reads is ideally >45 Million
- Mean S-read length should be consistent with expected 16S amplicon length (~1.5 kb)
- % of reads with full arrays is ideally >80%
- Mean array size is ideally >11.0 segments

Raw Base Yield	1,222 Gb
Mean Polymerase Read Length	69.0
P0	27%
P1	70%
P2	3%

Example sequencing metrics for a Kinnex 16S rRNA library sample run on a Revio system with Revio polymerase kit / 130 pM on-plate loading concentration (OPLC) / 24-hrs movie time.

HiFi Reads	6.1 M
HiFi Base Yield	114.21 Gb
Mean HiFi Read Length	18.78 kb
Median HiFi Read Quality	Q32
HiFi Read Mean # of Passes	7

For Kinnex 16S rRNA libraries, per-Revio SMRT Cell HiFi read counts were typically ~4 – 6 Million depending on the final library insert size and P1 loading performance.

Input HiFi Reads	6,050,730
Segmented reads (S-reads)	71,720,714
Mean length of S-reads	1,560 bp
Percent of reads with full arrays	95.03%
Mean array size (concentration factor)	11.85

For Kinnex 16S rRNA libraries, per-Revio SMRT Cell segmentation read counts were typically ~45 – 60 Million.

Improving sequencing performance of “difficult” 16S samples

Performing AMPure PB bead size-selection on Kinnex full-length 16S rRNA libraries can help improve *P1* loading of challenging metagenomic samples

Sample Name	P1 %	Gb Yield	Mean Length	Mean QV
16S_collaborator_SOP	26	46	18,813 bp	Q29
16S_collaborator_3.1X AMPure	80	87	18,851 bp	Q28

Some bacterial 16S samples may have carry-over contaminants present leading to low *P1* loading on Revo and Sequel II/IIe systems

→ Using AMPure size-selection (3.1X 35% AMPure PB beads) can help mitigate this issue




Kinnex 16S rRNA data analysis workflow overview

Kinnex 16S rRNA bioinformatics workflow overview

SMRT Link Read Segmentation processes HiFi reads generated from Kinnex 16S rRNA libraries to produce individual segmented reads (S-reads) that are compatible with tertiary analysis tools

SMRT Link read segmentation and barcode demultiplexing workflow¹



Read segmentation

- Use SMRT Link Read Segmentation data utility to split arrayed Kinnex library HiFi reads into individual segmented reads (**S-reads**) that represent the original barcoded 16S amplicon sequences

Demultiplex barcodes

- Use SMRT Link Demultiplex barcodes data utility to separate sequence reads by barcode

Tertiary analysis

- Use Github or other third-party tools (e.g., [PacBio GitHub](#), [DADA2](#), [QIIME2](#), [microbiomehelper](#), [One Codex](#), [EZBiome](#), [BugSeq](#))

Example Kinnex 16S rRNA data set

- ZymoBIOMICS Fecal Reference with TruMatrix Technology (human) [[Link](#)]
- Mixture: ZymoBIOMICS Gut Microbiome Standard, ZymoBIOMICS Fecal Reference with TruMatrix™ Technology, ATCC 20 Strain Even Mix Genomic Material, ATCC 20 Strain Staggered Mix Genomic Material [[Link](#)]

SMRT Link Read Segmentation and Demultiplex Barcodes video demonstration

Video demonstration of SMRT Link Read Segmentation and Demultiplex Barcodes workflow for analysis of Kinnex 16S rRNA samples

Data Utility Required

Read Segmentation

Import Analysis Settings Export

Associated Inputs

Segmentation Adapter Set

MAS-Seq Adapter v2 (MAS12)

Advanced Parameters

Analysis Name

test-ReadSeg


Analysis Datasets

Displaying rows 1 to 1 out of 1

ID	Name
21...	20231020-4_84028_13...

[Demo video](#) for Read Segmentation and Demultiplex Barcodes workflow (SMRT Link v13.0+)

- Demo video for analyzing Kinnex 16S rRNA data generated using Kinnex 16S rRNA kit in SMRT Link v13.0 and up



SMRT Link Read Segmentation setup

Specify **Read Segmentation** data utility in SMRT Link¹

Read Segmentation data utility splits arrayed Kinnex HiFi library reads at segmentation adapter positions, generating segmented reads (S-reads) which are the comprising fragments

- Accepts **HiFi reads** (BAM format) as input
 - HiFi reads are reads generated with CCS analysis whose quality value is equal to or greater than 20
- HiFi reads should be generated using the Kinnex 16S rRNA library preparation protocol ([103-238-800](#))
 - If the library is a regular (non-Kinnex) 16S monomer library without MAS-Seq concatenation, Read Segmentation does not need to be performed.

- For Kinnex 16S rRNA libraries, specify Segmentation Adapter Set = **MAS-Seq Adapter v2 (MAS12)**
- Specify a FASTA file, provided by PacBio, containing segmentation adapters. If you need a custom segmentation adapter set, click Advanced Parameters and use a custom FASTA file formatted as described in the SMRT Link User Guide

Example SMRT Link Read Segmentation data utility processing results¹ for Kinnex 16S rRNA library prepared from mock microbial community gDNA

SMRT Link Read Segmentation data utility job report – Summary Metrics and Segmentation Statistics

Summary Metrics


Value	Analysis Metric
6,050,730	Reads
71,720,714	Segmented reads (S-reads)
1,560	Mean length of S-reads
95.03 %	Percent of reads with full arrays
11.85	Mean array size (concatenation factor)

Example Revio system data shown.

For Revio system Kinnex 16S rRNA samples

- Yield of HiFi reads is ideally >4.5 Million
- Yield of segmented reads is ideally >45 Million
- Mean S-read length should be consistent with expected 16S amplicon length (~1.5 kb)
- % of reads with full arrays is ideally >80%
- Mean array size is ideally >11.0 segments

Segmentation Statistics



• Plot should ideally show a majority (>80%) of reads with array lengths = 12 segments

- **Reads:** Number of input arrayed HiFi reads
- **Segmented reads (S-reads):** Number of generated S-reads
- **Mean length of S-reads:** Mean read length of generated S-reads
- **Percent of reads with full arrays:** Percentage of input HiFi reads containing all adapter sequences in the order listed in the segmentation adapter FASTA file
- **Mean array size:** Mean number of fragments (or S-reads) found in input reads


Histogram distribution of the number of S-reads per HiFi read. (Example Revio system data shown.)

¹ HiFi read lengths, reads/data per SMRT Cell and other sequencing performance results can vary depending on DNA sample quality, insert size, P1 loading performance & movie time. For Sequel IIe systems, we recommend aiming for ~60 – 80% P1 loading. For Revio system, we recommend aiming for ~50 – 70% P1 loading.

Example SMRT Link Read Segmentation data utility processing results¹ for Kinnex 16S rRNA library prepared from mock microbial community gDNA (cont.)


SMRT Link Read Segmentation data utility job report – Length of Reads and S-read Length Distribution

Length of Reads



Histogram distribution of the number of HiFi reads by read length, in base pairs. (Example Revio system data shown.)

S-read Length Distribution



Histogram distribution of the number of S-reads by HiFi read length, in base pairs. (Example Revio system data shown.)



Technical documentation & applications support resources

Technical resources for Kinnex full-length 16S rRNA library preparation, sequencing & data analysis

Metagenomic DNA sample preparation resources

- Based on prior PacBio experience, the PowerFecal Pro DNA Kit (QIAGEN) extracts DNA of sufficient quality for the Kinnex full-length 16S rRNA library preparation workflow.

Kinnex 16S rRNA library preparation literature & other resources

- Application note – Kinnex 16S rRNA kit for full-length 16S sequencing ([102-326-601](#))
- Brochure – Scalable, cost-effective RNA sequencing with PacBio Kinnex kits ([102-326-597](#))
- Procedure & checklist – Preparing Kinnex libraries using Kinnex 16S rRNA kit ([103-238-700](#))
- Technical overview – Kinnex kits for single-cell RNA, full-length RNA and 16S rRNA sequencing ([103-343-700](#))
- Technical overview – Kinnex library preparation using Kinnex 16S rRNA kit ([103-344-800](#))
- Video tutorial – SMRT Link Sample Setup and Run Design setup procedure for Kinnex kits [[Link](#)]

Data analysis resources

- SMRT Link MAS-Seq troubleshooting guide ([102-994-400](#))
- SMRT Link software installation guide [[Link](#)]
- SMRT Link user guide [[Link](#)]
- SMRT Tools reference guide [[Link](#)]
- Video tutorial – Analyzing Kinnex 16S rRNA data in SMRT Link [[Link](#)]

Technical resources for Kinnex full-length 16S rRNA library preparation, sequencing & data analysis

Posters & publications

- ASM Microbe Poster (2023) – Increasing throughput of full-length 16S sequencing using concatenation [[Link](#)]

Webinars

- PacBio webinar (2023) – Creating a reliable microbiome testing service for companion animal health using PacBio sequencing [[Link](#)]

Example PacBio data sets

Application	Dataset	Data type	PacBio system
Kinnex 16S rRNA sequencing	ZymoBIOMICS Fecal Reference with TruMatrix Technology (human) [Link]	HiFi long read	Sequel II & Revio systems
	Mixture: ZymoBIOMICS Gut Microbiome Standard, ZymoBIOMICS Fecal Reference with TruMatrix™ Technology, ATCC 20 Strain Even Mix Genomic Material, ATCC 20 Strain Staggered Mix Genomic Material [Link]	HiFi long read	Sequel II & Revio systems



www.pacb.com

Research use only. Not for use in diagnostic procedures. © 2024 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, and PureTarget are trademarks of PacBio.