

The background of the slide is a blurred laboratory setting. In the upper right, a pipette tip is shown dispensing a drop of red liquid into a well of a multi-well plate. The plate contains several other wells, some of which are filled with red liquid. The overall scene is brightly lit, with a soft focus on the laboratory equipment.

**PacBio**

# **Technical overview – Kinnex library preparation using Kinnex full-length RNA kit**

Sequel II and IIe systems ICS v11.0

Revio system ICS v13.0+

SMRT Link v13.0+

PN 103-344-700 Rev 01 | March 2024


# Kinnex library preparation using Kinnex full-length RNA kit



## Technical Overview

1. Kinnex full-length RNA method overview
2. Kinnex full-length RNA library preparation workflow details
3. Kinnex full-length RNA sequencing preparation workflow details
4. Kinnex full-length RNA example sequencing performance data
5. Kinnex full-length RNA data analysis workflow overview
6. Technical documentation & applications support resources

# Kinnex library preparation using Kinnex full-length RNA kit: Getting started



**Application note**

## Kinnex full-length RNA kit for isoform sequencing

**Introduction**

Alternative splicing (AS) in eukaryotic species generates functional diversity by expressing different combinations of exons in the same gene. Accurate characterization of full-length transcript isoforms generated by AS is critical for biological and disease studies. Bulk RNA-Seq using short reads cannot fully resolve isoform structures, as the complex nature of AS prohibits unambiguous transcript assembly with even the most sophisticated computational tools [1, 2]. Long-read RNA-Seq using PacBio® technology (the Iso-Seq® method) eliminates the need for transcript assembly by sequencing full-length cDNAs and enables new discoveries across many applications (Figure 1).

The Kinnex™ full-length RNA kit takes total RNA as input and outputs a sequencing-ready library that results in an 8-fold throughput increase compared to typical Iso-Seq libraries. Combined with the Iso-Seq analysis in SMRT™ Link software, PacBio offers cost-effective isoform sequencing that does not require orthogonal sequencing methods. SMRT Link software produces an isoform classification report with abundance information that can be used by tertiary analysis tools [3].

© 2023 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures. 102-326-591 REV02 NOV2023

## Application note – Kinnex full-length RNA kit for isoform sequencing (102-326-591)

Summary overview of application-specific library preparation and data analysis workflow recommendations.

**Preparing Kinnex™ libraries using the Kinnex full-length RNA kit**

Procedure & checklist

**Before you begin**

This procedure describes the workflow for constructing Kinnex full-length RNA libraries from total RNA samples for sequencing on PacBio® Sequel II, Sequel IIe, and Revio™ systems.

Overview	
Samples	1–24
Workflow time	8 hours (for up to 24 samples)
Number of SMRT™ Cells per Kinnex Library Prep	Up to 2 SMRT Cells for Revio system Up to 4 SMRT Cells for Sequel II/Ile systems

RNA input	
Quality/size distribution	RIN (RNA integrity number) >7.0
Quantity	300 ng per library (minimum concentration 43 ng/μL per library)

© 2023 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures. 103-238-700 REV01 OCT 2023

## Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit (103-238-700)

Technical documentation containing application-specific library preparation protocol details.

**Technical overview – Kinnex library preparation using Kinnex full-length RNA kit**

Sequel II and Ile systems ICS v11.0  
Revio system ICS v13.0  
SMRT Link v13.0

PN 103-344-700 Rev 01 | December 2023

**Example sequencing performance for Kinnex full-length RNA libraries prepared with human cDNA**

Revio system example data\*

Raw Data Report	HiFi Read Length	Read Segmentation Metrics
<p>Mean HiFi Read Length: 18.41 kb</p> <p>Mean HiFi Read Quality: 99.9%</p>	<p>HiFi Reads: 8.3 M</p> <p>Mean HiFi Read Length: 18.41 kb</p> <p>Mean HiFi Read Quality: 99.9%</p>	<p>Input HiFi Reads: 8,000,154</p> <p>Mean length of reads: 17,216 bp</p> <p>Mean length of reads: 2,200 bp</p> <p>Percent of reads with full entry: 83.6%</p> <p>Mean entry size (concatenation depth): 7.40</p>

© 2023 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures. 103-344-700 REV01 OCT 2023

## Technical Overview – Kinnex library preparation using Kinnex full-length RNA kit (103-344-700)

Technical overview presentations describe sample preparation details for constructing Kinnex HiFi libraries for specific applications. Example sequencing performance data for a given application are also summarized.

**cDNA synthesis & amplification**  
(Iso-Seq express 2.0 kit)

300 ng input total RNA per sample  
RIN (RNA integrity number) ≥7.0  
Sample multiplexing options available

**Kinnex library preparation**  
(Kinnex PCR 8-fold kit + Kinnex concatenation kit)

Use amplified cDNA to generate Kinnex library containing 8-segment array  
Sample multiplexing up to 48-plex supported

**SMRT sequencing**  
(Sequel II/Ile & Revio systems)

Perform ABC\* and sequence Kinnex libraries on PacBio long-read systems

**Data analysis**  
(SMRT Link)


Use SMRT Link Read Segmentation data utility to split arrayed transcript HiFi reads  
Use SMRT Link Iso-Seq analysis application to identify novel genes and isoforms with abundance information



# Kinnex full-length RNA method overview

# Kinnex full-length RNA method overview

Use Kinnex full-length RNA kit to perform high-accuracy, full-length isoform sequencing with PacBio long-read systems




- Input 300 ng total RNA, RIN  $\geq 7$
- Generate up to 12-plex barcoded cDNA using **Iso-Seq express 2.0 kit (103-071-500)**
- 2-day Kinnex library preparation using **Kinnex full-length RNA kit (103-072-000)**
- SMRT Link Run Design support for 'Kinnex full-length RNA' application type with auto-analysis (read segmentation + isoform analysis)<sup>1</sup>
- SMRT Link Iso-Seq isoform-classification software to identify novel genes and isoforms with abundance information




# Kinnex full-length RNA method overview (cont.)

## Full-length cDNA synthesis & amplification



## Kinnex library prep, sequencing & analysis



<sup>1</sup> Twelve barcoded Iso-Seq primers (Iso-Seq primer bc01–12) are available for cDNA amplification step.

<sup>2</sup> Kinnex adapter barcode sequences can be downloaded from [SMRT Link](#) Data Management module.

# Kinnex full-length RNA library preparation procedure description

Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit ([103-238-700](#)) describes the workflow for constructing Kinnex libraries from total RNA samples using the **Iso-Seq express 2.0 kit** and **Kinnex full-length RNA kit** for sequencing on PacBio Sequel II, Sequel IIe, and Revio systems

Overview	
Samples	1–24
Workflow time	8 hours (for up to 24 samples)
Number of SMRT® Cells per Kinnex library Prep	Up to 2 SMRT Cells for Revio system Up to 4 SMRT Cells for Sequel II/IIe systems
RNA input	
Quality/size distribution	RIN (RNA integrity number) ≥7.0
Quantity	300 ng per library (minimum concentration 43 ng/μL per library)



**PacBio**

Preparing Kinnex™ libraries using the Kinnex full-length RNA kit

Procedure & checklist

---

Before you begin

This procedure describes the workflow for constructing Kinnex full-length RNA libraries from total RNA samples for sequencing on PacBio® Sequel® II, Sequel IIe, and Revio™ systems.

Overview	
Samples	1–24
Workflow time	8 hours (for up to 24 samples)
Number of SMRT® Cells per Kinnex library Prep	Up to 2 SMRT Cells for Revio system Up to 4 SMRT Cells for Sequel II/IIe systems
RNA input	
Quality/size distribution	RIN (RNA integrity number) ≥7.0
Quantity	300 ng per library (minimum concentration 43 ng/μL per library)

© 2023 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.  
103-238-700 REV 01 OCT 2023 **PacBio**

PacBio [Documentation](#) ([103-238-700](#))

- Kinnex full-length RNA library prep protocol uses **Kinnex full-length RNA kit** and **Iso-Seq express 2.0 kit**  
→ **Do not use** SMRTbell prep kit 3.0 with this protocol








# Kinnex full-length RNA kit bundle and Iso-Seq express 2.0 kit components

Kinnex full-length RNA kit bundle and Iso-Seq express 2.0 kit provide full support for Kinnex library prep workflow

## Iso-Seq express 2.0 kit (103-071-500)

Includes Iso-Seq Express template switching oligo, barcoded cDNA PCR Primers, and other reagents needed for performing 1<sup>st</sup>-strand cDNA synthesis and PCR amplification of cDNA products generated from input total RNA.





### Iso-Seq express 2.0 kit components

Component	Description
1 	<b>Iso-Seq RT buffer</b> <ul style="list-style-type: none"> <li>For 1<sup>st</sup>-strand cDNA synthesis</li> </ul>
2 	<b>Iso-Seq RT primer mix</b> <ul style="list-style-type: none"> <li>For 1<sup>st</sup>-strand cDNA synthesis</li> </ul>
3 	<b>Iso-Seq RT enzyme mix</b> <ul style="list-style-type: none"> <li>For 1<sup>st</sup>-strand cDNA synthesis</li> </ul>
4 	<b>Iso-Seq template switch oligo</b> <ul style="list-style-type: none"> <li>For 1<sup>st</sup>-strand cDNA synthesis</li> </ul>
5 	<b>Iso-Seq cDNA PCR mix</b> <ul style="list-style-type: none"> <li>Enzyme nucleotide mix for cDNA amplification</li> </ul>
6 	<b>Iso-Seq cDNA amplification primer</b> <ul style="list-style-type: none"> <li>Reverse primer for cDNA amplification</li> </ul>
7 	<b>Iso-Seq primers (bc01 – bc12)</b> <ul style="list-style-type: none"> <li>Barcoded forward primers (bc01 – bc04) for cDNA amplification</li> </ul>

## Kinnex full-length RNA kit bundle (103-072-000)

Includes Kinnex PCR kit, Kinnex concatenation and ancillary DNA cleanup reagents needed for incorporation of Kinnex segmentation adapters and Kinnex array formation for generating Kinnex full-length RNA libraries.

### Kinnex full-length RNA kit components

Component	Description
1 	<b>Kinnex PCR 8-fold kit (12 rxn)</b> <ul style="list-style-type: none"> <li>Contains reagents for Kinnex PCR to incorporate segmentation adapters</li> </ul>
2 	<b>Kinnex concatenation kit (12 rxn)</b> <ul style="list-style-type: none"> <li>Contains reagents for Kinnex array formation and SMRTbell template construction</li> <li>Includes barcoded Kinnex adapter mixes (bcM0001 – bcM0004)</li> </ul>
3 	<b>SMRTbell cleanup beads</b> <ul style="list-style-type: none"> <li>For DNA cleanup</li> </ul>
4 	<b>Elution buffer</b> <ul style="list-style-type: none"> <li>For DNA cleanup</li> </ul>




# Kinnex full-length RNA library barcoding options for sample multiplexing

Kinnex full-length RNA library preparation procedure supports up to 48-plex sample multiplexing

Kinnex full-length RNA library preparation procedure supports **up to 48-plex** sample multiplexing through combined use of:

- 12 different barcoded cDNA amplification PCR primers (bc01 – bc12)
- 4 different barcoded Kinnex terminal SMRTbell adapters (bc01 – bc04)



# Kinnex full-length RNA experimental design considerations

## Kinnex full-length RNA application use case recommendations for PacBio systems

Example application	Human genetics disease studies	Biopharma for identifying highly expressed targets	Plant & animal whole genome annotation
Experimental goal	Isoform discovery and quantification of moderate-to-rare transcripts	Isoform discovery of high expressed transcripts	Comprehensive transcript annotation in a species
Example study design	Disease vs. normal tissues with multiple replicates	Disease cohort with >20+ samples	Plant or animal with multiple tissue types
Target depth of coverage per sample	10 M reads per sample	5 M reads per sample	≤5 M reads per tissue (of same species)
Sample multiplexing <sup>1</sup>	<b>Sequel II/IIe system:</b> Up to 2 samples per SMRT Cell 8M (2-plex)	<b>Sequel II/IIe system:</b> Up to 3 samples per SMRT Cell 8M (3-plex)	<b>Sequel II/IIe system:</b> Up to 3 tissue types per SMRT Cell 8M (3-plex)
	<b>Revio system:</b> Up to 4 samples per Revio SMRT Cell (4-plex)	<b>Revio system:</b> Up to 8 samples per Revio SMRT Cell (8-plex)	<b>Revio system:</b> Up to 8 tissue types per Revio SMRT Cell (8-plex) <sup>2</sup>
Expected data throughput (per SMRT Cell)	<b>Sequel II/IIe system:</b> 15 M reads per SMRT Cell 8M divided by <i>N</i> samples <b>Revio system:</b> 40 M reads divided by <i>N</i> samples		
Kinnex library prep protocol	Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit ( <a href="https://www.pacb.com/products-services/kinnex/">103-238-700</a> )		
Total RNA input into Kinnex library prep workflow	300 ng total RNA (RIN ≥7) for 1 <sup>st</sup> -strand cDNA synthesis		
SMRT Link data analysis workflows	Read Segmentation and Iso-Seq analysis application with option to “pool reads and cluster together” to get a master isoform classification file with per-sample full-length read counts		
Community data analysis tools	Annotation & quantification: PIGEON, SQANTI3, Differential analysis: TappAS, Fusion calling: pbfusion, Visualization: SWAN		

<sup>1</sup> Kinnex concatenation kit (103-071-800) can support up to 48-plex sample multiplexing through the combined use of 12 different barcoded cDNA amplification primers and 4 different barcoded Kinnex terminal SMRTbell adapters during Kinnex full-length RNA library construction.

<sup>2</sup> If targeting <5 M transcripts reads per sample → can multiplex up to 12 tissues types per Revio SMRT Cell.



# **Kinnex full-length RNA library preparation workflow details**

# Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit (103-238-700)

Procedure & checklist [103-238-700](#) describes the workflow for constructing Kinnex libraries from total RNA samples using the **Iso-Seq express 2.0 kit** and **Kinnex full-length RNA kit** for sequencing on PacBio Sequel II, Sequel IIe, and Revio systems<sup>1</sup>

## Procedure & checklist contents

1. Total RNA input QC recommendations and general best practices for reagent & sample handling.
2. Enzymatic workflow steps for cDNA synthesis and amplification.
3. Enzymatic workflow steps for construction of 8-segment Kinnex arrays from amplified cDNA.
4. Enzymatic workflow steps for DNA damage repair & nuclease treatment of Kinnex libraries.
5. Workflow steps for final cleanup of Kinnex SMRTbell libraries using SMRTbell cleanup beads.

Preparing Kinnex™ libraries using the Kinnex full-length RNA kit **PacBio**

Procedure & checklist

### Before you begin

This procedure describes the workflow for constructing Kinnex full-length RNA libraries from total RNA samples for sequencing on PacBio® Sequel® II, Sequel IIe, and Revio™ systems.

Overview	
Samples	1–24
Workflow time	1.5 days (for up to 24 samples)
Number of SMRT® Cells per Kinnex library Prep	>2 SMRT Cells for Revio system >4 SMRT Cells for Sequel II/IIe systems


RNA input	
Quality/size distribution	RIN (RNA integrity number) ≥7.0
Quantity	300 ng per library (minimum concentration 43 ng/μL per library)

© 2024 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.  
103-238-700 REV 03 MAR2024 **PacBio**

PacBio [Documentation](#) ([103-238-700](#))

# Kinnex full-length RNA library construction workflow overview

Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit (103-238-700)



\* Excludes hands-on time for setting up enzymatic reaction steps and additional time required to perform DNA sizing QC and DNA concentration QC.




# General best practices recommendations for preparing Kinnex full-length RNA libraries

## Reagent and sample handling



- Take care to accurately pipette SMRTbell cleanup beads because small changes in volume can significantly alter the size distribution of your sample.
- Equilibrate the SMRTbell cleanup beads at room temperature for 30 mins prior to use.
- In cDNA amplification and Kinnex PCR, **keep sample(s) on ice until thermal cycler lid has reached 105°C** to avoid digestion of primers by polymerase exonuclease activity.
- This workflow takes ~8 hrs to complete.
  - If a stop is necessary, refer to the workflow for **safe stopping points**.

# Input RNA quality control

Input RNA quality control is highly recommended before proceeding to the MAS-Seq library prep workflow




- 300 ng of total RNA per library (minimum concentration 43 ng/μL per library) is required for this procedure
- Sample QC of input total RNA samples should be assessed by measuring RNA Integrity Number (RIN) using a Bioanalyzer 2100 instrument (Agilent Technology) with RNA 6000 Nano kit<sup>1</sup>
  - RIN ≥7.0 (ideally ≥8.0) is sufficient for Kinnex full-length RNA protocol
  - Samples with RIN <7.0 can be processed, but risk of significant underperformance or even failure is greatly increased
- RIN score (1 to 10) is related to ratio of the area under 28s and 18s fragment peaks and also takes into account signal intensity above baseline in the Inter-Region and Fast Region since this is where degradation products appear
- Higher RIN numbers are correlated with better overall sample quality and lower degradation




**Left:** Bioanalyzer electropherogram detailing the regions that are indicative of RNA quality. **Right:** Sample electro-pherograms corresponding to different RNA Integrity Number (RIN) scores. Samples range from intact (RIN 10), to degraded (RIN 2). Images from Agilent Application Note: RNA Integrity Number (RIN) – Standardization of RNA Quality Control ([5989-1165EN](https://www.agilent.com/chem/5989-1165EN))

# cDNA synthesis + cleanup

In this step, total RNA samples are converted to first-strand cDNA products




- Minimum total RNA input for cDNA synthesis = 300 ng



# cDNA synthesis + cleanup (cont.)

## Procedural notes



### 2.1 Thawing reagents for first-strand synthesis

Step	Instructions
2.1.1	Quick-spin the Iso-Seq RT enzyme mix in the microcentrifuge to collect liquid, then place on ice. Thaw the following components at room temperature, briefly vortex to mix, then quick-spin to collect liquid and place on ice.
Reagent	
	Iso-Seq RT primer mix (103-104-000)

- After thawing specified reagents **on ice**, perform a quick spin to collect liquid, then **place on ice**

### 2.2 Primer annealing for first-strand synthesis

Step	Instructions
2.2.1	For each RNA sample to be processed, prepare reagent mix 1 on ice by adding the following components to each tube in the PCR strip tube.
Components	Volume
Total RNA (300 ng)	<7 $\mu$ L

- Set up primer annealing reaction mix **on ice** and then transfer to thermal cycler for incubation

### 2.3 Reverse transcription and template switching

Step	Instructions
2.3.1	For each RNA sample, prepare reagent mix 2 on ice by adding the following components in the order and volume listed below. Adjust component volumes for the number of samples being prepared, plus 10% overage.
Components	Volume
Iso-Seq RT buffer (vortex briefly before use)	5 $\mu$ L

- Set up reverse transcription reaction mix **on ice** and then transfer to thermal cycler for incubation


### 2.4 1.3X SMRTbell bead cleanup

Step	Instructions
2.4.1	For each sample, add 29 $\mu$ L of elution buffer to the 21 $\mu$ L reverse transcription and template switching reaction ( <a href="#">Section 2.3</a> ) for a total volume of 50 $\mu$ L.
2.4.2	Add 65 $\mu$ L of resuspended, room-temperature SMRTbell cleanup beads.
2.4.3	Mix beads by pipetting 10 times or until evenly distributed.

- Perform **1.3X** SMRTbell bead cleanup

# cDNA amplification + cleanup

In this step, first-strand cDNA products are PCR-amplified and barcoded using barcoded Iso-Seq primers



Can barcode up to 12 samples during cDNA amplification step using PacBio barcoded Forward Iso-Seq primers<sup>1</sup>




2<sup>nd</sup>-strand synthesis & cDNA amplification

Forward (F) PCR Primer  
(Barcoded Iso-Seq primer bcXX<sup>1</sup>)



Reverse (R) PCR Primer  
(Iso-Seq cDNA amplification primer)

0.9X SMRTbell bead cleanup



Pool barcoded cDNA  
(if multiplexing)




Pooled barcoded full-length cDNA for input into Kinnex library construction



# cDNA amplification + cleanup

## Procedural notes



### 3.1 cDNA amplification

Step	Instructions														
3.1.1	For each sample, prepare reaction mix 3 on ice by adding the following components in the order and volume listed below. Adjust component volumes for the number of samples being prepared, plus 10% overage. Pipette mix master mix. Iso-Seq primer bc01–12 will be added to each sample individually and should not be added to the master mix.														
<table border="1"><thead><tr><th>Components</th><th>Volume</th></tr></thead><tbody><tr><td>Iso-Seq cDNA PCR mix</td><td>25 µL</td></tr><tr><td>Iso-Seq cDNA amplification primer</td><td>2 µL</td></tr><tr><td>Total volume</td><td>27 µL</td></tr></tbody></table>		Components	Volume	Iso-Seq cDNA PCR mix	25 µL	Iso-Seq cDNA amplification primer	2 µL	Total volume	27 µL						
Components	Volume														
Iso-Seq cDNA PCR mix	25 µL														
Iso-Seq cDNA amplification primer	2 µL														
Total volume	27 µL														
3.1.2	On ice, add 27 µL of reaction mix 3 to the 21 µL of the eluted cDNA (from previous <a href="#">Section 2.4</a> ). Add 2 µL of Iso-Seq primer bc01–12 for a total volume of 50 µL.														
3.1.3	Thoroughly mix by pipetting up and down 10 times and then quick spin to collect all liquid.														
3.1.4	Run the thermal cycler program below with the lid temperature set to 105°C. Keep sample on ice until thermal cycler lid has heated to 105°C.														
<table border="1"><thead><tr><th colspan="2">PCR program</th></tr></thead><tbody><tr><td>45 seconds at 98°C</td><td>1 cycle</td></tr><tr><td>10 seconds at 98°C</td><td></td></tr><tr><td>15 seconds at 60°C</td><td>10 cycles</td></tr><tr><td>3 minutes at 72°C</td><td></td></tr><tr><td>5 minutes at 72°C</td><td></td></tr><tr><td>Hold at 4°C</td><td></td></tr></tbody></table>		PCR program		45 seconds at 98°C	1 cycle	10 seconds at 98°C		15 seconds at 60°C	10 cycles	3 minutes at 72°C		5 minutes at 72°C		Hold at 4°C	
PCR program															
45 seconds at 98°C	1 cycle														
10 seconds at 98°C															
15 seconds at 60°C	10 cycles														
3 minutes at 72°C															
5 minutes at 72°C															
Hold at 4°C															

SAFE STOPPING POINT -- Store at 4°C or -20°C for long-term storage


- Set up cDNA amplification reaction mix on ice
- Do not add Barcoded Forward Iso-Seq primers to master mix reaction

- Add desired Barcoded Forward Iso-Seq primer (select one of bc01-bc12) to each individual cDNA amplification reaction

- After setting up reactions on ice, add PCR reactions to thermal cycler after the lid has preheated to 105°C to avoid digestion of primers by polymerase exonuclease activity

# cDNA amplification + cleanup (cont.)

## Procedural notes



### 3.2 Cleanup of amplified cDNA using 0.9X SMRTbell Cleanup beads

Step	Instructions
3.2.1	Add 45 $\mu$ L (0.9x) of resuspended, room-temperature SMRTbell cleanup beads to the 50 $\mu$ L of cDNA amplified reaction from <a href="#">Section 3.1</a> . The correct ratio of beads to sample is critical at this step.
3.2.2	Mix beads by pipetting 10 times or until evenly distributed.
3.2.3	Quick-spin strip tubes in a microcentrifuge to collect liquid.
3.2.4	Leave at room temperature for 10 minutes to allow DNA to bind beads.
3.2.14	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer 24 $\mu$ L of the supernatant to a new strip tube. Discard the old strip tube with beads. <b>Recommended: Measure concentration and size distribution of each cDNA sample.</b>
3.2.15	<ul style="list-style-type: none"> <li>Take a 1 <math>\mu</math>L aliquot from each strip tube. Dilute each aliquot with 4 <math>\mu</math>L of elution buffer.</li> <li>Measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit.</li> <li>Dilute 1:4 dilution further to 1.5 ng/<math>\mu</math>L based on the Qubit reading if needed.</li> <li>Run 1 <math>\mu</math>L on an Agilent Bioanalyzer using a High Sensitivity DNA kit.</li> </ul>
3.2.16	The expected recovery after cDNA amplification SMRTbell clean-up is >100 ng. A minimum of 55 ng of total cDNA is recommended to proceed with Kinnex PCR ( <a href="#">Step 4</a> ). If less than 55 ng but more than 25 ng is recovered, proceed with Kinnex PCR but expect lower yields. Do not proceed with less than 25 ng.

• Perform **0.9X** SMRTbell bead cleanup<sup>1</sup>


• Perform DNA concentration QC using Qubit ds DNA HS assay and DNA sizing QC using Bioanalyzer

• Expected yield of purified cDNA product is **>100 ng**  
 • Minimum cDNA amount needed to proceed with Kinnex PCR = **55 ng** (do not proceed with <25 ng)

### 3.3 Pooling barcoded cDNA (skip if not multiplexing)

Step	Instructions
3.3.1	Using the concentration reading from the Qubit fluorometer, pool an equal mass of each barcoded cDNA sample for a total mass of 55 ng. Store any remaining purified, amplified barcoded cDNA at 4°C for future use.
3.3.2	Quick-spin the tube strip in a microcentrifuge to collect liquid.
3.3.3	Proceed to next step of the protocol.


Example Bioanalyzer DNA sizing QC analysis results for amplified full-length cDNA generated from a universal human RNA reference (UHRR) total RNA sample.




• If multiplexing, perform **equal-mass pooling** of each barcoded cDNA sample after 0.9X SMRTbell bead cleanup

# Kinnex PCR

In this step, incorporate programmable Kinnex segmentation adapter sequences into amplified cDNA products




Set up 8 parallel Kinnex PCR reactions per sample with premixed Kinnex primers to generate amplified cDNA products containing programmable sequences at both ends



Example: Amplified cDNA products from Kinnex PCR 3 contain flanking segmentation adapter sequences that are **complementary** to the ends of cDNA products from Kinnex PCR 2 & Kinnex PCR 4

# Kinnex PCR (cont.)

## Procedural notes



### 4.1 Kinnex PCR

Step	Instructions																																				
4.1.1	Thaw primers. The entire volume of primers can be transferred to an 8-tube strip for ease of use with a multi-channel pipette.																																				
	<table border="1"> <thead> <tr> <th></th> <th>8X concatenation</th> <th>Tube color</th> <th>P/N</th> </tr> </thead> <tbody> <tr><td>1</td><td>Kinnex primer mix A</td><td>Orange</td><td>103-107-800</td></tr> <tr><td>2</td><td>Kinnex primer mix B</td><td>Orange</td><td>103-107-900</td></tr> <tr><td>3</td><td>Kinnex primer mix C</td><td>Orange</td><td>103-108-000</td></tr> <tr><td>4</td><td>Kinnex primer mix D</td><td>Orange</td><td>103-108-100</td></tr> <tr><td>5</td><td>Kinnex primer mix E</td><td>Orange</td><td>103-108-200</td></tr> <tr><td>6</td><td>Kinnex primer mix F</td><td>Orange</td><td>103-108-300</td></tr> <tr><td>7</td><td>Kinnex primer mix G</td><td>Orange</td><td>103-108-400</td></tr> <tr><td>8</td><td>Kinnex primer mix HQ</td><td>Orange</td><td>103-108-500</td></tr> </tbody> </table>		8X concatenation	Tube color	P/N	1	Kinnex primer mix A	Orange	103-107-800	2	Kinnex primer mix B	Orange	103-107-900	3	Kinnex primer mix C	Orange	103-108-000	4	Kinnex primer mix D	Orange	103-108-100	5	Kinnex primer mix E	Orange	103-108-200	6	Kinnex primer mix F	Orange	103-108-300	7	Kinnex primer mix G	Orange	103-108-400	8	Kinnex primer mix HQ	Orange	103-108-500
	8X concatenation	Tube color	P/N																																		
1	Kinnex primer mix A	Orange	103-107-800																																		
2	Kinnex primer mix B	Orange	103-107-900																																		
3	Kinnex primer mix C	Orange	103-108-000																																		
4	Kinnex primer mix D	Orange	103-108-100																																		
5	Kinnex primer mix E	Orange	103-108-200																																		
6	Kinnex primer mix F	Orange	103-108-300																																		
7	Kinnex primer mix G	Orange	103-108-400																																		
8	Kinnex primer mix HQ	Orange	103-108-500																																		
4.1.2	Briefly vortex to mix, then quick-spin to collect liquid and place the primer mixes on ice.																																				
4.1.3	Thaw the following components, briefly vortex to mix, then quick-spin to collect liquid and place on ice. Add the components on ice in a LoBind tube.																																				
	<table border="1"> <thead> <tr> <th>Master mix components</th> <th>Volume for 8X concatenation*</th> <th></th> </tr> </thead> <tbody> <tr> <td>PCR-grade water</td> <td>88-X <math>\mu</math>L</td> <td rowspan="4">X = 55 (ng)/purified pooled DNA concentration from <a href="#">step 3.2.16</a> (single-plex) or <a href="#">step 3.2.3</a> (multiplex) *10% overage included</td> </tr> <tr> <td>Kinnex PCR mix (103-107-700)</td> <td>110 <math>\mu</math>L</td> </tr> <tr> <td>55 ng of amplified cDNA from <a href="#">Step 3.2.16</a></td> <td>X <math>\mu</math>L</td> </tr> <tr> <td>Total volume</td> <td>198 <math>\mu</math>L</td> </tr> </tbody> </table>	Master mix components	Volume for 8X concatenation*		PCR-grade water	88-X $\mu$ L	X = 55 (ng)/purified pooled DNA concentration from <a href="#">step 3.2.16</a> (single-plex) or <a href="#">step 3.2.3</a> (multiplex) *10% overage included	Kinnex PCR mix (103-107-700)	110 $\mu$ L	55 ng of amplified cDNA from <a href="#">Step 3.2.16</a>	X $\mu$ L	Total volume	198 $\mu$ L																								
Master mix components	Volume for 8X concatenation*																																				
PCR-grade water	88-X $\mu$ L	X = 55 (ng)/purified pooled DNA concentration from <a href="#">step 3.2.16</a> (single-plex) or <a href="#">step 3.2.3</a> (multiplex) *10% overage included																																			
Kinnex PCR mix (103-107-700)	110 $\mu$ L																																				
55 ng of amplified cDNA from <a href="#">Step 3.2.16</a>	X $\mu$ L																																				
Total volume	198 $\mu$ L																																				
4.1.4.	Aliquot 22.5 $\mu$ L of Master Mix 1 into each of the 8 PCR tubes (for 8X concatenation).																																				
4.1.5	Add 2.5 $\mu$ L of Kinnex primer mix into each of 8 PCR tubes from Step 4.4.																																				
	Set up the thermal cycler program listed below with the lid set to 105°C. Keep sample(s) on ice until the lid is heated to 105°C.																																				
	The duration of PCR is approximately 1 hour.																																				
4.1.6	<table border="1"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Duration</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Initial denaturation</td> <td>98°C</td> <td>3 min</td> <td>1</td> </tr> <tr> <td>Denaturation</td> <td>98°C</td> <td>20 s</td> <td rowspan="3">9</td> </tr> <tr> <td>Annealing</td> <td>68°C</td> <td>30 s</td> </tr> <tr> <td>Extension</td> <td>72°C</td> <td>4 min</td> </tr> </tbody> </table>	Step	Temperature	Duration	Cycles	Initial denaturation	98°C	3 min	1	Denaturation	98°C	20 s	9	Annealing	68°C	30 s	Extension	72°C	4 min																		
Step	Temperature	Duration	Cycles																																		
Initial denaturation	98°C	3 min	1																																		
Denaturation	98°C	20 s	9																																		
Annealing	68°C	30 s																																			
Extension	72°C	4 min																																			

• Can transfer entire volume of primers to PCR tubes for ease of use with multi-channel pipettes (8 primer mix tubes)




• Set up Kinnex PCR reactions **ON ICE**  
 • PCR polymerase 3'→5' exonuclease activity negatively impacts amplification yield if prepared at room temp.

• **Critical step! Correct setup of all 8 Kinnex PCR reactions is required** – any missing/incorrect MAS primer pairs will result in no/low SMRTbell yield

• Set up on ice and add PCR reaction to thermal cycler after lid has preheated to 105°C to avoid digestion of primers by polymerase exonuclease activity

# Kinnex PCR (cont.)

## Procedural notes



### 4.2 Pooling of 8 Kinnex PCR products and 1.05X SMRTbell cleanup

✓ Step	Instructions
4.2.1	Add 23 $\mu\text{L}$ from each of the 8 PCR reactions into a 1.5 mL tube for a total volume of 184 $\mu\text{L}$ . An equal volume of each PCR product is necessary for efficient array assembly.
4.2.2	Add 193 $\mu\text{L}$ (1.05X v/v) of resuspended, room-temperature SMRTbell cleanup beads to a tube of pooled Kinnex PCR amplicon. The correct ratio of beads to pooled sample is critical at this step.
4.2.3	Pipette-mix the beads until evenly distributed.
4.2.4	Quick-spin the tube in a microcentrifuge to collect liquid.
4.2.5	Leave at room temperature for 10 minutes to allow the DNA to bind beads
4.2.6	Place the tube in a magnetic separation rack until the beads separate fully from the solution.
4.2.7	Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
4.2.8	Slowly dispense 200 $\mu\text{L}$ , or enough to cover the beads, of freshly prepared 80% ethanol into the tube. After 30 seconds, pipette off the 80% ethanol and discard.
4.2.9	Repeat the previous step.
	Remove residual 80% ethanol:
4.2.10	<ul style="list-style-type: none"><li>Remove the tube from the magnetic separation rack.</li><li>Quick-spin the tube in a microcentrifuge.</li><li>Place the tube back in the magnetic separation rack until the beads separate fully from the solution.</li><li>Pipette off residual 80% ethanol and discard.</li></ul>
4.2.11	Remove the tube from the magnetic rack. Immediately add 40 $\mu\text{L}$ of elution buffer to the tube and resuspend the beads by pipetting 10 times or until evenly distributed.
4.2.12	Quick-spin the tube in a microcentrifuge to collect liquid.
4.2.13	Leave at room temperature for 5 minutes to elute DNA.
4.2.14	Place tube in a magnetic separation rack until beads separate fully from the solution.
4.2.15	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new LoBind tube. Discard old tube with beads.
4.2.16	Make a 1:10 dilution of the sample in elution buffer and measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit. Typical yield is 6–12 $\mu\text{g}$ .

- Pool exactly 23  $\mu\text{L}$  from each Kinnex PCR reaction in a clean 1.5 mL DNA LoBind tube<sup>1</sup> for a total combined volume of 184  $\mu\text{L}$

- Add exactly 193  $\mu\text{L}$  of SMRTbell cleanup beads (1.05X)
- Kinnex PCR mix significantly increases stringency of SMRTbell clean up beads, so accurate pipetting is critical

- Perform DNA concentration QC to verify there is sufficient yield of Kinnex PCR products (min. 4  $\mu\text{g}$ ) to proceed to Kinnex array formation step




# Kinnex array formation


In this step, assemble cDNA transcripts (“segments”) containing programmable ends into a linear array

- 1 Input RNA quality control
- ↓ Safe stop
- 2 cDNA synthesis + cleanup
- ↓
- 3 cDNA amplification + cleanup
- ↓ Safe stop
- 4 Kinnex PCR + cleanup
- ↓
- 5 Kinnex array formation + cleanup**
- ↓ Safe stop
- 6 Nuclease treatment + cleanup

Treat pooled Kinnex PCR products with Kinnex enzyme to create single-stranded overhangs to enable subsequent directional assembly of cDNA transcripts into a linear, ordered array



Barcoded Kinnex terminal adapters<sup>1</sup> are ligated to specific overhang sequences at array ends




Kinnex library prep workflow supports incorporation of PacBio barcodes at the array formation step to enable up to 4-plex sample multiplexing at the SMRTbell library level

# Kinnex array formation (cont.)

Kinnex terminal adapters incorporate barcode sequences to enable up to 4-plex sample multiplexing at the library level


- Kinnex adapters contain **barcode sequences<sup>1</sup>** to enable (optional) sample multiplexing at the SMRTbell library level (**up to 4-plex**)
  - Forward and reverse Kinnex adapter pairs are pre-mixed in Kinnex concatenation kits
  - Kinnex concatenation kits contain a total of **4 barcoded Kinnex adapter mixes (bcM0001-bcM0004)** to enable multiplexing of up to **4 samples per SMRT Cell**



# Kinnex array formation (cont.)

Kinnex terminal adapters use a new design that enables improved SMRT sequencing performance

- Kinnex adapters enable:
  - Longer polymerase read length → Improved HiFi conversion rate (HiFi reads/Total *P1* reads)
  - Improved *P1* loading efficiency




- New Kinnex adapter design requires a **different sequencing primer (Kinnex sequencing primer 103-179-000)**



# Kinnex array formation (cont.)

## Procedural notes



### 5.1 Kinnex array formation

Step	Instructions									
5.1.1	<p>In a 0.2 mL PCR tube, add 5 µg of sample from <a href="#">Step 4.2.15</a>, in 39 µL of volume (128 ng/µL). Dilute with elution buffer <b>going into this step</b> if the sample is too concentrated.</p> <p>Add 2 µL of Kinnex adapter bc01–04 (select one barcode per library preparation).</p> <p>Note: if combining multiple barcoded Kinnex libraries for sequencing, make sure each library uses one of the 4 different Kinnex barcoded adapters.</p> <table border="1"><thead><tr><th>Tube color</th><th>Components</th><th>Volume</th></tr></thead><tbody><tr><td>Blue</td><td>Kinnex adapter bc</td><td>2 µL</td></tr></tbody></table>	Tube color	Components	Volume	Blue	Kinnex adapter bc	2 µL			
Tube color	Components	Volume								
Blue	Kinnex adapter bc	2 µL								
5.1.2	<p>Add the following components in the listed order.</p> <p>If processing multiple samples, make a master mix with 10% overage. Pipette mix master mix.</p> <table border="1"><thead><tr><th>Tube color</th><th>Components</th><th>Volume</th></tr></thead><tbody><tr><td>Blue</td><td>Kinnex adapter bc</td><td>2 µL</td></tr></tbody></table>	Tube color	Components	Volume	Blue	Kinnex adapter bc	2 µL			
Tube color	Components	Volume								
Blue	Kinnex adapter bc	2 µL								
5.1.4	<p>After running the Kinnex primer digestion/ligation program, add 2 µL of DNA repair mix directly to the Kinnex primer digestion/ligation sample.</p>									
5.1.5	<p>Thoroughly mix by pipetting up and down 10 times and then quick spin to collect all liquid.</p> <p>Run the DNA Damage Repair Program with the lid set to &gt;55°C.</p>									
5.1.6	<table border="1"><thead><tr><th>Step</th><th>Temperature</th><th>Duration</th></tr></thead><tbody><tr><td>1</td><td>45°C</td><td>30 min</td></tr><tr><td>2</td><td>4°C</td><td>Hold</td></tr></tbody></table>	Step	Temperature	Duration	1	45°C	30 min	2	4°C	Hold
Step	Temperature	Duration								
1	45°C	30 min								
2	4°C	Hold								

- Recommended input amount to proceed with Kinnex array formation is **5 µg** of Kinnex PCR amplicons (from Step 4)
- Proceeding with <3 µg is **not recommended** since lower input amounts may lead to insufficient final library yields to enable optimal sequencing results

- **IMPORTANT:** If combining multiple barcoded Kinnex libraries for sequencing, make sure each library uses one of the 4 different Kinnex barcoded adapters

- Perform **DNA Damage Repair** step to repair nicked / damaged DNA sites within newly formed Kinnex array products

### 5.2 1X SMRTbell bead cleanup


Step	Instructions
5.2.1	<p>Add 1X v/v (60 µL) of resuspended, room temperature SMRTbell cleanup beads to each sample.</p>
5.2.2	<p>Pipette-mix the beads until evenly distributed and quick-spin in a microcentrifuge to collect liquid.</p>

- Perform **1X SMRTbell bead cleanup** at room temp.

# Nuclease treatment & cleanup

Perform nuclease treatment and final SMRTbell bead cleanup to remove incomplete SMRTbell template arrays


- 1 Input RNA quality control
- 2 cDNA synthesis + cleanup
- 3 cDNA amplification + cleanup
- 4 Kinnex PCR + cleanup
- 5 Kinnex array formation + cleanup
- 6 Nuclease treatment + cleanup



• After nuclease treatment, most remaining SMRTbell templates are complete (full-length) 8-segment transcript array molecules capped with Kinnex adapters

# Nuclease treatment & cleanup (cont.)

## Procedural notes



### 6.1 Nuclease treatment


Step	Instructions															
✓	Add the following components listed below to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus 10% overage. Pipette mix master mix. For individual preps, add components directly to each sample from the previous step in the order and volume listed below.															
6.1.1	<table border="1"><thead><tr><th colspan="3">Nuclease master mix</th></tr><tr><th>Tube</th><th>Component</th><th>Volume</th></tr></thead><tbody><tr><td>Light purple</td><td>Nuclease buffer (103-110-200)</td><td>5 µL</td></tr><tr><td>Light green</td><td>Nuclease mix (103-110-100)</td><td>5 µL</td></tr><tr><td colspan="2">Total volume</td><td>10 µL</td></tr></tbody></table>	Nuclease master mix			Tube	Component	Volume	Light purple	Nuclease buffer (103-110-200)	5 µL	Light green	Nuclease mix (103-110-100)	5 µL	Total volume		10 µL
Nuclease master mix																
Tube	Component	Volume														
Light purple	Nuclease buffer (103-110-200)	5 µL														
Light green	Nuclease mix (103-110-100)	5 µL														
Total volume		10 µL														
6.1.2	Add 10 µL of Nuclease Master mix to each sample. The total volume should be 50 µL.															
6.1.3	Thoroughly mix by pipetting up and down 10 times and then quick spin to collect all liquid.															
6.1.4	Run the nuclease treatment program with the lid set to >47°C. <table border="1"><thead><tr><th>Step</th><th>Temperature</th><th>Duration</th></tr></thead><tbody><tr><td>1</td><td>37°C</td><td>15 min</td></tr><tr><td>2</td><td>4°C</td><td>Hold</td></tr></tbody></table>	Step	Temperature	Duration	1	37°C	15 min	2	4°C	Hold						
Step	Temperature	Duration														
1	37°C	15 min														
2	4°C	Hold														

• Perform nuclease treatment for 15 min



# Nuclease treatment & cleanup (cont.)

## Procedural notes



Final Kinnex library yield is typically sufficient to load **≥2 SMRT Cells**

### 6.2 Final cleanup with 1X SMRTbell cleanup beads

Step	Instructions
6.2.1	Add 50 $\mu\text{L}$ SMRTbell cleanup beads to each sample from the previous step. Pipette-mix the beads until evenly distributed.
6.2.2	Quick-spin the tube strip in a microcentrifuge to collect all liquid.
6.2.3	Leave at room temperature for 10 minutes to allow DNA to bind the beads.
6.2.4	Place the tube strip in a magnetic separation rack until the beads separate fully from the solution.
6.2.5	Slowly pipette off the cleared supernatant without disturbing the beads and discard the supernatant.
6.2.6	Slowly dispense 200 $\mu\text{L}$ , or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, pipette off the 80% ethanol and discard.
6.2.7	Repeat the previous step.
Remove residual 80% ethanol:	
6.2.8	<ul style="list-style-type: none"> <li>Remove the tube strip from the magnetic separation rack.</li> <li>Quick-spin the tube strip in a microcentrifuge.</li> <li>Place the tube strip back in a magnetic separation rack until beads separate fully from</li> </ul>
6.2.13	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new PCR tube strip. Discard old tube strip with beads.
6.2.14	Take a 1 $\mu\text{L}$ aliquot from each tube. Make a 1:5 dilution of the sample in elution buffer and measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit. Calculate the total mass. Expect 10-25% recovery of the starting Kinnex-PCR product.
	Recommended: Further dilute each aliquot to 250 $\text{pg}/\mu\text{L}$ with the Femto Pulse dilution buffer. Measure the final SMRTbell library size distribution with a Femto Pulse system.
6.2.15	Proceed to SMRT Link Sample Setup to prepare the SMRTbell library for sequencing. DNA concentration must be less than 60 $\text{ng}/\mu\text{L}$ to go into ABC. Using a concentration above 60 $\text{ng}/\mu\text{L}$ will result in lower loading during sequencing.
6.2.16	Store SMRTbell libraries at 4°C if sequencing within the week. Long-term storage should be at -20°C. Minimize freeze-thaw cycles when handling SMRTbell libraries.

PROTOCOL COMPLETE



- Perform **DNA concentration QC** on final purified Kinnex RNA library using a Qubit dsDNA HS assay
  - Typical final SMRTbell library yield from 5  $\mu\text{g}$  of input DNA into Kinnex array formation is **~10 – 25%** – a much higher observed yield might suggest incomplete digestion of partial SMRTbell templates
  - Troubleshooting tip:** If SMRTbell library yield is higher than expected and *P1* loading is lower than expected, consider repeating the nuclease treatment step



- Perform **DNA sizing QC** on final purified Kinnex full-length RNA library using a Femto Pulse system (expected final library insert size is **~12 – 16 kb**)



- Kinnex full-length RNA final SMRTbell library concentration must be **≤60  $\text{ng}/\mu\text{L}$**  to proceed with SMRT Link sample setup (ABC<sup>1</sup>)
  - Using a concentration above 60  $\text{ng}/\mu\text{L}$  will result in lower loading during sequencing

# Kinnex full-length RNA library prep inputs & expected step yields

Final Kinnex library yield is typically sufficient to load  $\geq 2$  SMRT Cells

Kinnex FL RNA

Preparing Kinnex libraries using the Kinnex full-length RNA kit (103-238-700)



Example Femto pulse DNA sizing QC analysis results for final Kinnex full-length RNA library prepared with human universal human reference RNA (UHRR) total RNA sample.

## Example Kinnex full-length RNA library prep yields

Total RNA input for cDNA synthesis	300 ng
cDNA input for Kinnex array formation	5900 ng
Post-nuclease treatment & final library cleanup yield (%) <sup>1</sup>	1460 ng (24.7%)

<sup>1</sup> Post-nuclease treatment & final cleanup yields typically ranged from ~10% to ~25% when using UHRR total RNA samples for Kinnex full-length RNA library construction.


Final Kinnex library yield is typically sufficient to load  $\geq 2$  SMRT Cells



# Kinnex full-length RNA sequencing preparation workflow details

# Sample Setup & Run Design recommendations for Kinnex full-length RNA libraries

SMRT Link supports Kinnex full-length RNA sequencing preparation & analysis workflow for PacBio systems<sup>1</sup>



SMRT Link module	Key setup parameters For Kinnex libraries	Sequel II/IIe system recommended settings for Kinnex libraries	Revio system recommended settings for Kinnex libraries
Sample setup	Library type	Kinnex	
	Primer	Kinnex sequencing primer	
	Binding/Polymerase kit <sup>1</sup>	Sequel II binding kit 3.2 (includes Kinnex sequencing primer)	Revio polymerase kit (includes Kinnex sequencing primer)
	Concentration on plate	40 – 60 pM	100 – 150 pM
Runs → Run design	Adapter / Library type	SMRTbell Adapter Design = SMRTbell Kinnex Prep Kit	Library type = Kinnex
	Movie collection time	30 hrs	24 hrs
	Use adaptive loading	YES	
	On-instrument CCS	CCS Analysis Output - Include Low Quality Reads = NO CCS Analysis Output - Include Kinetics Information = NO	Consensus Mode = MOLECULE

<sup>1</sup> Kinnex full-length RNA kit requires SMRT Link v13.0 or higher.

# SMRT Link Sample Setup and Run Design for Kinnex kits video demonstration


Video demonstration of SMRT Link Sample Setup and Run Design setup procedure for Kinnex kits supporting full-length RNA sequencing, single-cell RNA sequencing and full-length 16S rRNA sequencing

Sample Setup / Sample Calculation  
Sequel II binding kit 3.1/3.2, Revio polymerase kit

Conversion Calculator  
Autosaved at 2023-11-20, 09:23:31 AM

+ Add Sample Group

< Sample group >	
	Copy Remove Lock Download CSV
Name	My Batch of Samples
Application	Kinnex full-length RNA
Library type	Kinnex
Polymerase / Binding kit	Revio polymerase kit
Number of samples	1 samples
SMRT Cells per sample	1 cells
Available volume per sample	20 uL
Insert size	16000 bp
Sample concentration	40 ng/uL
Cleanup anticipated yield	75 %
Concentration on plate	130 pM Recommended: 100-150 pM
Minimum pipetting volume	1 uL
Comment	

 YouTube

## [Demo video for Sample Setup and Run Design for Kinnex kits \(SMRT Link v13.0+\)](#)

- Demo video for Sample Setup and Run Design for Kinnex kits in SMRT Link v13.0+
- Kinnex kits support full-length RNA sequencing (Kinnex full-length RNA kit), full-length 16S rRNA sequencing (Kinnex 16S rRNA kit) and full-length single-cell RNA sequencing (Kinnex single-cell RNA kit)

# SMRT Link Sample Setup procedure for Kinnex full-length RNA libraries



## Revo system



## Sequel II and Ile systems

< Sample group >		< Sample group >	
<input type="button" value="Copy"/> <input type="button" value="Remove"/> <input type="button" value="Lock"/> <input type="button" value="Download CSV"/>		<input type="button" value="Copy"/> <input type="button" value="Remove"/> <input type="button" value="Lock"/> <input type="button" value="Download CSV"/>	
Name	Kinnex full-length RNA library demo	Name	Kinnex full-length RNA library demo
Application	Kinnex full-length RNA	Application	Kinnex full-length RNA
Library type	Kinnex	Library type	Kinnex
Polymerase / Binding kit	Revo polymerase kit	Polymerase / Binding kit	Sequel II Binding Kit 3.2
Number of samples	1 samples	Number of samples	1 samples
SMRT Cells per sample	1 cells	SMRT Cells per sample	1 cells
Available volume per sample	20 uL	Available volume per sample	20 uL
Insert size	16000 bp	Insert size	16000 bp
Sample concentration	40 ng/uL	Sample concentration	40 ng/uL
Cleanup anticipated yield	75 %	Cleanup anticipated yield	75 %
Concentration on plate	130 pM Recommended: 100-150 pM	Concentration on plate	50 pM Recommended: 40-60 pM
Minimum pipetting volume	1 uL	Minimum pipetting volume	1 uL
Comment	Kinnex library containing array of 8 FL cDNA segments	Comment	Kinnex library containing array of 8 FL cDNA segments

- Select **application type** to autofill fields in green

### IMPORTANT: Specify Library type = Kinnex

- Library type field determines sequencing primer type to use for annealing step
  - Kinnex libraries require use of **Kinnex sequencing primer**<sup>1</sup>

- Select **Revo polymerase kit** for Revo system and **Sequel II Binding Kit 3.2** for Sequel II/Ile systems

- Recommended Kinnex full-length RNA library input concentration for sample setup is **20 – 60 ng/μL**

- Recommended OPLC range is **100 – 150 pM** for Revo system and **40 – 60 pM** for Sequel II/Ile systems

### Recommended target P1 loading range

- Revo system: **~50 – 70%**
- Sequel II and Ile systems: **~60 – 80%**



# SMRT Link Run Design procedure for **Revio system**

## Sample and run information

**Kinnex full-length RNA**

▼ **Plate 1, Well A01: Kinnex full-length RNA library demo**

Application Required	Kinnex full-length RNA
Plate Well Required	Plate 1, Well A01
Well Name Required	Kinnex full-length RNA library demo
Well Comment	
Library Type Required	Kinnex
Insert Size (bp) Required	16000
Polymerase Kit Required	Revio polymerase kit
Movie Acquisition Time (hours)	24

**Use Adaptive Loading**  
 YES  NO



Forward and reverse standard terminal adapters have the same structure



Forward and reverse Kinnex terminal adapters have different structures

**Library Type** field determines which adapter finding algorithm is used during post-primary analysis<sup>1</sup>

<sup>1</sup> **Note:** When sequencing a Kinnex library sample, if 'Standard' library type is mistakenly selected instead of 'Kinnex' then a higher missing adapter rate (> 95%) and a slight degradation in barcode demultiplexing performance (~93-96% barcoded HiFi read yield) will be observed.

# SMRT Link Run Design procedure for **Revio system** (cont.)

## Sample indexing (barcoding) information



### Kinnex full-length RNA

Default = YES for Sample is indexed

Samples

Sample is indexed  YES  NO

Indexes Required MAS SMRTbell barcoded adapters (v2)


Same Barcodes on Both Ends of Sequence  YES  NO

Biosample names Required Interactively From a File

Specify Indexes FASTA = MAS SMRTbell barcoded adapters (v2)

Specify YES for Same barcodes on both ends of sequences

Example complete Kinnex library molecule containing barcoded Kinnex terminal MAS adapters<sup>1</sup> at both ends



Example interactive biosample name specification for a multiplexed Kinnex library sample

Barcode Selector and Sample Name Editor

Barcode ID	Bio Sample ID
<input type="checkbox"/> bcM0001--bcM0001	
<input type="checkbox"/> bcM0002--bcM0002	
<input type="checkbox"/> bcM0003--bcM0003	
<input type="checkbox"/> bcM0004--bcM0004	

Barcode Selector and Sample Name Editor

Barcode ID	Bio Sample ID
<input type="checkbox"/> bcM0003--bcM0003	
<input type="checkbox"/> bcM0004--bcM0004	
<input type="checkbox"/> bcM0001--bcM0001	Kinnex adapter-barcoded library 1
<input type="checkbox"/> bcM0002--bcM0002	Kinnex adapter-barcoded library 2

SMRT Link



Data Management

MAS SMRTbell barcoded adapter indexes

```
>bcM0001
ACAGTC
>bcM0002
ATGACG
>bcM0003
CACGTG
>bcM0004
CATCGC
```

# SMRT Link Run Design procedure for **Revio system** (cont.)

## Run options and data options



### **Kinnex full-length RNA**

▼ **Run Options**

Library Concentration (pM)   
Required

On-plate loading concentration is required for Revio samples

▼ **Data Options**

Include Base Kinetics  YES  NO

Consensus Mode  MOLECULE  STRAND

Assign Data To Project ⓘ

Default = NO for Include Base Kinetics

Default Consensus Mode = MOLECULE<sup>1</sup>

Can leave Include Base Kinetics and Consensus Mode fields at their default settings for Kinnex library samples

# SMRT Link Run Design procedure for Sequel II/Ile systems

## Sample information and run information



### Kinnex full-length RNA

- Select desired **Kinnex application** from the **Application** field drop-down menu
- The following fields are **auto-populated** with default recommended values and high-lighted in **green**:

**SMRTbell Adapter Design**

→ SMRTbell Kinnex Prep Kit

**Binding Kit**

→ Sequel II Binding Kit 3.2

**Sequencing Kit**

→ Sequel II Sequencing Plate 2.0 (4 rxn or 1 rxn)

**DNA Control Complex**

→ Sequel II DNA Internal Control Complex 3.2

**Movie Time per SMRT Cell**

→ 30 hrs

**Pre-Extension Time**

→ 2 hrs

**SMRTbell Adapter Design** field determines which adapter finding algorithm is used during post-primary analysis<sup>1</sup>

Default SMRTbell adapter design for Kinnex samples is **SMRTbell Kinnex Prep Kit**

Recommended OPLC for Sequel II/Ile Kinnex library samples is **40 – 60 pM**

Recommended movie time = **30 hrs**

Select desired Kinnex application type from drop-down menu

SAMPLE 1: Kinnex full-length RNA library demo , A01, 30 hour movie, 16000 bp insert

Copy Delete

Import from Sample Setup Select Sample

Application Required Kinnex full-length RNA

Well Sample Name Required Kinnex full-length RNA library demo

Bio Sample Name

Sample Comment

Sample Well A01

SMRTbell Adapter Design Required SMRTbell® Kinnex Prep Kit

Binding Kit Required Sequel® II Binding Kit 3.2

Sequencing Kit Required Sequel® II Sequencing Plate 2.0 (4 rxn)

DNA Control Complex Sequel® II DNA Internal Control Complex 3.2

Insert Size (bp) Required 16000

Recommended Concentration on Plate (pM) 40 – 60 pM

On-Plate Loading Concentration (pM) Required 45

Movie Time per SMRT Cell (hours) 30

Use Pre-Extension YES NO

Pre-Extension Time (hours) 2

CCS Analysis will be performed on-instrument to produce HiFi .bam files.

Example sample information entered into a Sequel Ile system run design worksheet for a Kinnex full-length RNA library sample.

# SMRT Link Run Design procedure for Sequel II/Ile systems (cont.)

## Advanced options



### Kinnex full-length RNA

- For all Kinnex library samples, leave the following **Advanced Options** fields at their **default settings**
  - Use Adaptive Loading**  
→ YES
  - Loading Target (P1 + P2)**  
→ 0.85
  - Maximum Loading Time**  
→ 2 hours
  - CCS Analysis Output - Include Low Quality Reads**  
→ NO
  - CCS Analysis Output - Include Kinetics Information**  
→ NO
  - Pre-Extension Time**  
→ 2 hrs
- If desired, specify to use an alternative project folder for the **Add Data to Project** field

**Advanced Options**

Use Adaptive Loading  YES  NO

Loading Target (P1 + P2)

Maximum Loading Time (hours)

CCS Analysis Output - Include Low Quality Reads  YES  NO

CCS Analysis Output - Include Kinetics Information  YES  NO

Add Data to Project

Leave these Advanced Options fields at their **default values**

Can specify to use a different Project folder

Example default Advanced Options settings entered into a Sequel IIe system run design worksheet for a Kinnex full-length RNA library sample.

# SMRT Link Run Design procedure for Sequel II/Ile systems (cont.)

## Barcoded sample options



### Kinnex full-length RNA

- For Kinnex library samples, can leave most **Barcoded Sample Options** fields at their **default settings**

Specify Bio Sample Names, either interactively or by downloading a CSV file (**Interactively** or **From a file**)

If desired, specify to perform barcode demultiplexing on-instrument or in SMRT Link (default = **On-instrument** for Sequel IIe system)

Barcoded Sample Options

Sample Is Barcoded  YES  NO

Barcode Set Required MAS SMRTbell barcoded adapters (v2)

Same Barcodes on Both Ends of Sequence i  YES  NO

Assign Bio Sample Names to Barcodes i Required


Demultiplex Barcodes  ON INSTRUMENT  IN SMRT LINK  
 DO NOT GENERATE

Can leave most of these fields at their default values

Example default Barcoded Sample Options settings entered into a Sequel IIe system run design worksheet for a Kinnex full-length RNA library sample.



# SMRT Link Run Design analysis options for **Revio system** and **Sequel II/IIe systems**



**Kinnex full-length RNA**

**Analysis Options** Default = YES for Add Analysis

Add Analysis  YES  NO

Analysis Name Required Kinnex\_Full-Length\_RNA\_Demo\_Analysis\_Job\_Name

Select Analysis Workflow Required Read Segmentation and Iso-Seq

Segmentation Adapter Set MAS-Seq Adapter v3 (MAS8)

Primer Set Required Iso-Seq v2 Barcoded cDNA Primers


Reference Set Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples ? Pool reads and cluster together

Advanced Parameters

Analysis Workflow is automatically filled in (Default = Read Segmentation and Iso-Seq)

## Read Segmentation and Iso-Seq



Perform isoform-classification analysis to **identify novel genes & isoforms with abundance information** (raw counts & normalized counts per million)

# SMRT Link Run Design analysis options for **Revio system** and **Sequel II/IIe systems** (cont.)



## Kinnex full-length RNA

**Analysis Options**

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required

Segmentation Adapter Set


Primer Set Required

Reference Set

Cluster of Barcoded Samples

[Advanced Parameters](#)

Specify **Segmentation Adapter Set** that corresponds to the Kinnex library concatenation method used  
→ For Kinnex full-length RNA samples, specify **MAS-Seq Adapter v3 (MAS8)**



# SMRT Link Run Design analysis options for **Revio system** and **Sequel II/IIe systems** (cont.)



## Kinnex full-length RNA

Specify **Primer Set** used for cDNA amplification

Specify primer sequence file in FASTA format to **identify cDNA primers for removal** (include the 5' and 3' cDNA primers)



For Kinnex full-length RNA analysis, select 'Iso-Seq v2 Barcoded cDNA primers'

**Analysis Options**

Add Analysis  YES  NO

Analysis Name Required: Kinnex\_Full-Length\_RNA\_Demo\_Analysis\_Job\_Name

Select Analysis Workflow Required: Read Segmentation and Iso-Seq

Segmentation Adapter Set: MAS-Seq Adapter v3 (MAS8)


Primer Set Required: Iso-Seq v2 Barcoded cDNA Primers

Reference Set: Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples ⓘ: Pool reads and cluster together

Advanced Parameters

Specify reference genome & annotation sets to **align high quality isoforms to, and to collapse isoforms mapped to the same genomic loci.**



Specify Reference Set. Default sets are:

- Human Genome hg38, with Gencode v39 annotations
- Mouse Genome mm39, with Gencode vM28 annotations

```
>IsoSeqX_bc01_5p
CTACACGACGCTCTCCGATCTACTACACGCAATGAAGTCGCAGGGTTGGG
>IsoSeqX_bc02_5p
CTACACGACGCTCTCCGATCTACTAGTAGCAATGAAGTCGCAGGGTTGGG
>IsoSeqX_bc03_5p
CTACACGACGCTCTCCGATCTAGTGACGCAATGAAGTCGCAGGGTTGGG
⋮
>IsoSeqX_bc11_5p
CTACACGACGCTCTCCGATCTGTAGTAGCAATGAAGTCGCAGGGTTGGG
>IsoSeqX_bc12_5p
CTACACGACGCTCTCCGATCTGTATGACGCAATGAAGTCGCAGGGTTGGG
>IsoSeqX_3p
AAGCAGTGGTATCAACGCAGAGTAC
```

# SMRT Link Run Design analysis options for **Revio system** and **Sequel II/IIe systems** (cont.)



**Analysis Options**

Add Analysis  YES  NO

Analysis Name Required: Kinnex\_Full-Length\_RNA\_Demo\_Analysis\_Job\_Name

Select Analysis Workflow Required: Read Segmentation and Iso-Seq

Segmentation Adapter Set: MAS-Seq Adapter v3 (MAS8)

Primer Set Required: Iso-Seq v2 Barcoded cDNA Primers

Reference Set: Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples  Pool reads and cluster together

Advanced Parameters

Specify how to perform read clustering for barcoded samples

Pool reads and cluster together





-- select --

Cluster reads separately

Pool reads and cluster together

Specification of **Cluster of Barcoded Samples** setting determines whether all FLNC reads will be pooled for clustering (Does not apply to non-barcoded samples.)

## Clustering options for Kinnex full-length RNA samples



Analyze pooled samples (default)

Analyze by sample






# Kinnex full-length RNA example sequencing performance data

# Example Kinnex full-length RNA library preparation QC results


Kinnex full-length RNA library prepared with human UHRR total RNA

## Amplified full-length cDNA QC



Example Bioanalyzer DNA sizing QC analysis results for amplified full-length cDNA generated from a universal human RNA reference (UHRR) total RNA sample.

## Final Kinnex full-length RNA library QC



Example Femto Pulse DNA sizing QC analysis results for final Kinnex full-length RNA library.

Final Kinnex library yield is typically sufficient to load  $\geq 2$  SMRT Cells

Total RNA input for cDNA synthesis	300 ng
cDNA input for Kinnex array formation	5900 ng
Post-nuclease treatment & final library cleanup yield (%) <sup>1</sup>	1460 ng (24.7%)


<sup>1</sup> Post-nuclease treatment & final cleanup yields typically ranged from ~10% to ~25% when using UHRR total RNA samples for Kinnex full-length RNA library construction.




# Example sequencing performance for Kinnex full-length RNA libraries prepared with human cDNA

Sequel IIe system example data<sup>1</sup>


## Raw Data Report



## HiFi Read Length



## Read Segmentation Metrics



**For Sequel II/IIe system Kinnex FL RNA samples**

- Yield of HiFi reads is ideally >1.5 Million
- Yield of segmented reads is ideally >15 Million
- Mean S-read length will depend on size distribution of cDNA sample (ideally >1,500 bp)
- % of reads with full arrays is ideally >80%
- Mean array size is ideally >7.0 segments

Raw Base Yield	652 Gb
Mean Polymerase Read Length	115.36 kb
P0	28%
P1	71%
P2	1%

Example sequencing metrics for a Universal Human Reference RNA (UHRR) Kinnex full-length RNA library sample run on a Sequel IIe system with Binding Kit 3.2 (Polymerase 2.2) / 80 pM on-plate loading concentration (OPLC) / 30-hrs movie time / 2-hrs pre-extension time.

HiFi Reads	2.3 M
HiFi Base Yield	43.2 Gb
Mean HiFi Read Length	19.1 kb
Median HiFi Read Quality	Q32
HiFi Read Mean # of Passes	10

For UHRR Kinnex full-length RNA libraries, per-SMRT Cell 8M HiFi read counts typically ranged from ~2 – 3 Million depending on the final library insert size.


Reads	2,260,039
Segmented reads (S-reads)	17,213,165
Mean length of S-reads	2,420 bp
Percent of reads with full arrays	91.07%
Mean array size (concentration factor)	7.62

For UHRR Kinnex libraries, per-SMRT Cell 8M segmentation read counts were typically ~15 – 20 Million.


# Example sequencing performance for Kinnex full-length RNA libraries prepared with human cDNA

Revio system example data<sup>1</sup>


## Raw Data Report



## HiFi Read Length



## Read Segmentation Metrics



**For Revio system Kinnex FL RNA samples**

- Yield of HiFi reads is ideally >4.5 Million
- Yield of segmented reads is ideally >30 Million
- Mean S-read length will depend on size distribution of cDNA sample (ideally >1,500 bp)
- % of reads with full arrays is ideally >80%
- Mean array size is ideally >7.0 segments

Raw Base Yield	1,168 Gb
Mean Polymerase Read Length	58.6 kb
P0	16%
P1	79%
P2	5%

Example sequencing metrics for a Universal Human Reference RNA (UHRR) Kinnex full-length RNA library sample run on a Revio system with Revio polymerase kit / 130 pM on-plate loading concentration (OPLC) / 24-hrs movie time.

HiFi Reads	5.1 M
HiFi Base Yield	93.47 Gb
Mean HiFi Read Length	18.41 kb
Median HiFi Read Quality	Q28
HiFi Read Mean # of Passes	7

For UHRR Kinnex full-length RNA libraries, per-Revio SMRT Cell HiFi read counts were typically ~5 – 6 Million depending on the final library insert size and P1 loading performance.

Input HiFi Reads	5,027,154
Segmented reads (S-reads)	37,216,151
Mean length of S-reads	2,393 bp
Percent of reads with full arrays	85.84%
Mean array size (concentration factor)	7.40


For UHRR Kinnex libraries, per-Revio SMRT Cell segmentation read counts were typically ~30 – 45 Million.



# **Kinnex full-length RNA data analysis workflow overview**

# Kinnex full-length RNA bioinformatics workflow overview

SMRT Link Read segmentation and Iso-Seq workflow processes HiFi reads generated from Kinnex full-length RNA libraries to produce classified isoforms with read counts that are compatible with tertiary analysis tools



## Read segmentation

- HiFi reads are segmented into individual segmented reads (**S-reads**) that represent the original cDNA sequences

## Primer extraction

- Primers and polyA tails are removed, but also used to orient the read into 5' → 3' orientation

## Isoform clustering

- FLNC reads are clustered by their sequencing similarity to produce isoform consensus sequences
- This step is the last step of Iso-Seq analysis if no genome is provided

## Mapping

- If a genome is provided, isoform consensus sequences from the previous step are mapped and further collapsed by their exonic structures to produce isoforms as GFF files for visualization

## Transcript classification<sup>2</sup>

- If an annotation (e.g., Gencode) is provided, isoforms are classified against it using pigeon (the PacBio implementation of SQANTI3) to identify known and novel genes/isoforms
- The Iso-Seq workflow can jointly analyze pooled sample reads to produce a unified isoform annotation with per-sample read counts, both raw and normalized as counts per million (CPM)

# SMRT Link Read Segmentation and Iso-Seq analysis video demonstration

Video demonstration of SMRT Link Read Segmentation and Iso-Seq application workflow for analysis of Kinnex full-length RNA samples

Analysis Application Required

Read Segmentation and Iso-Seq

Import Analysis Settings Export

Associated Inputs

Segmentation Adapter Set

MAS-Seq Adapter v3 (MAS8)

Primer Set Required

Iso-Seq v2 Barcoded cDNA Primers

Reference Set

Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples

Pool reads and cluster together

Advanced Parameters

Analysis Name

test


Analysis Datasets

Displaying rows 1 to 1 out of 1

ID	Name
21...	3230211_KPoS_64007_...

[Demo video for Read Segmentation and Iso-Seq workflow \(SMRT Link v13.0+\)](#)

- Workflow supports full-length isoform analysis for data generated on PacBio Sequel II/IIe and Revio systems using **Kinnex full-length RNA kit**
- End-to-end workflow begins with HiFi reads and **outputs full-length isoform classifications with supporting read count information**



# Kinnex full-length RNA bioinformatics workflow recommendations

SMRT Link Read Segmentation and Iso-Seq workflow common considerations and recommendations for analysis of Kinnex full-length RNA data<sup>1</sup>

## Analysis recommendations for Iso-Seq data based on reference genome and annotation availability<sup>1</sup>

- With SMRT Link v13, the **Read segmentation and Iso-Seq workflow** analysis application supports **human and mouse** reference genomes and annotations to produce classified isoforms with read counts.
- If working with **other organisms**, see table below for analysis recommendations

Available reference or annotation	Analysis workflow recommendation
Human or mouse	<ul style="list-style-type: none"><li>• Use the Iso-Seq workflow with preloaded human / mouse annotation to get mapped, unique isoforms with classifications and read count information (FASTA, GFF, TXT).</li></ul>
Model organism with good annotation	<ul style="list-style-type: none"><li>• Run Iso-Seq workflow with uploaded reference genome to get mapped, unique isoforms (FASTA, GFF)</li><li>• Generate pigeon-compliant annotation and use the command line for isoform classification with read count information (TXT)</li></ul>
Non-model organism with genome	<ul style="list-style-type: none"><li>• Run Iso-Seq workflow with uploaded reference genome to get mapped, unique isoforms (FASTA, GFF)</li></ul>
No genome	<ul style="list-style-type: none"><li>• Run Iso-Seq workflow without reference genome to get unique isoforms (FASTA)</li></ul>



# Kinnex full-length RNA bioinformatics workflow recommendations (cont.)

SMRT Link Read Segmentation and Iso-Seq workflow common considerations and recommendations for analysis of Kinnex full-length RNA data<sup>1</sup>

## Sequencing depth recommendations for Iso-Seq data based on experimental goals and study design

Example application	Human genetics disease studies	Biopharma for identifying highly expressed targets	Plant & animal whole genome annotation
Experimental goal	Isoform discovery and quantification of moderate-to-rare transcripts	Isoform discovery of high expressed transcripts	Comprehensive transcript annotation in a species
Example study design	Disease vs. normal tissues with multiple replicates	Disease cohort with >20+ samples	Plant or animal with multiple tissue types
Target depth of coverage per sample	10 M reads per sample	5 M reads per sample	≤5 M reads per tissue (of same species)
Sample multiplexing <sup>1</sup>	<b>Sequel II/IIe system:</b> Up to 2 samples per SMRT Cell 8M (2-plex)	<b>Sequel II/IIe system:</b> Up to 3 samples per SMRT Cell 8M (3-plex)	<b>Sequel II/IIe system:</b> Up to 3 tissue types per SMRT Cell 8M (3-plex)
	<b>Revio system:</b> Up to 4 samples per Revio SMRT Cell (4-plex)	<b>Revio system:</b> Up to 8 samples per Revio SMRT Cell (8-plex)	<b>Revio system:</b> Up to 8 tissue types per Revio SMRT Cell (8-plex) <sup>2</sup>
SMRT Link data analysis workflows	Read Segmentation and Iso-Seq analysis application with option to “pool reads and cluster together” to get a master isoform classification file with per-sample full-length read counts		



# SMRT Link Read Segmentation and Iso-Seq analysis application setup

Specify **Read Segmentation and Iso-Seq** analysis application in SMRT Link<sup>1</sup>

• **Read Segmentation and Iso-Seq** analysis application processes HiFi reads generated from a Kinnex full-length RNA library to produce classified isoforms with read counts that are compatible with tertiary analysis tools

- Accepts **HiFi reads** (BAM format) as input
  - HiFi reads are reads generated with CCS analysis whose quality value is equal to or greater than 20
- HiFi reads should be generated using the Kinnex full-length RNA library preparation protocol ([103-238-700](#))
  - If the library is a regular (non-Kinnex) Iso-Seq monomer library without MAS-Seq concatenation, use the SMRT Link Iso-Seq Analysis workflow instead

**Note on barcoded libraries**

The **Read Segmentation and Iso-Seq** workflow will only process barcoded libraries at the cDNA level (such as using **Iso-Seq v2 Barcoded cDNA Primers** as part of the MAS-Seq for bulk Iso-Seq kit)

→ Demultiplexing of **barcoded adapters** (also part of the Kinnex full-length RNA kit) should first be performed by running the **Demultiplex Barcodes** data utility workflow in SMRT Link.

# SMRT Link Read Segmentation and Iso-Seq analysis application setup (cont.)

Specify Read Segmentation and Single-Cell Iso-Seq analysis application required associated inputs<sup>1</sup>

PacBio SMRT Analysis

SMRT Analysis / Create New Analysis

1. Select Data 2. Select Analysis

Analysis Application Required

Read Segmentation and Iso-Seq

Import Analysis Settings Export

Associated Inputs

1 Segmentation Adapter Set  
MAS-Seq Adapter v3 (MAS8)

2 Primer Set Required  
Iso-Seq v2 Barcoded cDNA Primers

3 Reference Set  
Human Genome hg38, with Gencode v39 annotatic

Cluster of Barcoded Samples  
Pool reads and cluster together

Advanced Parameters

## 1. Segmentation Adapter Set (Default = MAS-Seq Adapter v3 (MAS8))

- Specify a FASTA file, provided by PacBio, containing segmentation adapters. If you need a custom segmentation adapter set, click Advanced Parameters and use a custom FASTA file formatted as described in the SMRT Link User Guide [documentation](#).

## 2. Primer Set (Required) (Default = Iso-Seq v2 Barcoded cDNA Primers)

- Specify a primer sequence file in FASTA format to identify cDNA primers for removal. The primer sequence includes the 5' and 3' cDNA primers
- Primer IDs must be specified using the suffix `_5p` to indicate 5' cDNA primers and the suffix `_3p` to indicate 3' cDNA primers. The 3' cDNA primer should not include the `Ts` and is written in reverse complement. (See the [SMRT Link User Guide](#) for example Iso-Seq v2 Barcoded cDNA Primer IDs and sequences)
- Each primer sequence must be unique

## 3. Reference Set (Required)

- Specify one of two default reference genome and annotation sets to align high quality isoforms to, and to collapse isoforms mapped to the same genomic loci. The default sets are `Human_hg38_Gencode_v39` and `Mouse_mm39_Gencode_vM28` annotations
- Alternatively, choose other reference genomes (but not with annotations) that were custom-uploaded to SMRT Link
- The Reference Set can be left blank. If blank, the workflow will stop after the isoform clustering step (`isoseq cluster`)

# SMRT Link Read Segmentation and Iso-Seq analysis application setup (cont.)

Specify **Read Segmentation** and **Single-Cell Iso-Seq** analysis application required associated inputs<sup>1</sup>

**PacBio** SMRT Analysis

SMRT Analysis / Create New Analysis

1. Select Data 2. Select Analysis

Analysis Application Required

Read Segmentation and Iso-Seq

Import Analysis Settings Export

Associated Inputs

Segmentation Adapter Set

MAS-Seq Adapter v3 (MAS8)

Primer Set Required

Iso-Seq v2 Barcoded cDNA Primers

Reference Set

Human Genome hg38, with Gencode v39 annotatic


**4** Cluster of Barcoded Samples

Pool reads and cluster together

Advanced Parameters

## 4. Cluster of Barcoded Samples

- This option specifies barcoded samples that were barcoded at the cDNA level, where the (barcoded) cDNA primers are specified in the Primer Set option. This option does **not** address libraries that were barcoded using barcoded adapters
- Specify whether all FLNC reads will be pooled for clustering, then demultiplexed based on pooled result. **Note:** This setting does **not** apply to non-barcoded samples
- Specify **Pool reads and cluster together** if barcoded samples are from the **same** species, but different tissues, or samples of the same genes but different individuals. The samples are clustered with **all** barcodes pooled
- Specify **Cluster reads separately** if barcoded samples are from different species. The samples are clustered **separately** by barcode
- In either case, the samples on the results page are automatically named BioSample\_1 through BioSample\_N



# Example SMRT Link Read Segmentation data utility processing results<sup>1</sup> for Kinnex full-length RNA library prepared with human UHRR sample

## SMRT Link Read Segmentation data utility job report – Summary Metrics and Segmentation Statistics

### Summary Metrics


Value	Analysis Metric
6,400,421	Reads
48,911,387	Segmented reads (S-reads)
2,111	Mean length of S-reads
91.89 %	Percent of reads with full arrays
7.64	Mean array size (concatenation factor)

Example Revo system data shown.

**For Revo system with optimal sample P1 loading:**

- Yield of HiFi reads is ideally >4.5 Million
- Yield of segmented reads is ideally >30 Million
- Mean S-read length will depend on size distribution of cDNA sample (ideally >1,500 bp)
- % of reads with full arrays is ideally >80%
- Mean array size is ideally >7.0 segments

### Segmentation Statistics



• Plot should ideally show a majority (>80%) of reads with array lengths = 8 segments

- **Reads:** Number of input arrayed HiFi reads
- **Segmented reads (S-reads):** Number of generated S-reads
- **Mean length of S-reads:** Mean read length of generated S-reads
- **Percent of reads with full arrays:** Percentage of input HiFi reads containing all adapter sequences in the order listed in the segmentation adapter FASTA file
- **Mean array size:** Mean number of fragments (or S-reads) found in input reads


Histogram distribution of the number of S-reads per HiFi read. (Example Revo system data shown.)

<sup>1</sup> HiFi read lengths, reads/data per SMRT Cell and other sequencing performance results can vary depending on DNA sample quality, insert size, P1 loading performance & movie time. Note: Refer to **SMRT Link MAS-Seq troubleshooting guide (102-994-400)** for example performance metrics typically achievable with Kinnex libraries under optimal P1 loading conditions. For Sequel IIe systems, we recommend aiming for ~60 – 80% P1 loading. For Revo system, we recommend aiming for ~50 – 70% P1 loading.

# Example SMRT Link Read Segmentation data utility processing results<sup>1</sup> for Kinnex full-length RNA library prepared with human UHRR sample (cont.)


SMRT Link Read Segmentation data utility job report – Length of Reads and S-read Length Distribution

### Length of Reads



Histogram distribution of the number of HiFi reads by read length, in base pairs. (Example Revio system data shown.)

### S-read Length Distribution



Histogram distribution of the number of S-reads by HiFi read length, in base pairs. (Example Revio system data shown.)

# Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample

## SMRT Link Iso-Seq analysis job report – Read Classification statistics

### Summary Metrics

Value	Analysis Metric
46,396,697	Reads
44,875,144	Reads with 5' and 3' Primers
44,784,675	Non-Concatamer Reads with 5' and 3' Primers
44,739,994	Non-Concatamer Reads with 5' and 3' Primers and Poly-A Tail (FLNC Reads)
1,981	Mean Length of FLNC Reads
12	Unique Primers
3,739,595	Mean Reads per Primer
4,379,484	Max. Reads per Primer
3,143,785	Min. Reads per Primer
1,521,553	Reads without Primers
96.22%	Percent Bases in Reads with Primers
96.72%	Percent Reads with Primers

Example Revio system data shown.

- **Reads:** Total number of CCS reads
- **Reads with 5' and 3' Primers:** Number of CCS reads with 5' and 3' cDNA primers detected
- **Non-Concatamer Reads with 5' and 3' Primers:** Number of nonconcatemer CCS reads with 5' and 3' primers detected
- **Non-Concatamer Reads with 5' and 3' Primers and Poly-A Tail (FLNC Reads):** Number of non-concatemer CCS reads with 5' and 3' primers and polyA tails detected. This is usually the number for full-length, nonconcatemer (FLNC) reads, unless polyA tails are not present in the sample
- **Mean Length of FLNC Reads:** Mean length of the non-concatemer CCS reads with 5' and 3' primers and polyA tails detected
- **Unique Primers:** Number of unique primers in the sequence
- **Mean Reads per Primer:** Mean number of CCS reads per primer
- **Max. Reads per Primer:** Maximum number of CCS reads per primer
- **Min. Reads per Primer:** Minimum number of CCS reads per primer
- **Reads without Primers:** Number of CCS reads without a primer
- **Percent Bases in Reads with Primers:** Percentage of bases in CCS reads in the sequence data that contain primers
- **Percent Reads with Primers:** Percentage of CCS reads in the sequence data that contain primers

# Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)

## SMRT Link Iso-Seq analysis job report – Read Classification statistics

### Primer Data

Bio Sample Name	Primer Name	CCS Reads	Mean Primer Quality	Reads with 5' and 3' Primers	Non-Concatamer Reads with 5' and 3' Primers	Non-Concatamer Reads with 5' and 3' Primers and Poly-A Tail
BioSample_1	IsoSeqX_bc01_5p--IsoSeqX_3p	3,648,886	99.6	3,648,886	3,641,475	3,638,289
BioSample_2	IsoSeqX_bc02_5p--IsoSeqX_3p	3,826,584	99.7	3,826,584	3,818,994	3,815,554
BioSample_3	IsoSeqX_bc03_5p--IsoSeqX_3p	4,296,994	99.7	4,296,994	4,288,437	4,284,413
BioSample_4	IsoSeqX_bc04_5p--IsoSeqX_3p	3,588,147	99.7	3,588,147	3,580,832	3,577,622
BioSample_5	IsoSeqX_bc05_5p--IsoSeqX_3p	4,543,731	99.7	4,543,731	4,534,661	4,530,441
BioSample_6	IsoSeqX_bc06_5p--IsoSeqX_3p	4,606,161	99.7	4,606,161	4,596,642	4,592,382
BioSample_7	IsoSeqX_bc07_5p--IsoSeqX_3p	4,009,556	99.7	4,009,556	4,001,451	3,997,812
BioSample_8	IsoSeqX_bc08_5p--IsoSeqX_3p	3,651,440	99.7	3,651,440	3,644,176	3,641,097
BioSample_9	IsoSeqX_bc09_5p--IsoSeqX_3p	4,390,535	99.7	4,390,535	4,381,510	4,377,425
BioSample_10	IsoSeqX_bc10_5p--IsoSeqX_3p	4,049,234	99.7	4,049,234	4,041,043	4,037,382
BioSample_11	IsoSeqX_bc11_5p--IsoSeqX_3p	3,323,541	99.7	3,323,541	3,316,702	3,313,724
BioSample_12	IsoSeqX_bc12_5p--IsoSeqX_3p	3,372,840	99.7	3,372,840	3,365,834	3,362,745
Bio Sample 4	No Primer	1,603,738	0.0	0	0	0

Example Revio system data shown.


- **Bio Sample Name:** Name of the biological sample associated with the primer
- **Primer Name:** A string containing the pair of primer indices associated with this biological sample
- **CCS Reads:** Number of CCS reads associated with the primer
- **Mean Primer Quality:** Mean primer quality associated with the primer
- **Reads with 5' and 3' Primers:** Number of CCS reads with 5' and 3' cDNA primers detected
- **Non-Concatemer Reads with 5' and 3' Primers:** Number of non-concatemer CCS reads with 5' and 3' primers detected
- **Non-Concatemer Reads with 5' and 3' Primers and Poly-A Tail:** Number of non-concatemer CCS reads with 5' and 3' primers and polyA tails detected. This is usually the number for full-length, non-concatemer (FLNC) reads, unless polyA tails are not present in the sample.




# Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)

SMRT Link Iso-Seq analysis job report – Read Classification statistics


## Primer Read Statistics



**Number Of Reads Per Primer:** Maps the number of reads per primer, sorted by primer ranking



**Primer Frequency Distribution:** Maps the number of samples with primers by the number of reads with primers




**Mean Read Length Distribution:** Maps the read mean length against the number of samples with primers

# Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)



SMRT Link Iso-Seq analysis job report – Read Classification statistics

## Primer Quality Scores



**Primer Quality Score Distribution:** Histogram of primer scores

## Length of Full-length Non-Concatemer Reads




**Length of Full-Length Non-Concatemer Reads:** Per-sample histograms of the read length distribution of non-concatemer CCS reads with 5' and 3' primers and polyA tails detected



# Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)

## SMRT Link Iso-Seq analysis job report – Transcript Clustering statistics

### Summary Metrics

Sample Name ↓	Number of High-Quality Isoforms ↓
BioSample_1	650,495
BioSample_2	666,648
BioSample_3	712,671
	
BioSample_11	617,944
BioSample_12	621,931
All Samples	2,001,226

### Length of Consensus Isoforms



- **Sample Name:** Sample name for which the following metrics apply
- **Number of High-Quality Isoforms:** Number of consensus isoforms that have an estimated accuracy **above** the specified threshold

- **Length of Consensus Isoforms:** Per-sample histograms of the consensus isoform lengths and the distribution of isoforms exceeding a read length cutoff. Also includes a single histogram plot for all samples.

# Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)

## SMRT Link Iso-Seq analysis job report – Transcript Mapping and Classification statistics

### Summary Metrics (All samples)

Sample Name ↕	Total Unique Genes ↕	Total Unique Genes, filtered ↕	Total Unique Isoforms ↕	Total Unique Isoforms, filtered ↕
BioSample_1	127,560	18,576	317,862	127,517
BioSample_2	131,678	18,833	325,744	129,550
BioSample_3	143,056	19,089	348,790	135,483
BioSample_4	125,866	18,529	315,326	126,923
BioSample_5	148,422	19,451	360,453	138,918
BioSample_6	149,269	19,495	361,022	138,742
BioSample_7	136,195	18,867	336,140	132,612
BioSample_8	128,393	18,608	317,667	126,641
BioSample_9	145,921	19,291	353,923	136,903
BioSample_10	137,496	18,999	338,295	133,063
BioSample_11	120,079	18,308	301,435	122,466
BioSample_12	120,440	18,355	303,142	123,454

Example Revio system data shown.

- **Sample Name:** Sample name for which the following metrics apply
- **Total unique genes:** The total number of unique genes across all cells.
- **Total unique genes, filtered:** The total number of unique genes, after filtering out reads based on the SQANTI transcript filtering criteria.
- **Total unique isoforms:** The total number of unique isoforms across all cells
- **Total unique isoforms, filtered:** The total number of unique isoforms across all cells, after filtering out reads based on the SQANTI transcript filtering criteria.

# Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)

## SMRT Link Iso-Seq analysis job report – Transcript Mapping and Classification statistics

### Transcript Classification, filtered (All samples)

Category ††	Count ††	CAGE Detected ††	CAGE Detected, (%) ††	polyA Detected ††	polyA Detected, (%) ††
FSM	171662	83485	48.63%	92562	53.92%
ISM	257444	32719	12.70%	158888	61.71%
NIC	149131	97650	65.47%	78854	52.87%
NNC	106396	67432	63.37%	57440	53.98%
Antisense	1501	324	21.58%	923	61.49%
Fusion	3446	2021	58.64%	1901	55.16%
More junctions	83	49	59.03%	55	66.26%
Genic intron	0	0	0.00%	0	0.00%
Genic genomic	1166	524	44.93%	704	60.37%
Intergenic	3073	285	9.27%	2441	79.43%

Example Revio system data shown.<sup>1</sup>

- **Category:** Transcript classification<sup>2</sup> assigned by the classification and filtering tool `pigeon`, based on the [SQANTI3](#) software
- **Count:** The number of transcripts, after filtering out reads based on the SQANTI filtering criteria, in a specific classification
- **CAGE Detected:** The number of transcripts where the transcription start site falls within 50 bp of an annotated CAGE (Cap Analysis of Gene Expression) peak site
- **CAGE Detected, (%):** The percentage of transcripts where the transcription start site falls within 50 bp of an annotated CAGE peak site
- **polyA Motif Detected:** The number of transcripts where a known polyA motif is detected upstream of the transcription end site
- **polyA Motif Detected, (%):** The percentage of transcripts where a known polyA motif is detected upstream of the transcription end site

# Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)

SMRT Link Iso-Seq analysis job report – Transcript Mapping and Classification statistics


## Transcript Classification Plots

## Transcript Classification Plots, Filtered


### Isoform distributions across structural categories:

- Distribution of the % of isoforms by structural categories

Example Revo system data shown.




Filter out reads based on the SQANTI3 transcript filtering criteria




### Structural categories by isoform lengths:

- Histogram display of the number of isoforms by their length in kb and their structural category

Example Revo system data shown.



Filter out reads based on the SQANTI3 transcript filtering criteria




# Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)


SMRT Link Iso-Seq analysis job report – Transcript Mapping and Classification statistics

## Gene Saturation

Gene Saturation, all genes, filtered (All Samples)



Gene Saturation, known genes only, filtered (All Samples)



**Gene Saturation, all genes, filtered:** Saturation plot showing the level of gene saturation for **all genes**, after filtering out reads based on the SQANTI transcript filtering criteria

**Gene Saturation, known genes only, filtered:** Saturation plot showing the level of gene saturation, for unique **known genes only** (genes annotated in the reference annotation) per cell, after filtering out reads based on the SQANTI transcript filtering criteria



# Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)

## File Downloads tab

Edit Output File Name Prefix Example:analysis-Bio Sample 4-2110

File ↑	Size ↓	Type ↓
Non-passing reads, unaligned	4 GB	bam
Report read_segmentation	3 KB	JsonReport
SMRT Link Log	13 KB	log
Segmented Reads, passing, unaligned	65 GB	bam
UHRR_verification_DL_bc4 (demux1) (Bio Sample 4) Segmented Reads	20 KB	ConsensusReadSet
Unique mapped transcripts, GFF (All Samples)	895 MB	gff
Unique mapped transcripts, classification TXT (All Samples)	728 MB	txt
Unique mapped transcripts, filtered, GFF (All Samples)	599 MB	gff
Unique mapped transcripts, filtered, classification TXT (All Samples)	307 MB	txt
Unique mapped transcripts, filtered, junctions TXT (All Samples)	624 MB	txt
Unique mapped transcripts, junctions TXT (All Samples)	771 MB	txt

Refer to [SMRT Link user guide](#) for descriptions of downloadable output files

- These files are useful for **visualizing** isoform structures in Integrative Genomics Viewer (IGV) / UCSC genome browser and enable understanding of why an isoform is novel/known, etc.
  - GFF file containing unique mapped transcripts after filtering
  - Text file containing unique mapped transcript classifications against annotations, after filtering
  - Text file containing information about unique mapped transcript junctions, after filtering

Files shown in the File Downloads tab are available on the analysis results page. Additional files are also available on the SMRT Link server in the analysis output directory.



# Technical documentation & applications support resources

# Technical resources for Kinnex full-length RNA library preparation, sequencing & data analysis

## RNA sample preparation resources

- Note: The products below have not been tested or validated by PacBio but are listed here as examples of third-party kits used by other PacBio customers for isolating total RNA for PacBio RNA sequencing (Iso-Seq) applications

Kit type	Product name
mRNA isolation	Ambion Poly(A) Purist MAG Kit [ <a href="#">Link</a> ]
	Qiagen RNeasy Plus Kits [ <a href="#">Link</a> ]
Total RNA isolation	Sigma Spectrum Plant Total RNA Kit [ <a href="#">Link</a> ]
	iNtRON Easy Spin Total RNA [ <a href="#">Link</a> ]
	TRIzol Reagent can be used to isolate total RNA from tissues or cells, including lipid-rich and difficult samples [ <a href="#">Link</a> ]
RNA stabilization & storage	RNALater is an aqueous, nontoxic tissue storage reagent that rapidly permeates tissues to stabilize and protect cellular RNA [ <a href="#">Link</a> ]

# Technical resources for Kinnex full-length RNA library preparation, sequencing & data analysis (cont.)

## Kinnex full-length RNA library preparation literature & other resources

- Application note – Kinnex full-length RNA kit for isoform sequencing ([102-326-591](#))
- Procedure & checklist – Preparing Kinnex libraries using Kinnex full-length RNA kit ([103-238-700](#))
- Technical overview – Kinnex kits for single-cell RNA, full-length RNA and 16S rRNA sequencing ([103-343-700](#))
- Technical overview – Kinnex library preparation using Kinnex full-length RNA kit ([103-344-700](#))
- Video tutorial – SMRT Link Sample Setup and Run Design setup procedure for Kinnex kits [ [Link](#) ]
- Whitepaper – Bulk and single-cell isoform sequencing for human disease research ([102-326-576](#))

## Data analysis resources

- Application note – Bioinformatics tools for full length isoform sequencing ([102-326-593](#))
- SMRT Link MAS-Seq troubleshooting guide ([102-994-400](#))
- SMRT Link software installation guide [ [Link](#) ]
- SMRT Link user guide [ [Link](#) ]
- SMRT Tools reference guide [ [Link](#) ]
- Video tutorial – Read Segmentation and Iso-Seq workflow in SMRT Link [ [Link](#) ]

# Technical resources for Kinnex full-length RNA library preparation, sequencing & data analysis (cont.)

## Publications

- Schertzer, M.D. et al. (2023) Cas13d-mediated isoform-specific RNA knockdown with a unified computational and experimental toolbox. BioRxiv preprint [ [Link](#) ]
- Al'Khafaji, A.M. et al. (2023) High-throughput RNA isoform sequencing using programmable cDNA concatenation. Nature biotechnology. [ [Link](#) ]

## Webinars

- PacBio Iso-Seq social club webinar (2022) – Introduction to Iso-Seq method [ [Link](#) ]
- PacBio Iso-Seq social club webinar (2022) – SQANTI3 for isoform classification and annotation [ [Link](#) ]
- PacBio Iso-Seq social club webinar (2022) – TappAS for isoform differential expression analysis [ [Link](#) ]

## Example PacBio data sets

Application	Dataset	Data type	PacBio system
Kinnex full-length RNA sequencing	Homo sapiens – universal human reference RNA (UHRR) [ <a href="#">Link</a> ]	HiFi long read	Sequel II & Revio systems
	Homo sapiens – HG002 [ <a href="#">Link</a> ]	HiFi long read	Revio system
	Homo sapiens – Heart [ <a href="#">Link</a> ]	HiFi long read	Revio system
	Homo sapiens – Cerebellum [ <a href="#">Link</a> ]	HiFi long read	Revio system



[www.pacb.com](http://www.pacb.com)

Research use only. Not for use in diagnostic procedures. © 2024 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](http://pacb.com/license). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, and PureTarget are trademarks of PacBio.