



## Introduction

The number of useable bases per SMRT Cell can exceed 1 Gb with the latest P6-C4 chemistry and 6-hour movies. For applications such as microbial sequencing, targeted sequencing, Iso-Seq full-length isoform sequencing and Roche NimbleGen's target enrichment method, current PacBio RS II SMRT Cell yields could be in excess relative to project requirements.


To this end, barcoding is a viable option for multiplexing samples. For microbial sequencing, multiplexing can be accomplished by tagging sheared genomic DNA during library construction with SMRTbell barcoded adapters. Results from 2-, 4-, 6- and 8-plex will be presented.

For HLA typing, full-length HLA genes as large as 5 kb may be barcoded during amplification or during SMRTbell library construction. The preferred barcoding strategy depends on the user's existing workflow and flexibility for changing and/or updating existing workflows. Five Class I/II genes (3.3 – 5.8 kb) x 96 patients can be multiplexed and typed in up to 4 SMRT Cells.

For Iso-Seq analysis, barcodes are incorporated during first-strand synthesis using oligo-dT tailed with 16-bp barcodes. Six maize RNA samples were barcoded, pooled and constructed into a single SMRTbell library.

Finally for target enrichment using NimbleGen's SeqCap EZ Target Enrichment method, linear barcodes are incorporated before the capture step. Results from a 12-plex target enrichment experiment will be presented.


## Multiplexing for Microbial Assembly



**Figure 1:**


- Workflow for multiplexing microbial genomes
- Barcodes are added to the SMRTbell adapters.
- Recommended barcoded adapters available from PacBio.

## MICROBIAL WHOLE GENOME ASSEMBLY OF MULTIPLEXED *H. PYLORI* STRAINS



**Figure 2:** N50 contig vs Multiplex level of sequenced *H. pylori* strains using 10 kb library preparations, P6-C4 chemistry, and sequenced on the PacBio RS II. Increasing P1 loading helped achieve higher multiplex level to 6, possibly 8.


## Targeted Sequencing Multiplexing



**Figure 3:** Barcoding options for targeted sequencing. [www.pacb.com/wp-content/uploads/2015/09/ProductNote-Barcoded-Adapters-Barcoded-Universal-Primers.pdf](http://www.pacb.com/wp-content/uploads/2015/09/ProductNote-Barcoded-Adapters-Barcoded-Universal-Primers.pdf)

- Barcoded Universal Primers:** Barcode can be incorporated into the amplicon via a two-step tailed primer approach. Barcodes are commercially available from PacBio.
- Barcoded Adapters:** Barcodes are incorporated during ligation with barcoded adapters. Barcodes are commercially available from PacBio.
- Locus-specific primers tailed with barcodes.** Primers may be ordered from any oligo synthesis providers. The first 96 barcodes out of 384 available sequences are available: [www.pacb.com/wp-content/uploads/PacBio-PCR-Primer-Barcodes-0001-to-0096-IDT-Template.xlsx](http://www.pacb.com/wp-content/uploads/PacBio-PCR-Primer-Barcodes-0001-to-0096-IDT-Template.xlsx).


## MULTIPLEXING FOR HLA SEQUENCING



**Figure 4:**

- GenDx HLA-A, HLA-B, HLA-C, HLA-DRB1 and HLA-DQB1 are amplified and pooled for barcoding. All 96 wells are subsequently pooled into a single tube for SMRTbell library construction.
- Recommendations for multiplexing NGS-go amplicons. [www.pacb.com/applications/targeted-sequencing/hla/](http://www.pacb.com/applications/targeted-sequencing/hla/)

## Iso-Seq Application Multiplexing




**Figure 5:** Barcodes are incorporated into the cDNA during first-strand synthesis. Double-stranded cDNAs are pooled during SMRTbell library construction.

## Isoform Sequencing Multiplexing – cont.

### GENOME ANNOTATION OF MULTIPLEXED MAIZE TISSUES


**Figure 6:** Bioanalyzer traces of barcoded cDNA from 6 Maize tissues prior to pooling.



**Figure 7:**

- Full-length read distribution from pooled maize library.
- Full-length yield summary per fraction per sample.


## Target Enrichment Multiplexing




**Figure 8:** Target enrichment multiplex workflow. Sheared DNA samples are tagged with linear barcoded adapters. Workflow in development.

For information: [www.pacb.com/products-and-services/consumables/target-enrichment/](http://www.pacb.com/products-and-services/consumables/target-enrichment/)

## MULTIPLEXING WITH NIMBLEGEN NEUROLOGY PANEL



**Figure 9:** Multiplex of 12 samples. Number of mapped subreads of 5 kb target enrichment using Nimblegen's neurology panel.



**Figure 10:** A 60 kb fragment showing intron-exon regions of the NPC1 gene.

## Conclusions

- Variety of barcoding methods available for multiplexing samples for microbial whole genome sequencing, HLA typing, isoform sequencing, and targeted sequencing applications with PacBio Systems.
- Barcoding increases efficiency and sample throughput thus reducing the cost of library preparation and sequencing.

*The authors would like to acknowledge and thank everyone Involved with the data generation for this poster.*