

Description

The Magbead kit v2 contains reagent supplies to capture DNA polymerase and template complex for loading onto the PacBio systems. This kit consists of the following components: Magbeads (100-125-900), Magbead binding buffer kit (101-046-400), and Magbead wash buffer kit (101-046-500).

Kit configuration

Configuration

The Magbead kit v2 (100-676-500) includes:

- Magbeads (P/N 100-125-900)

#	Component	Part number	Qty	Color	Volume
1	Magbeads	100-125-900	1	N/A	2.88 mL

- Magbead binding buffer kit v2 (P/N 101-046-400)

#	Component	Part number	Qty	Color	Volume
1	Magbead binding buffer v2	100-672-100	9	red	1.7 mL

- Magbead wash buffer kit v2 (P/N 101-046-500)

#	Component	Part number	Qty	Color	Volume
1	Magbead wash buffer v2	100-672-200	9	yellow	1.7 mL

Quantities

Varies – refer to kit configuration section

Storage and handling requirements

Magbeads (100-125-900) should be stored refrigerated between 2 to 8 °C upon arrival. DO NOT FREEZE.

Magbead binding buffer kit v2 (101-046-400) and Magbead wash buffer kit v2 (101-046-500) should be stored frozen between -15 to -25 °C.

Refer to protocols for additional recommendations for reagent handling.

Shelf life

Refer to expiration date on kit label.

Intended use

Research use only. Not for use in diagnostic procedures.

Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, and SBB are trademarks of PacBio. All other trademarks are the sole property of their respective owners.