

Description

Onso library conversion kit supports reagents to enable conversion of Illumina-style P5/P7 libraries into Onso-compatible libraries. This kit contains 2x PCR master mix and conversion primers.

Kit configuration & layout

Configuration

#	Component	Part number	Qty	Color	Volume
1	Conversion primer mix	102-554-700	1	Light Blue	550µL
2	PCR master mix (2x)	102-552-700	1	Green	1200 µL
3	PCR master mix (2x)	102-552-700	1	Green	1200 µL

Layout



Quantities

96 reactions

Storage and handling requirements

Product should be stored frozen between -15 to -25 °C.

Refer to protocols for additional recommendations for reagent handling.

Shelf life

Refer to expiration date on kit label.

Intended use

Research use only. Not for use in diagnostic procedures.

Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, and Onso are trademarks of PacBio.