Description

The Onso library quant kit reagents enable accurate quantification of A/P1-tailed libraries. This kit contains qPCR mix, standards, dilution buffer, ROX dye (hi/lo) and PCR primers.

Kit configuration & layout

Configuration

#	Component	Part number	Qty	Color	Volume
1	Library quant primer mix	102-553-100	1	Light Blue	500 μL
2	Rox (high)	102-552-900	1	Green	100 μL
3	Rox (low)	102-553-000	1	Red	100 μL
4	Library quant DNA standard 1	102-553-300	1	Red	120 µL
5	Library quant DNA standard 2	102-553-400	1	Orange	120 µL
6	Library quant DNA standard 3	102-553-500	1	Yellow	120 µL
7	Library quant DNA standard 4	102-553-600	1	Green	120 µL
8	Library quant DNA standard 5	102-553-700	1	Blue	120 µL
9	Library quant DNA standard 6	102-553-800	1	Purple	120 µL
10	Library quant master mix	102-647-800	1	Clear	7500 µL
11	10x dilution buffer	102-647-900	1	Clear	7500 µL
12	10x dilution buffer	102-647-900	1	Clear	7500 µL

Storage and handling requirements

Product should be stored frozen between -15 to -25 °C.

Refer to protocols for additional recommendations for reagent handling.

Shelf life

Refer to expiration date on kit label.

Intended use

Research use only. Not for use in diagnostic procedures.

Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our website.

Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, and Onso are trademarks of PacBio.

Layout



Quantities

500 reactions



Kit insert PN: 102-851-900 Rev: 02