

Description

The Revio SPRQ polymerase kit 96 contains reagents for binding sequencing polymerase to SMRTbell[®] libraries. The kit includes 96 reactions, each supporting one Revio SMRT[®] Cell.

Kit configuration & layout

Configuration

#	Component	Part number	Qty	Color	Volume
1	Annealing buffer 96	103-287-400	1	light blue	1450 µL
2	Standard sequencing primer 96	103-287-500	1	light green	1450 µL
3-5	Polymerase buffer 96	103-287-600	3	yellow	1700 µL
6	Sequencing polymerase 96	103-512-900	1	purple	357 µL
7	Dilution buffer 96	103-288-000	1	blue	1505 µL
8	Sequencing control	103-508-800	1	red	27 µL
9-12	Loading buffer	103-485-400	4	green	1460 µL

Layout



Quantities

96 reactions

Storage and handling requirements

Product should be stored frozen between -15 to -25°C.

Refer to protocols for additional recommendations for product handling.

Shelf life

Refer to expiration date on kit label.

Intended use

Research use only. Not for use in diagnostic procedures.

Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, SPRQ, and Vega are trademarks of PacBio.