

Description

The Revio SPRQ sequencing plate – 1 rxn contains reagents for sequencing 1 Revio SMRT[®] Cell on the Revio system. Foil seals to prevent sample evaporation are included.

Kit configuration

Configuration

#	Component	Part number	Qty	Color	Volume
1	Revio SPRQ sequencing plate	103-496-800	1	N/A	N/A
2	Sample well seals	102-449-500	2	N/A	N/A

Quantities

1 reaction

Storage and handling requirements

Product should be stored frozen between -15 to -25 °C.

Refer to protocols for additional recommendations for product handling.

Shelf life

Refer to expiration date on kit label.

Intended use

Research use only. Not for use in diagnostic procedures.

Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, SPRQ, and Vega are trademarks of PacBio.