

Description

The Revio SPRQ sequencing plate contains reagents for sequencing 4 Revio SMRT[®] Cells on the Revio system. Foil seals to prevent sample evaporation are included.

Kit configuration

Configuration

| # | Component | Part number | Qty | Color | Volume |
|---|-----------------------------|-------------|-----|-------|--------|
| 1 | Revio SPRQ sequencing plate | 103-496-700 | 1 | N/A | N/A |
| 2 | Sample well seals | 102-449-500 | 2 | N/A | N/A |

Quantities

4 reactions

Storage and handling requirements

Product should be stored frozen between -15 to -25 °C.

Refer to protocols for additional recommendations for product handling.

Shelf life

Refer to expiration date on kit label.

Intended use

Research use only. Not for use in diagnostic procedures.

Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, SPRQ, and Vega are trademarks of PacBio.