

## Description

The SMRT Cell 8M single-use tray contains 1 SMRT Cell 8M for sequencing on the Sequel<sup>®</sup> II and IIe systems.

## Kit configuration

### Configuration

#	Component	Part number	Qty	Color	Volume
1	SMRT Cell 8M	N/A	1	NA	N/A

### Quantities

1 SMRT Cell tray containing 1 SMRT Cell

## Storage and handling requirements

Product should be stored at ambient temperature between 8 to 40°C.

Refer to protocols for additional recommendations for reagent handling.

## Shelf life

Refer to expiration date on kit label.

## Intended use

Research use only. Not for use in diagnostic procedures.

## Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

## Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](http://pacb.com/license). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, and Onso are trademarks of PacBio.