# SMRTbell<sup>®</sup> enzyme cleanup kit 2.0

#### Description

The SMRTbell enzyme cleanup kit 2.0 contains reagents for enzymatic removal of incomplete or damaged SMRTbell molecules (18 reactions).

# Kit configuration & layout

#### Configuration

#	Component	Part number	Qty	Color	Volume
1	SMRTbell enzyme clean up buffer 2.0	101-932-700	1	white	130 µL
2	SMRTbell enzyme clean up mix 2.0	101-932-900	1	blue	130 µL

#### Layout



### Quantities

18 reactions

### Storage and handling requirements

Product should be stored frozen between -15 to -25 °C.

Refer to protocols for additional recommendations for reagent handling.

## Shelf life

Refer to expiration date on kit label.

#### Intended use

Research use only. Not for use in diagnostic procedures.

#### Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our <u>website</u>.

#### Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at <u>pacb.com/license</u>. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, and Kinnex are trademarks of PacBio.

