

Description


The SMRTbell enzyme cleanup kit 1.0 contains 4 enzymes and buffer for enzymatic removal of incomplete or damaged SMRTbell molecules (18 reactions).

Kit configuration & layout

Configuration

#	Component	Part number	Qty	Color	Volume
1	10X clean up buffer	101-741-600	1	white	1000 µL
2	Enzyme A	101-741-100	1	yellow	80 µL
3	Enzyme B	101-741-700	1	blue	25 µL
4	Enzyme C	101-741-400	1	red	25 µL
5	Enzyme D	101-741-500	1	orange	40 µL

Layout



Quantities

18 reactions

Storage and handling requirements

Product should be stored frozen between -15 to -25 °C.

Refer to protocols for additional recommendations for reagent handling.

Shelf life

Refer to expiration date on kit label.

Intended use

Research use only. Not for use in diagnostic procedures.

Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, and Kinnex are trademarks of PacBio.