

Description

Reagents to complete genomic DNA size selection for HiFi sequencing by depleting DNA <10 kb in length (24 reactions).

Kit configuration

Configuration

#	Component	Part number	Qty	Color	Volume
1	Buffer SRE	102-280-900	1	clear	1.7 mL
2	Buffer LTE	102-277-100	1	clear	5 mL

Quantities

24 rxn

Storage and handling requirements

Product should be stored ambient between 15 to 30 °C.

Refer to protocols for additional recommendations for reagent handling.

Shelf life

Refer to expiration date on kit label.

Intended use

Research use only. Not for use in diagnostic procedures.

Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, and PureTarget are trademarks of PacBio.