

Description

The Sequencing primer v3 is used to anneal to SMRTbell® templates. Sequencing primer v3 is compatible with SMRTbell templates and should only be used for SMRTbell libraries for Magbead loading on the Sequel® systems.

Kit configuration

Configuration

#	Component	Part number	Qty	Color	Volume
1	Sequencing primer v3	101-374-000	1	blue	25 µL

Quantities

10 reactions

Storage and handling requirements

Product should be stored frozen between -15 to -25 °C.

Refer to protocols for additional recommendations for reagent handling.

Shelf life

Refer to expiration date on kit label.

Intended use

Research use only. Not for use in diagnostic procedures.

Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](https://www.pacb.com/license). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, and SBB are trademarks of PacBio. All other trademarks are the sole property of their respective owners.